



**See Some Important Members of
AuPreP Family**

AuPreP Citations

AuPreP Extra Mile First Strand cDNA System

AuPreP Extra Mile First Strand cDNA System is a high quality, robust, stable and optimized premixed system with RNase H minus Reverse Transcriptase including the requisite components of a cDNA assembly system. The synthesized cDNA is of very high quality and length, for your onward applications viz. PCR, cDNA libraries

Cat. No. : AUP-EXM-50

Storage: - 20°C

Shipping Conditions: Dry Ice

System Details:

AuPreP Extra Mile First Strand cDNA System is a stable, high quality system with RNase H minus Reverse Transcriptase (point mutant lacking RNase H activity). Using this powerful system, RNA is reverse transcribed from RNA into single stranded cDNA, which further, can be amplified by the Polymerase Chain Reaction (PCR) to analyze gene expression. The site of the new cDNA strand synthesized by Reverse Transcriptase (RT) enzyme is determined by the type of primer used by the user viz. Oligo (dT) primer, random primer, or a sequence-specific primer. The synthesized cDNA is of very high quality and length and can then be used as a template for PCR.

The Extra Mile mix supplied in the system contains buffer, Point mutant RNase H minus reverse transcriptase, and RNase inhibitor sufficient for 50 reactions.

Salient Features:

- **AuPreP Extra Mile First Strand cDNA System produces very high quality and long cDNA fragments, which shows excellent performance in RT-PCR experiments and cDNA library construction.**
- **AuPreP Extra Mile First Strand cDNA System shows excellent performance with sequence specific primers, Oligo (dT) primers as well as with random primers.**
- **AuPreP Extra Mile First Strand cDNA System shows superior performance due to RNase H minus Reverse Transcriptase and very high quality ultra pure**



LIFE TECHNOLOGIES (INDIA) PVT. LTD.

your molecular & cell technology partner

Bringing zest to your research



System Components

The following components are supplied with the AuPreP Extra Mile First Strand cDNA System:

Extra Mile Mix Contains; Reverse transcriptase, RNase inhibitor and >99% pure dNTP's in buffer solution	400µl
DEPC-Treated Water	1.5ml
Random Hexamer Primer , 40µM	50µl
Primer Mix Human G3PDH amplimers, 10µM each	50µl
Oligo (dT)₂₀ Primer , 40µM	50µl
DTT Solution , 100mM	100µl
Control RNA Human, Total RNA, 1µg/µl	25µl

NB: Material not supplied with the system include PCR components viz. Amplification primers specific to your target cDNA, DNA Polymerase for PCR, Thermal Cycler, Microcentrifuge, 42°C and 70°C water baths or heat blocks.

Important Notes

Please read the notes below before starting the procedure.

- *Wear gloves and use sterile pipette tips to avoid RNase contamination and degradation of RNA.*
- *High quality RNA preparation is critical for the synthesis of cDNA for the synthesis of cDNA for PCR. RNA should have a A₂₆₀/A₂₈₀ ratio of 1.7 or higher and the integrity and purity should be evaluated. The RNA should be stored at -70°C or below.*
- *cDNA Priming: AuPreP Extra Mile First Strand cDNA System allows the user to choose the desired primer for cDNA synthesis.*
 1. *Oligo (dT)₂₀ or random primer: the entire population of mRNA molecules is converted into cDNA by priming with Oligo (dT) or random primer. Both primers are provided in the kit, but the use of Oligo (dT) primer is recommended. The random priming of cDNA may be beneficial when the reverse transcriptase fails to fully transcribe an mRNA template or if secondary structures exist.*


Contact Us: **Life Technologies (India) Pvt. Ltd .**, 306, Agarwal City Mall, opposite M2K Pitampura , Delhi-110034 (India) Tel # +91-11-4220-8000;4220-8111;4220-8222 Fax# +91-11-4220-8444,Mobile# +91-98105-21400, Email - customerservice@lifetechindia.com Web: www.lifetechindia.com

Last Updated on: 04/04/2009



2. *Gene Specific Primer: Primer design, concentration and annealing temperature should be evaluated before use.*

Flow Chart Procedure for cDNA Synthesis





LIFE TECHNOLOGIES (INDIA) PVT. LTD.

your molecular & cell technology partner

Bringing zest to your research



Subsequent Suggested Protocol

1. For PCR Amplification:

- Dilute the reaction mixture to a final volume of 100µl by adding 80µl of DEPC-Treated Water. Vortex gently and centrifuge.
Note: After thawing frozen samples, vortex and spin briefly before use.
- The PCR amplification parameters will vary depending on the specific primers, template DNA and thermal cycler used.
For each 50µl PCR reaction, use 5-10µl of the diluted cDNA.
- Analyze the RT-PCR product by electrophoresis on an agarose gel. Expected results will appear as a single band of size determined by the PCR primers used.

2. Controls:

- AuPreP Extra Mile First Strand cDNA System includes Control RNA and PCR Primers. For cDNA synthesis, use 1µl of Control RNA (human, total RNA, 1µg/µl). Then use the Primer Mix (G3PDH) for the PCR amplification.
- Use the following PCR protocol for control primers and cDNA template.

- **PCR reaction**

Sterile H ₂ O -----	33.6µl
10X PCR Buffer -----	5µl
dNTP's Mix (2.5mM each) -----	4µl
Primer Mix (G3PDH) -----	2µl
Taq DNA Polymerase (5u/µl) -----	0.4µl
Diluted cDNA -----	5µl

Total:	50µl
---------------	-------------

- **Amplification Parameters**

Initial Denaturation-----94°C-----	2 minutes	} 30 cycles
Denaturation-----94°C-----	45 seconds	
Annealing-----60°C-----	45 seconds	
Extension-----72°C-----	2 minutes	
Final Extension-----72°C-----	7 minutes	

- Upon gel electrophoresis (1.8% agarose gel), a 983 bp fragment should be observed.



<u>Other AuPreP™ DNA/RNA Kits</u>	<u>Other Related Products</u>
AuPreP™ Plasmid Maxi Kit	AuPreP Oligos (High Affinity Purified Oligo synthesis available in different scales, purifications & modifications)
AuPreP™ Plasmid Midi Kit	AuPreP TaQ DNA Polymerase (Ultrapure, Ultra-stable & Ultra-sensitive Taq DNA Polymerase)
AuPreP™ SPIN™ SPIN Miniprep Kit	AuPreP Hotstart TaQ DNA Polymerase (Robust Polymerase for Hotstart PCR assays)
AuPreP™ Blood Genomic DNA Maxi	AuPreP Super Fidelity TaQ DNA Polymerase (High fidelity Polymerase produces blunt ended amplicons upto 5Kb)
AuPreP™ Blood Genomic DNA Extraction Midi Kit	PCR Doctor - (PCR enhancer for AuPreP Hotstart Taq or Super Fidelity Taq especially designed for GC/AT/Dirty/Difficult Templates)
AuPreP™ GEN ^{bt} DNA Extraction Kit	AuPreP Longjump Polymerase (Robust Long Polymerase for templates > 4kb to 18kb+ for challenging PCRs)
AuPreP™ DNA easy Plant Maxi kit	AuPreP Red PCR Master Mix (2x Master mix with Red Dye without Enhancer)
AuPreP™ DNA easy Plant Mini Kit	AuPreP DIAMOND MASTER-MIX (2x Mastermix with PCR Enhancer & Stabilizer without tracking dyes)
AuPreP™ PCR Purification Kit	AuPreP DIAMOND DOUBLE DYE MASTERMIX (2x Mastermix with PCR Enhancer, Stabilizer & tracking dyes)
AuPreP™ Plant RNA Maxi Kit	AuPreP DNA Extraction System (A fast Reagent for pure genomic DNA isolation for down stream applications)
AuPreP™ Plasmid Maxi Kit	AuPreP RNA Extraction System (for Purest & High Quality RNA extraction with simple cost effective protocol)
AuPreP™ RNA Easy Midi Kit	AuPreP Gold cDNA Synthesis Kit (Highly Cost effective cDNA Synthesis Kit using RT with reduce Rnase H activity)
AuPreP™ RNA ^m Mini Kit	AuPreP Gold RT-PCR Combo Kit (2 step RT-PCR protocol with tracking Dye)
AuPreP™ RNV ^m Viral RNA Extraction Miniprep Kit	AuPreP Extra Mile First Strand cDNA System (Premium cDNA Synthesis Kit using RT with point mutant Rnase H minus activity)
	Novascript III RNase H⁻ RT (Premium Ultra-stable Rnase H minus RT for long high quality cDNA construction)
	Novascript III single step RT-PCR System (Premium 1step RT-PCR system using Novascript & AuPreP Hotstart DNA Polymerase)
	AuPreP Random Primer labeling Mix System (Premixed solution for the labeling of DNA with radiolabeled dCTP using random sequence oligonucleotides)