

*n*Taq

 Cat.#
 Size
 Conc.

 P025B
 250 units
 5 units/µl

 P050B
 500 units
 5 units/µl

Store at -20℃

Supplied with: 10X nTaq Buffer (Mg2+ free)

dNTP Mixture 25 mM MgCl₂ Sterile water

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposité M2K Pitampura, Delhi – 110034 (INDIA).

Ph: +91-11-42208000, 42208111, 42208222

Mobile: +91-9810521400 Fax: +91-11-42208444

Email: customerservice@lifetechindia.com

Web: www.lifetechindia.com

Product description

A gene encoding the Thermus aquaticus (Taq) DNA polymerase was cloned and expressed in *E. coli*, and the enzyme was purified to homogeneity. The purified Taq polymerase (nTaq) has optimal activity at high temperatures (72° C), which helps amplify secondary-structured regions. nTaq DNA Polymerase has greatly reduced activity of 5° — 3° exonuclease activity. As a result, formation of erroneously amplified products caused by the 5° — 3° exonuclease activity of wild type polymerase is effectively prevented. In addition, high GC-content or secondary-structured regions can be efficiently amplified due to the genetic modifications in the catalytic domain of the polymerase.

Characteristics

- Molecular weight: 94 kDa
- Error rate: 2.4 X 10⁻⁵
- Thermal stability: Half life of 40 min at 95℃
- A-tail formation at 3' ends of amplified duplex DNA.

Applications

- Amplification of DNA fragments shorter than 3 kb (Suitable for general PCR analysis)
- Amplification of cDNA and genomic DNA,
- Primer extension
- Colony PCR

For Research Use Only, Not for use in diagnostic procedures.

ISO9001 [ISO14001] ISO13485

- Multiplex PCR
- Labeling of DNA fragments with radioactive-isotopes
- Nucleotide sequencing

Quality control

- Purity: >99% on SDS-PAGE
- Endonuclease-free
- Exonuclease-free
- RNase-free
- Inhibitor-free

Unit definition

One unit is defined as the amount of enzyme required to incorporate 10 nmol of dNTP into acid insoluble materials in 30 min at 74°c in a 50- μ l reaction mixture (20 mM Tris-HCl/pH 8.8, 50 mM KCl, 2.5 mM MgCl₂, 10 mM β -mercaptoethanol, 12.5 μ g of calf thymus DNA).

Storage buffer

20 mM Tris-HCl (pH 7.9), 100 mM KCl, 0.1 mM EDTA, 1 mM DTT, 0.5% NP-40, 0.5% Tween-20, 50% glycerol.

- Labeling of DNA fragments with radioactive-isotopes

10X nTaq Buffer (Mg2+ Free)

Mg²⁺ plus Buffer: containing 15 mM MgCl₂ Mg²⁺ free: MgCl₂ free

Life Technologies™

Life Technologies™

Standard PCR conditions

10X nTag Buffer (Mg2+ free)

dNTP Mixture (2 mM each)

Template DNAd (0.1~500 ng/µl)

nTag DNA Polymeraseb (5 units/µl)

^aAssemble the reaction mixture on ice

the lower Tm of the two primers used.

bAdd the PCR polymerase at the final step

cAdjust Ma2+ to 1.5~5 mM with 25 mM MaCla

dPlasmid DNA, 0.1 ng~30 ng; genomic DNA, 50 ng~500 ng

95℃

95℃

72℃

72℃

When cycles are over, keep the reaction mixture at 4°C; may

^aRecommended annealing temperatures is 5 to 10℃ below

add 10 mM EDTA until use to prevent DNA degradation.

55℃~65℃

25~35 times

- PCR mixture^a

MgCl₂c (25 mM)

Sterile water

- PCR cycle

Primer 1 (5 pmole/ul)

Primer 2 (5 pmole/µl)

Initial denaturation

Number of cycles

Final elongation

Denaturation

Annealing^a

Elongation

Your Molecular & Cell Technology Partner

2 ul

2 µl

Xμl

1 µl

1 ul

1 µl

2 min

30 sec

5 min

30~60 sec

1 min/kb

up to 20 µl

0.2 µl

Your Molecular & Cell Technology Partner

Standard PCR conditions

- PCR mixture^a

10X nTaq Buffer (Mg ²⁺ free)	2 µl
nTaq DNA Polymerase ^b (5 units/µl)	0.2 µl
dNTP Mixture (2 mM each)	2 µl
$MgCl_2^c$ (25 mM)	Xμl
Template DNAd (0.1~500 ng/µl)	1 µl
Primer 1 (5 pmole/µl)	1 µl
Primer 2 (5 pmole/µl)	1 µl
Sterile water	up to 20 µl
^a Assemble the reaction mixture on ice	
bAdd the PCR polymerase at the final step	
cAdjust Mg ²⁺ to 1.5~5 mM with 25 mM MgCl ₂	
dPlasmid DNA, 0.1 ng~30 ng; genomic DNA, 50	ng~500 ng

- PCR cvcle

Initial departuration

	iriitiai deriaturation	950	2 11111		
	Denaturation	95℃	30 sec		
	Annealing ^a	55℃~65℃	30~60 sec		
	Elongation	72℃	1 min/kb		
	Number of cycles	25~35 times			
	Final elongation	72℃	5 mir		
When cycles are over, keep the reaction mixture at 4°C; may					
add 10 mM EDTA until use to prevent DNA degradation.					
^a Recommended annealing temperatures is 5 to 10℃ below					
	the lower Tm of the two primers used.				

*n*Tag

enzynomics

 Cat.#
 Size
 Conc.

 P025B
 250 units
 5 units/µl

 P050B
 500 units
 5 units/µl

Store at -20℃

Supplied with: 10X nTaq Buffer (Mg2+ free)

dNTP Mixture 25 mM MgCl₂ Sterile water

India Contact:

Life Technologies (India) Pvt. Ltd. 306, Aggarwal City Mall, Opposite M2K Pitampura,

Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222

Mobile: +91-9810521400 Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Web: www.lifetechindia.com

Product description

A gene encoding the Thermus aquaticus (Taq) DNA polymerase was cloned and expressed in *E. coli*, and the enzyme was purified to homogeneity. The purified Taq polymerase (nTaq) has optimal activity at high temperatures (72°c), which helps amplify secondary-structured regions. nTaq DNA Polymerase has greatly reduced activity of 5'—3' exonuclease activity. As a result, formation of erroneously amplified products caused by the 5'—3' exonuclease activity of wild type polymerase is effectively prevented. In addition, high GC-content or secondary-structured regions can be efficiently amplified due to the genetic modifications in the catalytic domain of the polymerase.

Characteristics

- Molecular weight: 94 kDa
- Error rate: 2.4 X 10⁻⁵
- Thermal stability: Half life of 40 min at 95℃
- A-tail formation at 3' ends of amplified duplex DNA.

Applications

- Amplification of DNA fragments shorter than 3 kb
- (Suitable for general PCR analysis)
 Amplification of cDNA and genomic DNA.
- Primer extension
- Colony PCR

For Research Use Only, Not for use in diagnostic procedures, [SO9001] [SO14001] [SO13485]

Unit definition

- Multiplex PCR

Ouality control

- Endonuclease-free

- Exonuclease-free

- RNase-free

- Inhibitor-free

- Nucleotide sequencing

- Purity: >99% on SDS-PAGE

One unit is defined as the amount of enzyme required to incorporate 10 nmol of dNTP into acid insoluble materials in 30 min at 74°c in a 50-µl reaction mixture (20 mM Tris-HCl/pH 8.8, 50 mM KCl, 2.5 mM MgCl₂, 10 mM β -mercaptoethanol, 12.5 μg of calf thymus DNA).

Storage buffer

20 mM Tris-HCl (pH 7.9), 100 mM KCl, 0.1 mM EDTA, 1 mM DTT, 0.5% NP-40, 0.5% Tween-20, 50% glycerol.

10X nTag Buffer (Mg²⁺ Free)

Mg²⁺ plus Buffer: containing 15 mM MgCl₂ Mg²⁺ free: MgCl₂ free