



## 2X TOPsimple™ DyeMIX-nTaq

Cat.#	Size
P510T	1 ml (0.5 ml X 2)
P525T	2.5 ml (0.5 ml X 5)

Store at 20°C

### India Contact:

**Life Technologies (India) Pvt. Ltd.**  
 306, Aggarwal City Mall, Opposite M2K Pitampura,  
 Delhi – 110034 (INDIA).  
 Ph: +91-11-42208000, 42208111, 42208222  
 Mobile: +91-9810521400  
 Fax: +91-11-42208444  
 Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com)  
 Web: [www.lifetechindia.com](http://www.lifetechindia.com)

### Product description

2X TOPsimple™ DyeMIX-nTaq is similarly formulated to the 2X TOPsimple™ PreMIX-nTaq except that it contains loading dyes for further convenience of use. Thus, the reaction mixtures after PCR cycles are ready for agarose gel electrophoresis. This product includes nTaq DNA Polymerase (Cat.# P725, P750), and is optimized for high efficiency PCR amplification with a high success rate.

### Characteristics

- 2-Dye system: Easy for gel electrophoresis (Xylene cyanol and Orange G)
- A-tail formation at 3' ends of amplified duplex DNA

### Applications

- Amplification of DNA fragments shorter than 3 kb (Suitable for general PCR analysis)
- Amplification of cDNA and genomic DNA.
- Primer extension
- Colony PCR
- Multiplex PCR

### Components

- nTaq DNA Polymerase: 0.2 units/μl
- nTaq buffer (containing 3 mM Mg<sup>2+</sup>)
- dNTP mixture: 0.4 mM each (2X)
- Stabilizer
- Dye: Xylene cyanol and Orange G

### Migration of dyes in agarose gel

In an ordinary agarose gel, xylene cyanol co-migrate with 4-kb DNA fragments, and Orange G with 50-bp DNA fragments.

### Use of the product

- To carry out PCR amplification, add template DNA and a pair of primers dissolved in water (10 μl) to the same volume (10 μl) of 2X TOPsimple™ DyeMIX-nTaq
- Completely thaw the premix solution prior to use.
- Frequent cycles of freezing and thawing reduce activity.
- Recommend using the 2X TOPsimple™ DyeMIX-nTaq as soon as possible once thawed.

### Storage

Stable up to 18 months at -20°C or 3 months at 4°C (Storage at -20°C is recommended).



### Standard PCR conditions

- PCR mixture<sup>a</sup>

2X TOPsimple™ DyeMIX-nTaq	10 μl
Template DNA <sup>b</sup> (0.1–500 ng/μl)	1 μl
Primer 1 (5 pmole/μl)	1 μl
Primer 2 (5 pmole/μl)	1 μl
Distilled water	up to 20 μl

<sup>a</sup>Assemble the reaction mixture on ice

<sup>b</sup>Plasmid DNA, 0.1 ng–30 ng; genomic DNA, 50 ng–500 ng

- PCR cycles

Initial denaturation	95°C	2 min
Denaturation	95°C	30 sec
Annealing <sup>a</sup>	55°C~65°C	30–60 sec
Elongation	72°C	1 min/kb
Number of cycles	25–35 times	
Final elongation	72°C	5 min

When cycles are over, keep the reaction mixture at 4°C; may add 10 mM EDTA until use to prevent DNA degradation.

<sup>b</sup>Recommended annealing temperatures is 5 to 10°C below the lower T<sub>m</sub> of the two primers used.



## 2X TOPsimple™ DyeMIX-nTaq

Cat.#	Size
P510T	1 ml (0.5 ml X 2)
P525T	2.5 ml (0.5 ml X 5)

Store at 20°C

### India Contact:

**Life Technologies (India) Pvt. Ltd.**  
 306, Aggarwal City Mall, Opposite M2K Pitampura,  
 Delhi – 110034 (INDIA).  
 Ph: +91-11-42208000, 42208111, 42208222  
 Mobile: +91-9810521400  
 Fax: +91-11-42208444  
 Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com)  
 Web: [www.lifetechindia.com](http://www.lifetechindia.com)

### Product description

2X TOPsimple™ DyeMIX-nTaq is similarly formulated to the 2X TOPsimple™ PreMIX-nTaq except that it contains loading dyes for further convenience of use. Thus, the reaction mixtures after PCR cycles are ready for agarose gel electrophoresis. This product includes nTaq DNA Polymerase (Cat.# P725, P750), and is optimized for high efficiency PCR amplification with a high success rate.

### Characteristics

- 2-Dye system: Easy for gel electrophoresis (Xylene cyanol and Orange G)
- A-tail formation at 3' ends of amplified duplex DNA

### Applications

- Amplification of DNA fragments shorter than 3 kb (Suitable for general PCR analysis)
- Amplification of cDNA and genomic DNA.
- Primer extension
- Colony PCR
- Multiplex PCR

### Components

- nTaq DNA Polymerase: 0.2 units/μl
- nTaq buffer (containing 3 mM Mg<sup>2+</sup>)
- dNTP mixture: 0.4 mM each (2X)
- Stabilizer
- Dye: Xylene cyanol and Orange G

### Migration of dyes in agarose gel

In an ordinary agarose gel, xylene cyanol co-migrate with 4-kb DNA fragments, and Orange G with 50-bp DNA fragments.

### Use of the product

- To carry out PCR amplification, add template DNA and a pair of primers dissolved in water (10 μl) to the same volume (10 μl) of 2X TOPsimple™ DyeMIX-nTaq
- Completely thaw the premix solution prior to use.
- Frequent cycles of freezing and thawing reduce activity.
- Recommend using the 2X TOPsimple™ DyeMIX-nTaq as soon as possible once thawed.

### Storage

Stable up to 18 months at -20°C or 3 months at 4°C (Storage at -20°C is recommended).



### Standard PCR conditions

- PCR mixture<sup>a</sup>

2X TOPsimple™ DyeMIX-nTaq	10 μl
Template DNA <sup>b</sup> (0.1–500 ng/μl)	1 μl
Primer 1 (5 pmole/μl)	1 μl
Primer 2 (5 pmole/μl)	1 μl
Distilled water	up to 20 μl

<sup>a</sup>Assemble the reaction mixture on ice

<sup>b</sup>Plasmid DNA, 0.1 ng–30 ng; genomic DNA, 50 ng–500 ng

- PCR cycles

Initial denaturation	95°C	2 min
Denaturation	95°C	30 sec
Annealing <sup>a</sup>	55°C~65°C	30–60 sec
Elongation	72°C	1 min/kb
Number of cycles	25–35 times	
Final elongation	72°C	5 min

When cycles are over, keep the reaction mixture at 4°C; may add 10 mM EDTA until use to prevent DNA degradation.

<sup>b</sup>Recommended annealing temperatures is 5 to 10°C below the lower T<sub>m</sub> of the two primers used.