

M- MLV Reverse Transcriptase

Cat.#	Size	Conc.
RT001S	10,000 units	200 units/μl
RT001M	10,000 units	200 units/μl
RT001L	50,000 units	200 units/μl
RT001H	50,000 units	1000 units/μl

Store at -20°C

Supplied with: 10X M-MLV RT Buffer
Sterile water (RNase free)

India Contact:

Life Technologies (India) Pvt. Ltd.
306, Aggarwal City Mall, Opposite M2K Pitampura,
Delhi – 110034 (INDIA).
Ph: +91-11-42208000, 42208111, 42208222
Mobile: +91-9810521400
Fax: +91-11-42208444
Email: customerservice@lifetechindia.com
Web: www.lifetechindia.com

Product description

Moloney Murine Leukemia Virus Reverse Transcriptase (M-MLV RT) is a RNA-directed DNA polymerase that can synthesize a complementary DNA (cDNA) strand from a primer hybridized to either RNA or single stranded DNA as a template. M-MLV Reverse Transcriptase is capable of synthesizing DNA longer than 5 kb from messenger RNA. RNase H activity of M-MLV reverse transcriptase is weaker than that of Avian Myeloblastosis Virus (AMV) reverse transcriptase.

Characteristics

- Molecular weight: 71 kDa
- Reaction temperature: 37°C or 42°C
- Heat inactivation: 70°C, 10 min

Applications

- Synthesis of first-strand cDNA
- Primer extension

Quality control

- Purity: >99% on SDS-PAGE
- Endonuclease-free
- Exonuclease-free
- RNase-free

For Research Use Only. Not for use in diagnostic procedures.

ISO9001 ISO14001 ISO13485

M- MLV Reverse Transcriptase

Cat.#	Size	Conc.
RT001S	10,000 units	200 units/μl
RT001M	10,000 units	200 units/μl
RT001L	50,000 units	200 units/μl
RT001H	50,000 units	1000 units/μl

Store at -20°C

Supplied with: 10X M-MLV RT Buffer
Sterile water (RNase free)

India Contact:

Life Technologies (India) Pvt. Ltd.
306, Aggarwal City Mall, Opposite M2K Pitampura,
Delhi – 110034 (INDIA).
Ph: +91-11-42208000, 42208111, 42208222
Mobile: +91-9810521400
Fax: +91-11-42208444
Email: customerservice@lifetechindia.com
Web: www.lifetechindia.com

Product description

Moloney Murine Leukemia Virus Reverse Transcriptase (M-MLV RT) is a RNA-directed DNA polymerase that can synthesize a complementary DNA (cDNA) strand from a primer hybridized to either RNA or single stranded DNA as a template. M-MLV Reverse Transcriptase is capable of synthesizing DNA longer than 5 kb from messenger RNA. RNase H activity of M-MLV reverse transcriptase is weaker than that of Avian Myeloblastosis Virus (AMV) reverse transcriptase.

Characteristics

- Molecular weight: 71 kDa
- Reaction temperature: 37°C or 42°C
- Heat inactivation: 70°C, 10 min

Applications

- Synthesis of first-strand cDNA
- Primer extension

Quality control

- Purity: >99% on SDS-PAGE
- Endonuclease-free
- Exonuclease-free
- RNase-free

For Research Use Only. Not for use in diagnostic procedures.

ISO9001 ISO14001 ISO13485

Your Molecular & Cell Technology Partner

Standard reaction conditions

10X M-MLV RT Buffer	2 μl
M-MLV Reverse Transcriptase (200 units/μl)	1 μl
dNTP Mixture (2 mM each)	1 μl
^a)Template	X μl
^b)Primer	1 μl
RNase Inhibitor (40 units/μl)	0.5 μl
Distilled water	up to 20 μl

^a) Prepare one of the following RNA template.

- Total RNA: 1 ng~5 μg
- Messenger RNA (mRNA): 1 ng~250 ng
- Specific RNA: 0.01 pg~0.5 μg

^b) Prepare one of the following primers.

- Oligo (dT)₁₈: 50~100 μM
- Random hexamer: 50~100 μM
- Specific primer: 15~20 pmol

→ An additional annealing step is recommended:

- if using oligo(dT)₁₈, incubate at 42°C for 5 min.
- if using random hexamer, incubate at 25°C for 10 min.

→ Incubate at 42°C for 60 min.

→ Incubate at 95°C for 5 min to inactivate the reaction.

Your Molecular & Cell Technology Partner

Standard reaction conditions

10X M-MLV RT Buffer	2 μl
M-MLV Reverse Transcriptase (200 units/μl)	1 μl
dNTP Mixture (2 mM each)	1 μl
^a)Template	X μl
^b)Primer	1 μl
RNase Inhibitor (40 units/μl)	0.5 μl
Distilled water	up to 20 μl

^a) Prepare one of the following RNA template.

- Total RNA: 1 ng~5 μg
- Messenger RNA (mRNA): 1 ng~250 ng
- Specific RNA: 0.01 pg~0.5 μg

^b) Prepare one of the following primers.

- Oligo (dT)₁₈: 50~100 μM
- Random hexamer: 50~100 μM
- Specific primer: 15~20 pmol

→ An additional annealing step is recommended:

- if using oligo(dT)₁₈, incubate at 42°C for 5 min.
- if using random hexamer, incubate at 25°C for 10 min.

→ Incubate at 42°C for 60 min.

→ Incubate at 95°C for 5 min to inactivate the reaction.

Unit definition

One unit is the amount of enzyme required to incorporate 1 nmol of dTTP into acid-insoluble materials using 0.4 mM poly(rA)-oligo(dT) as substrate at 37°C in 10 min.

Storage buffer

50 mM Tris-HCl (pH 7.6), 150 mM NaCl, 1 mM DTT, 0.1 mM EDTA, 0.1% NP40, 50% glycerol.

10X M-MLV RT Buffer

500 mM Tris-HCl (pH 8.3), 30 mM MgCl₂, 100 mM DTT, 750 mM KCl

Unit definition

One unit is the amount of enzyme required to incorporate 1 nmol of dTTP into acid-insoluble materials using 0.4 mM poly(rA)-oligo(dT) as substrate at 37°C in 10 min.

Storage buffer

50 mM Tris-HCl (pH 7.6), 150 mM NaCl, 1 mM DTT, 0.1 mM EDTA, 0.1% NP40, 50% glycerol.

10X M-MLV RT Buffer

500 mM Tris-HCl (pH 8.3), 30 mM MgCl₂, 100 mM DTT, 750 mM KCl