

TOPscript™ Reverse Transcriptase

Cat.#	Size	Conc.
RT002S	10,000 units	200 units/μl
RT002M	20,000 units	200 units/μl
RT002L	50,000 units	200 units/μl
RT002H	50,000 units	1000 units/μl

Store at -20°C

Supplied with: 10X TOPscript™ RT Buffer
dNTP Mixture (2 mM each)
Sterile water (RNase free)

India Contact:

Life Technologies (India) Pvt. Ltd.
306, Aggarwal City Mall, Opposite M2K Pitampura,
Delhi – 110034 (INDIA).
Ph: +91-11-42208000, 42208111, 42208222
Mobile: +91-9810521400
Fax: +91-11-42208444
Email: customerservice@lifetechindia.com
Web: www.lifetechindia.com

TOPscript™ Reverse Transcriptase

Cat.#	Size	Conc.
RT002S	10,000 units	200 units/μl
RT002M	20,000 units	200 units/μl
RT002L	50,000 units	200 units/μl
RT002H	50,000 units	1000 units/μl

Store at -20°C

Supplied with: 10X TOPscript™ RT Buffer
dNTP Mixture (2 mM each)
Sterile water (RNase free)

India Contact:

Life Technologies (India) Pvt. Ltd.
306, Aggarwal City Mall, Opposite M2K Pitampura,
Delhi – 110034 (INDIA).
Ph: +91-11-42208000, 42208111, 42208222
Mobile: +91-9810521400
Fax: +91-11-42208444
Email: customerservice@lifetechindia.com
Web: www.lifetechindia.com

Product description

TOPscript™ Reverse Transcriptase is genetically engineered version of M-MLV RT which is highly thermostable, thus can synthesize cDNA at elevated temperatures up to 60°C. This property is very useful when RNA templates are long and have extensive secondary structures. TOPscript™ Reverse Transcriptase is capable of synthesizing cDNA longer than 20 kb from messenger RNA.

Characteristics

- Molecular weight: 71 kDa
- Broad reaction temperature: 37°C~60°C
- Synthesis of long cDNA
- Excellent sensitivity

Applications

- Synthesis of first-strand cDNA,
- Array labeling
- cDNA library construction
- 3' and 5' RACE, RT-PCR
- Primer extension

Quality control

- Purity: >99% on SDS-PAGE
- Endonuclease-free
- Exonuclease-free
- RNase-free
- Inhibitor-free
- Satisfactory yield and length of cDNA products

For Research Use Only. Not for use in diagnostic procedures.

ISO9001 ISO14001 ISO13485

Product description

TOPscript™ Reverse Transcriptase is genetically engineered version of M-MLV RT which is highly thermostable, thus can synthesize cDNA at elevated temperatures up to 60°C. This property is very useful when RNA templates are long and have extensive secondary structures. TOPscript™ Reverse Transcriptase is capable of synthesizing cDNA longer than 20 kb from messenger RNA.

Characteristics

- Molecular weight: 71 kDa
- Broad reaction temperature: 37°C~60°C
- Synthesis of long cDNA
- Excellent sensitivity

Applications

- Synthesis of first-strand cDNA,
- Array labeling
- cDNA library construction
- 3' and 5' RACE, RT-PCR
- Primer extension

Quality control

- Purity: >99% on SDS-PAGE
- Endonuclease-free
- Exonuclease-free
- RNase-free
- Inhibitor-free
- Satisfactory yield and length of cDNA products

Unit definition

One unit is the amount of enzyme required to incorporate 1 nmol of dTTP into acid-insoluble materials using 0.4 mM poly(rA)-oligo(dT) as substrate at 37°C in 10 min.

Storage buffer

20 mM Tris-Cl (pH 7.5), 100 mM NaCl, 1 mM DTT, 0.1 mM EDTA, 0.01% NP40, 50% glycerol

10X TOPscript™ RT Buffer

500 mM Tris-HCl (pH 8.3), 30 mM MgCl₂, 100 mM DTT, 750 mM KCl

Note

TOPscript™ Reverse Transcriptase performs optionally over the full range of 42°C-60°C. Typically, 50°C is a good starting point. For RNAs containing secondary structure and other challenging targets, a synthesis temperature of 60°C may be used without loss of performance.

Unit definition

One unit is the amount of enzyme required to incorporate 1 nmol of dTTP into acid-insoluble materials using 0.4 mM poly(rA)-oligo(dT) as substrate at 37°C in 10 min.

Storage buffer

20 mM Tris-Cl (pH 7.5), 100 mM NaCl, 1 mM DTT, 0.1 mM EDTA, 0.01% NP40, 50% glycerol

10X TOPscript™ RT Buffer

500 mM Tris-HCl (pH 8.3), 30 mM MgCl₂, 100 mM DTT, 750 mM KCl

Note

TOPscript™ Reverse Transcriptase performs optionally over the full range of 42°C-60°C. Typically, 50°C is a good starting point. For RNAs containing secondary structure and other challenging targets, a synthesis temperature of 60°C may be used without loss of performance.

Standard reaction conditions

10X TOPscript™ RT Buffer	2 μl
TOPscript™ Reverse Transcriptase (200 units/μl)	1 μl
dNTP Mixture (2 mM each)	1 μl
a) Template RNA	X μl
b) Primer	1 μl
RNase Inhibitor (40 units/μl)	0.5 μl
Distilled water	up to 20 μl

a) Prepare one of the following RNA template.

- Total RNA: 1 ng~5 μg
- Messenger RNA (mRNA): 1 ng~250 ng
- Specific RNA: 0.01 pg~0.5 μg

b) Prepare one of the following primers.

- Oligo (dT)₁₈: 50~100 μM
- Random hexamer: 50~100 μM
- Specific primer: 15~20 pmol

→ An additional annealing step is recommended:

- if using oligo(dT)₁₈, incubate at 42°C for 5 min.
- if using random hexamer, incubate at 25°C for 10 min.

→ Incubate at 42°C-60°C for 60 min.

→ Incubate at 95°C for 5 min to inactivate the reaction.

Standard reaction conditions

10X TOPscript™ RT Buffer	2 μl
TOPscript™ Reverse Transcriptase (200 units/μl)	1 μl
dNTP Mixture (2 mM each)	1 μl
a) Template RNA	X μl
b) Primer	1 μl
RNase Inhibitor (40 units/μl)	0.5 μl
Distilled water	up to 20 μl

a) Prepare one of the following RNA template.

- Total RNA: 1 ng~5 μg
- Messenger RNA (mRNA): 1 ng~250 ng
- Specific RNA: 0.01 pg~0.5 μg

b) Prepare one of the following primers.

- Oligo (dT)₁₈: 50~100 μM
- Random hexamer: 50~100 μM
- Specific primer: 15~20 pmol

→ An additional annealing step is recommended:

- if using oligo(dT)₁₈, incubate at 42°C for 5 min.
- if using random hexamer, incubate at 25°C for 10 min.

→ Incubate at 42°C-60°C for 60 min.

→ Incubate at 95°C for 5 min to inactivate the reaction.