

TOPscript™ II Reverse Transcriptase (RNase H-)

Cat.#	Size	Conc.
RT003S	10,000 units	200 units/μl
RT003M	20,000 units	200 units/μl
RT003L	50,000 units	200 units/μl
RT003H	50,000 units	1000 units/μl

Store at -20°C

Supplied with: 10X TOPscript™ RT Buffer
dNTP Mixture (2 mM each)
Sterile water (RNase free)

India Contact:

Life Technologies (India) Pvt. Ltd.
306, Aggarwal City Mall, Opposite M2K Pitampura,
Delhi – 110034 (INDIA).
Ph: +91-11-42208000, 42208111, 42208222
Mobile: +91-9810521400
Fax: +91-11-42208444
Email: customerservice@lifetechindia.com
Web: www.lifetechindia.com

TOPscript™ II Reverse Transcriptase (RNase H-)

Cat.#	Size	Conc.
RT003S	10,000 units	200 units/μl
RT003M	20,000 units	200 units/μl
RT003L	50,000 units	200 units/μl
RT003H	50,000 units	1000 units/μl

Store at -20°C

Supplied with: 10X TOPscript™ RT Buffer
dNTP Mixture (2 mM each)
Sterile water (RNase free)

India Contact:

Life Technologies (India) Pvt. Ltd.
306, Aggarwal City Mall, Opposite M2K Pitampura,
Delhi – 110034 (INDIA).
Ph: +91-11-42208000, 42208111, 42208222
Mobile: +91-9810521400
Fax: +91-11-42208444
Email: customerservice@lifetechindia.com
Web: www.lifetechindia.com

Product description

TOPscript™ II Reverse Transcriptase (*RNase H-*) is genetically engineered recombinant protein that possesses moderate thermostability (42~55°C) and lacks RNase H activity. TOPscript™ II Reverse Transcriptase (*RNase H-*) does not degrade RNA strands from RNA-DNA hybrids and therefore is capable of high yield of synthesizing cDNA. The optimum temperature of TOPscript™ II Reverse Transcriptase (*RNase H-*) is 50°C. TOPscript™ II Reverse Transcriptase (*RNase H-*) does not show inhibitory effect on one-step RT-PCR reaction.

Characteristics

- Themostable and RNase H variant of M-MLV RTase (*RTase H-*)
- Synthesize cDNA at 42~55 °C
- Especially the highest activity at 50 °C
- Synthesis of long length cDNA (< 14 kb)
- Efficiently synthesize cDNA by inhibiting the formation of secondary structure of RNA

Applications

- First-strand cDNA synthesis
- RT-PCR
- Real-Time PCR
- Primer extension analysis

For Research Use Only. Not for use in diagnostic procedures.

ISO9001 ISO14001 ISO13485

Product description

TOPscript™ II Reverse Transcriptase (*RNase H-*) is genetically engineered recombinant protein that possesses moderate thermostability (42~55°C) and lacks RNase H activity. TOPscript™ II Reverse Transcriptase (*RNase H-*) does not degrade RNA strands from RNA-DNA hybrids and therefore is capable of high yield of synthesizing cDNA. The optimum temperature of TOPscript™ II Reverse Transcriptase (*RNase H-*) is 50°C. TOPscript™ II Reverse Transcriptase (*RNase H-*) does not show inhibitory effect on one-step RT-PCR reaction.

Characteristics

- Themostable and RNase H variant of M-MLV RTase (*RTase H-*)
- Synthesize cDNA at 42~55 °C
- Especially the highest activity at 50 °C
- Synthesis of long length cDNA (< 14 kb)
- Efficiently synthesize cDNA by inhibiting the formation of secondary structure of RNA

Applications

- First-strand cDNA synthesis
- RT-PCR
- Real-Time PCR
- Primer extension analysis

Quality control

- Purity: >99% on SDS-PAGE
- Endonuclease-free
- Exonuclease-free
- RNase-free
- Inhibitor-free
- Yield and length of cDNA product

Unit definition

One unit is the amount of enzyme required to incorporate 1 nmol dTTP into acid-insoluble materials using 0.4 mM poly(rA)-oligo(dT) as substrate at 37°C in 10 min.

Storage buffer

20 mM Tris-Cl (pH 7.5), 100 mM NaCl, 1 mM DTT, 0.1 mM EDTA, 0.01% NP40, 50% glycerol

10X TOPscript™ RT Buffer

500 mM Tris-HCl (pH 8.3), 30 mM MgCl₂, 100 mM DTT, 750 mM KCl

Note

TOPscript™ Reverse Transcriptase performs optionally over the full range of 42°C-60°C. Typically, 50°C is a good starting point. For RNAs containing secondary structure and other challenging targets, a synthesis temperature of 60°C may be used without loss of performance.

Quality control

- Purity: >99% on SDS-PAGE
- Endonuclease-free
- Exonuclease-free
- RNase-free
- Inhibitor-free
- Yield and length of cDNA product

Unit definition

One unit is the amount of enzyme required to incorporate 1 nmol dTTP into acid-insoluble materials using 0.4 mM poly(rA)-oligo(dT) as substrate at 37°C in 10 min.

Storage buffer

20 mM Tris-Cl (pH 7.5), 100 mM NaCl, 1 mM DTT, 0.1 mM EDTA, 0.01% NP40, 50% glycerol

10X TOPscript™ RT Buffer

500 mM Tris-HCl (pH 8.3), 30 mM MgCl₂, 100 mM DTT, 750 mM KCl

Note

TOPscript™ Reverse Transcriptase performs optionally over the full range of 42°C-60°C. Typically, 50°C is a good starting point. For RNAs containing secondary structure and other challenging targets, a synthesis temperature of 60°C may be used without loss of performance.

Your Molecular & Cell Technology Partner

Standard reaction conditions

10X TOPscript™ RT Buffer	2 μl
TOPscript™ II Reverse Transcriptase (200 units/μl)	1 μl
dNTP mixture (2 mM each)	2 μl
^{a)} Template RNA	X μl
^{b)} Primer	1 μl
RNase Inhibitor (40 units/μl)	0.5 μl
Sterile water (RNase free)	up to 20 μl

^{a)}Prepare one of the following RNA template.

- Total RNA: 1 ng~5 μg
- Messenger RNA (mRNA): 1 ng~250 ng
- Specific RNA: 0.01 pg~0.5 μg

^{b)}Prepare one of the following primers.

- Oligo (dT)₁₈: 50~100 μM
- Random hexamer: 50~100 μM
- Specific primer: 15~20 pmol

→An additional annealing step is recommended:

- if using oligo(dT)₁₈, incubate at 42°C for 5 min.
- if using random hexamer, incubate at 25°C for 10 min.

→Incubate at 42°C-60°C for 60 min.

→Incubate at 95°C for 5 min to inactivate the reaction.

Your Molecular & Cell Technology Partner

Standard reaction conditions

10X TOPscript™ RT Buffer	2 μl
TOPscript™ II Reverse Transcriptase (200 units/μl)	1 μl
dNTP mixture (2 mM each)	2 μl
^{a)} Template RNA	X μl
^{b)} Primer	1 μl
RNase Inhibitor (40 units/μl)	0.5 μl
Sterile water (RNase free)	up to 20 μl

^{a)}Prepare one of the following RNA template.

- Total RNA: 1 ng~5 μg
- Messenger RNA (mRNA): 1 ng~250 ng
- Specific RNA: 0.01 pg~0.5 μg

^{b)}Prepare one of the following primers.

- Oligo (dT)₁₈: 50~100 μM
- Random hexamer: 50~100 μM
- Specific primer: 15~20 pmol

→An additional annealing step is recommended:

- if using oligo(dT)₁₈, incubate at 42°C for 5 min.
- if using random hexamer, incubate at 25°C for 10 min.

→Incubate at 42°C-60°C for 60 min.

→Incubate at 95°C for 5 min to inactivate the reaction.