

FIREScript® RT cDNA synthesis MIX

Convenient cDNA synthesis mix with fast and robust FIREScript® Reverse Transcriptase



Genetically modified MMLV-based reverse transcriptase with increased thermostability and improved performance at elevated temperatures

- high specificity and yield
- wide reaction temperature from 37°C to 60°C
- fast 15 min reaction time
- a convenient 3-vial kit with all reaction components
- RNase inhibitor and water included
- reaction set-up and shipment without ice

Ordering

Some applications of this product may require a license which is not provided by the purchase of this product.

For research use only.

Choose Product Size

- 100 rxn | 100 x 20 µl rxn
- 500 rxn | 500 x 20 µl rxn
- 20 rxn | 20 x 20 µl rxn **free**

sample

REQUEST FOR BULK SIZE

— Reagents Provided

Item	Pcs.	Vial size
FIREScript® Enzyme mix	1	100 rxn 150 µl
10x RT Reaction Premix with random primers	1	250 µl
Water, nuclease free	2	1.25 ml

Description


FIREScript® Reverse Transcriptase (RT) is a genetically engineered MMLV (Moloney Murine Leukemia Virus) based Reverse Transcriptase.

This is an RNA-directed DNA polymerase that can synthesize a complementary DNA strand from ssRNA or ssDNA and is active over a broad range of reaction temperatures from 37°C to 60°C.


FIREScript® RT is a robust enzyme for RNA detection and has enhanced stability at room temperature with no activity loss for up to 1 month. This RT contains a functional RNase H domain which can increase the sensitivity of RT-qPCR (quantitative reverse transcription PCR).

Applications

- First-strand cDNA synthesis
- RT-PCR
- RT-qPCR



Highly competitive enzyme



Exceptional stability

Properties

Sample type: RNA

Final product: cDNA (first strand)

Product format: Enzyme mix with all reagents included

Enzyme: FIREScript® Reverse Transcriptase

Source: Purified from an E.coli strain that carries an overproducing plasmid containing a FIREScript® Reverse Transcriptase gene.

Unit description: One unit is defined as the amount of enzyme that will incorporate 1 nmol of dTTP into acid-precipitable material in 10 minutes at 37°C using poly(rA)•oligo(dT) as template in a total reaction volume of 50 µl.

Storage and dilution buffer: 50% glycerol (v/v), 20 mM Tris-HCl pH 7.5 at 25°C, 100 mM KCl, 0.1 mM EDTA and stabilizers.

Reagents

FIREScript® Enzyme Mix (FIREScript® RT and RiboGrip RNase inhibitor)

10x RT Reaction Premix with Random primers

RT reaction buffer with DTT, dNTPs and Random primers

Water, nuclease free

	Edge RT primer	Random primers	Edge RT and Random primers	Gene specific primers
Recommended Reaction	5 µl	5 µl	2.5 µl each	0.1 - 1 µl
Reaction volume at 42°C for 5-10 min	-	-	-	-
Benefits	Full length cDNA	All RNA in sample is converted to cDNA	Reducing 3' bias	All reaction resources used for genes of interest (increased sensitivity)
Disadvantages	Potential 3' bias; Can't bind to RNA lacking poly A sequence; Can bind to long poly A sequences in the middle of RNA	Truncated cDNA; multiple potential binding sites per RNA molecule	Edge RT doesn't bind to RNA lacking poly A sequence	Only 1 specific gene of interest per cDNA synthesis run can be analyzed downstream
Target RNA	RNA containing poly A tail (eg. mature messenger RNA, which accounts for 15% of total RNA in the cell)	Prokaryotic total RNA; All eukaryotic RNA types (mRNA, rRNA, tRNA); Degraded RNA	Continued targets of both primers	Specific gene of interest

cDNA synthesis priming options

India Contact:
Life Technologies (India) Pvt. Ltd.
 306, Aggarwal City Mall, Opposite M2K Pitampura,
 Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222
 Mobile: +91-9810521400
 Fax: +91-11-42208444
 Email: customerservice@lifetechindia.com
 Web: www.lifetechindia.com