

**ELISA kits available from ADI (see details at the web site)**

#0010	Human Leptin		
#200-120-AGH	Human globular Adiponectin (gAcrp30)		
#0700	Human Sex Hormone Binding Glob (SHBG)		
#0900	Human IGF-Binding Protein 1 (IGFBP1)		
#1000	Human C-Reactive Protein (CRP)		
#100-110-RSH	Human Resistin /FIZZ3		
#100-140-ADH	Human Adiponectin (Acrp30)		
#100-160-ANH	Human Angiogenin		
#100-180-APH	Human Angiopoietin-2 (Ang-2)		
#100-190-B7H	Human Bone Morphogenic Protein 7 (BMP-7)		
#1190	Human Serum Albumin	#1200	Human Albumin (Urinary)
#1750	Human IgG (total)	#1760	Human IgM
#1800	Human IgE	#1810	Human Ferritin
#1210	Human Transferrin (Tf)	#0020	Beta-2 microglobulin
#1600	Human Growth Hormone (GH)		
#0060	Human Pancreatic Colorectal cancer (CA-242)		
#1820	Human Ovarian Cancer (CA125)	#1830	Human CA153
#1840	Human Pancreatic & GI Cancer (CA199)		
#1310	Human Pancreatic Lipase		
#1400	Human Prostatic Acid Phosphatase (PAP)		
#1500	Human Prostate Specific Antigen (PSA)	#1510	free PSA (fPSA)
#0500	Human Alpha Fetoprotein (AFP)		
#0050	Human Neuron Specific Enolase (NSE)		
#0030	Human Insulin	#0040	Human C-peptide
#0100	Human Luteinizing Hormone (LH)		
#0200	Human Folicle Stimulating Hormone (FSH)		
#0300	Human Prolactin (PRL)		
#0400	Human Chorionic Gonadotropin (HCG)	#0410	HCG-free beta
#0600	Human Thyroid Stimulating Hormone (TSH)		
#1100	Human Total Thyroxine (T4)	#1110	Human Free T4 (fT4)
#1650	Human free triiodothyronine (fT3)	#1700	Human T3 (total)
#1850	Human Cortisol	#1860	Human Progesterone
#1865	Human Pregnenolone	#1875	Human Aldosterone
#1880	Human Testosterone	#1885	Human free Testosterone
#1910	Human Androstenedione	#1920	Human Estradiol
#1925	Human Estrone	#1940	Dihydrotestosterone (DHT)
#1950	Human DHEA-sulphate (DHEA-S)		
#3400	Human serum Neopterin		
#3000	Human Rheumatoid Factors IgM (RF)		
#3100	Human anti-dsDNA		
#3200	Anti-Nuclear Antibodies (ANA)		

*Instruction Manual No. M-1310*

## Human Pancreatic Lipase

**ELISA KIT Cat. No. 1310**

**For Quantitative Determination of Lipase In Serum**

*For In Vitro Research Use Only*



**ALPHA DIAGNOSTIC  
INTERNATIONAL**

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034  
Ph: +91-11-42208000, 42208111, 42208222

Mobile: +91-9810521400

Fax: +91-11-42208444

Email: [customerservice@atzlabs.com](mailto:customerservice@atzlabs.com)

Web: [www.atzlabs.com](http://www.atzlabs.com)



## Pancreatic Lipase ELISA Kit Cat. No. 1310

For Quantitative Determination of Pancreatic Lipase In Serum

Kit Contents: (reagents for 96 tests)

<b>C o m p o n e n t s</b>	<b>C a t . N o .</b>
Anti-h pancreatic Lipase coated microwell strip plate (96 wells). Ready-to-use	1 3 1 1
Pancreatic Lipase Sample Diluent 11 ml	S D - 1 3 1 0
Pancreatic Lipase Std. A, 0.75 ml (0 U/L)	1 3 1 2
Pancreatic Lipase Std. B, 0.75 ml (10 U/L)	1 3 1 3
Pancreatic Lipase Std. C, 0.75 ml (50 U/L)	1 3 1 4
Pancreatic Lipase Std. D, 0.75 ml (100 U/L)	1 3 1 5
Pancreatic Lipase Std. E, 0.75 ml (200 U/L)	1 3 1 6
Pancreatic Lipase Std. F, 0.75 ml (400 U/L)	1 3 1 7
Pancreatic Lipase <b>control serum</b> , 0.75 ml; #C1310 (exact values printed on vial)	
Anti-human Pancreatic Lipase HRP Conjugate, 11 ml (Ready-to-use)	1 3 1 8
Wash buffer (100X), 10 ml, <b>dilute 1:100 with distilled water</b>	W - 1 0 0
HRP substrate solution (ready-to-use), 11 ml	T M B 1 3 1 0
Stop solution (ready-to-use), 10 ml	T - 1 0
Complete Instruction Manual	1 3 1 0

### Introduction

Lipase (triacylglycerol acylhydrolase EC 3.1.1.3) hydrolyzes preferentially glycerol esters of long chain fatty acid, an enzyme that is anticipated to be specific for pancreas. The serum lipase activity tends to be elevated at about the same time as the elevation of serum amylase in acute pancreatitis. Analyses of pancreatic lipase and amylase in serum have added a new dimension to the laboratory detection and differentiation of pancreatic diseases. The current methodologies for determination of lipase activity using turbidimetric methods are relatively complex. These methods are non-specific and have low sensitivity.

ADI's human pancreatic lipase ELISA kit provides a direct assay in which pancreatic lipase is specifically recognized both by a solid phase antibody coated on the ELISA plates and the HRP-conjugated antibody. This eliminates the interference of the other isoamylase and the problem derived from non-standardized substrate preparation and conditions, which employed in most other enzymatic assay.

### 3. RECOVERY

A known amount of pancreatic Lipase (10-100 U/L) was added to three patient sera (with original Lipase concentrations of 22, 39, and 83 U/L) and the total Lipase measured. The assay showed excellent mean recoveries of about 105% (range 98-105%).

### 4. SPECIFICITY

The test is specific for human pancreatic lipase. The interference of transferin, gamma globulins, bilirubin, triglycerides, hemoglobin, and ascorbic acid were studied were not unusual and expected of any other ELISA test.

### 5. HIGH DOSE HOOK EFFECT

High pancreatic Lipase concentrations of up to 13200 U/L did not cause any hook effect.

### 6. Species Reactivity

Human pancreatic lipase test is specific for human samples. We have not tested other species (mouse, rat, monkey etc).

### References

1. Lott JA et al (1986) Clin Chem. 32, 1290;  
Lott JA et al (1986) Clin Enzymology, 245
2. Tietz NW (1986) Clin Chem. 32, 301-307
3. Panteghini M et al (1989) Clin Chem. 35, 417
4. Tietz NW et al (1987) Clin Chem. 33, 1624
5. Verduin PA (1973) Clin Chi Acta 46, 11-19

**TEST PROCEDURE** (ALLOW ALL REAGENTS TO REACH ROOM TEMPERATURE BEFORE USE). **Dilute wash buffer (1:100) with distilled water (10 ml stock in total of 1-liter).**

Label or mark the microtiter well strips to be used on the plate.

1. Pipet **25 ul of standards**, control, and serum samples into appropriate wells in *duplicate*. **Add 100 ul of Anti-lipase-enzyme conjugate** into each well. Mix gently for 5-10 seconds.
2. Cover the plate and incubate at room temp. for **60 minutes**.
3. Aspirate and wash the wells 5 times with 300 ul of 1x wash buffer. We recommend using an automated ELISA plate washer for better consistency. Failure to wash the wells properly will lead to high blank or zero values. If washing manually, plate must be tapped over paper towel between washings to ensure proper washing.
4. Dispense **100 ul TMB substrate per well**. Cover the plate and incubate at room temp. for **30 minutes**. Blue color develops in standards and positive wells.
5. Stop the reaction by adding **50 ul of stop solution** to all wells. Mix gently for 5-10 seconds (blue color turns yellow).
6. Measure the **Absorbance at 450 nm** using an ELISA reader. Color is stable for at least one hr after stopping.

## NOTES

Read instructions carefully before the assay. Do not allow reagents to dry on the wells. Careful aspiration of the washing solution is essential for good assay precision.

Since timing of the incubation steps is important to the performance of the assay, pipet the samples without interruption and it should not exceed 5 minutes to avoid assay drift. If more than one plate is being used in one run, it is recommended to include a standard curve on each plate. The unused strips should be stored in a sealed bag at 4°C.

Addition of the HRP substrate solution starts a kinetic reaction, which is terminated by dispensing the stopping solution. Therefore, keep the incubation time for each wells the same by adding the reagents in identical sequence. Plate readers measure Absorbance vertically. Do not touch the bottom of the wells.

## CALCULATION OF RESULTS

Calculate the mean Absorbance for each duplicate. Subtract the Absorbance of the zero standard from the mean Absorbance values of standards, control, and samples. Draw the standard curve on log-log graph paper by plotting net Absorbance values of standards against appropriate pancreatic Lipase concentrations. Read off the pancreatic Lipase concentrations of the control and patient samples.

### Expected Values

It is recommended that each laboratory determine its own normal and abnormal range.

A clinical study of ADI ELISA kit on 36 apparently normal sample yielded values of 26-150 U/L.

**Reference range** : 56-239 U/L.

## PERFORMANCE CHARACTERISTICS

### 1. DETECTION LIMIT

Based on sixteen replicate determinations of the zero standard, the minimum pancreatic Lipase concentration detectable using this assay is 1.0 U/L. The detection limit is defined as the value deviating by 2 SD from the zero standard.

### 2. PRECISION

#### *Intra-assay precision:*

Three serum samples (mean Lipase concentrations: 44-100 U/L were run in sixteen replicates. The samples showed good intra-assay precision with %CV of 4.9 and S.D. 2.6 U/L.


#### *Inter-assay precision:*

Twelve serum samples were run in duplicate in sixteen independent assays. The samples showed good inter-assay precision (7-9 % CV). The actual values were: mean 42.33 U/L , 71.92 U/L, and 137 U/L

## WORKSHEET OF TYPICAL ASSAY

Wells	Stds/samples	Net Abs.	Calculated Conc'n (U/L)
A1, A2	Std. A (0 U/L)	-	
B1, B2	Std. B (10 U/L)	0.071	
C1, C2	Std. C (50 U/L)	0.302	
D1, D2	Std. D (100 U/L)	0.580	
E1, E2	Std. E (200 U/L)	1.156	
F1, F2	Std. F (400 U/L)	1.920	22
G1, G2	Sample 1	0.437	75

Note: These data are for demonstration purpose only. A complete standard curve must be run in every assay to determine sample values. Each laboratory should determine their own normal reference values.



A typical std. assay curve (do not use this for calculating sample values)

## PRINCIPLE OF THE TEST

Pancreatic Lipase ELISA kit is based on sequential binding of human pancreatic Lipase from samples to two antibodies, one immobilized on microtiter well plates, and other conjugated to the enzyme horseradish peroxidase. After a washing step, chromogenic substrate is added and colors developed. The enzymatic reaction (color) is directly proportional to the amount of Lipase present in the sample. Adding stopping solution terminates the reaction. Absorbance is measured on a microtiter well ELISA reader at 450 nm. The unknown sample values are then read-off the standard curve.

## MATERIALS AND EQUIPMENT REQUIRED

Adjustable micropipet (20-100  $\mu$ l) and multichannel pipet with disposable plastic tips. Reagent troughs, plate washer (recommended) and ELISA plate Reader.

## PRECAUTIONS

The ADI Pancreatic Lipase ELISA kit is intended for *in vitro research* use only. The reagents contain thimerosal as preservative; necessary care should be taken when disposing solutions. The Control Serum and Standards have been prepared from human sera shown to be negative for HBsAg and HIV antibodies. Nevertheless, such tests are unable to prove the complete absence of viruses, therefore, sera should be handled with appropriate precautions.

Applicable **MSDS**, if not already on file, for the following reagents can be obtained from ADI or the web site.

TMB (substrate), 2N HCl (stop solution), and Proclin-300 (0.1% v/v in standards, sample diluent and HRP-conjugates).

## SPECIMEN COLLECTION AND HANDLING

Collect blood by venipuncture, allow to clot, and separate the serum by centrifugation at room temperature. Do not heat inactivate the serum. If sera can not be immediately assayed, these should be stored at  $-20^{\circ}\text{C}$  for up to six months. Avoid repeated freezing and thawing of samples. No preservatives should be added to the serum.

## STORAGE AND STABILITY

The microtiter well plate and all other reagents are stable at  $2-8^{\circ}\text{C}$  until the expiration date printed on the label. The whole kit stability is usually 6 months from the date of shipping under appropriate storage conditions.

## Reagent Preparation:

Dilute wash buffer (1:100) with distilled water (10 ml stock in total of 1-liter). Store at  $4^{\circ}\text{C}$ .