

Phytohaemagglutinin (PHA)- Liquid form

Principle

Peripheral blood is the most commonly utilized tissue for determination of constitutional karyotype. Peripheral blood lymphocytes are not dividing in normal, healthy adult; they must be stimulated to divide by exposure to mitogen.

The most commonly used mitogen for Cytogenetic studies is PHA. It is a mucoprotein isolated from seeds of *Phaseolus vulgaris*. It has the property of selectively agglutinating and sedimenting erythrocytes and of inducing division in leukocytes. When the T-lymphocytes are exposed to this plant antigen (PHA), initiates immunological response through behaviors resembling an antigen. The mature T cells dedifferentiate to a T-lymphoblastic cell, which has the capacity to synthesize DNA to under go mitosis. During the first 24 hrs exposure to PHA, the T cells undergo this transformation and as the culture duration increases undergo several rounds of cell division.

Storage Conditions: Store frozen at -20°C

After thawing PHA is stable for at least one month at +2/+8°C.

Shelf Life : 24 months

Pack Size: 5ml

Recommended Use:

- † Store the product at the mentioned temperature condition.
- † Do not use the product after expiry.
- † Use proper PPE at the time of handling product.
- † The storage and work environmental should be free from humidity.
- † The product is intended to used in laboratory only. Do not use in therapy, human or veterinary application.

Final Concentration

† 2.5/3ml per 100ml of culture medium
(Add 2.5/3ml of PHA in 100ml of peripheral blood culture medium for Karyotype)

OR

† 0.1/0.2ml per 5ml of culture medium
(Add 0.1/0.2ml of PHA in 5ml of peripheral blood culture medium for Karyotype)

Methodology-

Working Medium Preparation

- 1) Add 200mg of sodium bicarbonate and 2.5/3ml of PHA in 100ml of peripheral blood plain culture medium (Preferably RPMI 1640).
- 2) Adjust the PH at 7.2 with 1N HCL/NAOH.

Cell culture

- 3) Label the 15ml centrifuge tube with patient name and lab ID
(Do not use culture flask, as it is close system culture)
- 4) Add 5ml of working medium and 1ml of FBS
- 5) Add 0.3/0.5ml of heparinised whole blood sample
- 6) Incubate the cell culture at 37°C for 72 hrs (No need of CO2 incubator)

Chromosome preparation

- 7) After completion of 72hrs add 0.2ml colcemid solution
- 8) Incubate the culture for another 30 minutes
- 9) Centrifuge the tube at 1000rpm for 10 minutes
- 10) Remove the supernatant and add 7-8ml hypotonic solution (0.075M KCL)
- 11) Incubate the tube at 37°C for 15mints
- 12) Add 5-6 drops of chilled fixative (3:1-Methanol: Glacial acetic acid), mix slowly
- 13) Centrifuge at 1000rpm for 10mints
- 14) Remove supernatant and add 7-8ml chilled fixative, mix it properly.
- 15) Repeat steps 13 and 14 for 3 times
- 16) In the final wash, will get white lymphocyte clear pellet

Slide preparation

- 17) Re-suspend the cell pellet in a small volume of chilled fixative approx. 0.4 to 0.8ml
- 18) Drop the cells on chilled glass slide and allow to dry
- 19) Metaphases can be observe under phase contrast microscope

Trouble Shooting

- † PHA may appear cloudy at +2/+8°C, but this turbidity has no effects on the activity of the products.
- † Do not add more then 4ml of PHA per 100ml of medium
- † Always run parallel culture before use of every lot

Ref

1. D.E. Rooney, Human Cytogenetics constitutional analysis, Third edition, PP. 33-44.
2. Marilyn S. Arsham, Margaret J. Barch, The AGT cytogenetics Laboratory manual, Fourth edition, PP 91-100.
3. Verma RS, Babu A., Human Chromosome- Manual of Basic Techniques, 1988, PP.9-10.

India Contact:

Life Technologies (India) Pvt. Ltd.

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com Mobile: +91-9810521400