


Cas9 Nuclease

Manual



Important!

-20°C Storage Required

- * Immediately inspect packages
- * Freeze upon receipt

FOR RESEARCH USE ONLY. NOT FOR HUMAN OR DIAGNOSTIC USE

Catalog #3273 and #3276

India Contact:
Life Technologies (India) Pvt. Ltd.
Ph: +91-11-42208000, 42208111, 42208222
Mobile: +91-9810521400
Email: customerservice@lifetechindia.com
www.lifetechindia.com




Description:

CRISPR-associated (Cas) systems


Clustered regularly interspaced short palindromic repeats (CRISPR)/CRISPR-associated (Cas) systems provide bacteria and archaea with adaptive immunity against viruses and plasmids by using CRISPR RNAs (crRNAs) to guide the silencing of invading nucleic acids (1). The CRISPR system consists of a short non-coding guide RNA (sgRNA) made up of a target complementary CRISPR RNA (crRNA) and an auxiliary transactivating crRNA (tracrRNA). The sgRNA guides the Cas9 endonuclease to a specific genomic locus via base pairing between the crRNA sequence and the target sequence, and cleaves the DNA to create a double-strand break. The location of the break is within the target sequence 3 bases from the NGG PAM (Protospacer Adjacent Motif) (1). The PAM sequence, NGG, must follow the targeted region on the opposite strand of the DNA with respect to the region complementary sgRNA sequence (Fig.1).

Fig. 1: Overview of the CRISPR-associated (Cas) systems.



Intact Genomics (ig®)Cas9 Nuclease is the purified recombinant *Streptococcus pyogenes* Cas9 enzyme containing a nuclear localization signal (NLS) at the C-terminal for targeting to the nucleus. This enzyme is designed to perform CRISPR/Cas9-mediated genome editing (1, 2). The physical purity of this enzyme is $\geq 98\%$ as assessed by SDS-PAGE with Coomassie® blue staining.

Fig. 2: Lane 1. Protein Marker
Lane 2. Cas9 Nuclease



Product Source

E. coli BL21 (DE3) strain expressing a Cas9 gene from *Streptococcus pyogenes* with an N-terminal 6xHis tag and C-terminal SV40 nuclear localization signal (NLS).

Components and Storage

Cas9 Nuclease Kits contains the below items. **Store all components at -20°C.**

- Cas9 Nuclease
- 10x Cas9 Nuclease Reaction Buffer
- **Storage Buffer**
 - 50 mM Tris-HCl, 50 mM KCl, 1 mM DTT, 0.1 mM EDTA, 50% Glycerol, pH 7.5 @ 25 °C
- **1x Cas9 Reaction Buffer**
 - 20 mM HEPES, 100 mM NaCl, 5 mM MgCl₂, 0.1 mM EDTA, pH 6.5 @ 25 °C

Quality Control Assays

Cas9 nuclease is free from detectable RNase, Endonuclease (nicking) and non-specific DNase activities.


Functional Testing

Cas9 Nuclease functional testing was done by in vitro DNA cleavage assay with the following protocol which gives more than 95% digestion of the substrate DNA as determined by agarose gel electrophoresis (Fig. 3).

1. Set up 30 µl reaction in a microcentrifuge tube on ice with the following combinations.

Target DNA	x µl (~100 ng)
sgRNA	x µl (~4000 ng)
10x Cas9 Reaction Buffer	3.0 µl
Cas9 Nuclease	1.0 µl (~160 ng)
Add H ₂ O up to	30.0 µl

2. Gently mix the reaction mixture and centrifuge briefly.
3. Incubate at 37 °C for 60 min.
4. Add 1 µl RNase (4 mg/ml)
5. Incubate at 37 °C for 20 min.
6. Run 0.7 to 1% agarose TBE gel



Related Products

- Taq DNA Polymerase(Cat.# 3243)
- Taq DNA Polymerase 2x Premix(Cat.# 3249)
- T4 DNA Ligase(Cat.# 3212)
- ig[®] 10B Chemically Competent Cells(Cat.# 1011-12)

References

Jinek M, Chylinski K, Fonfara I, Hauer M, Doudna JA, Charpentier E. (2012) A programmable dual-RNA-guided DNA endonuclease in adaptive bacterial immunity. *Science*. Aug 17;337(6096):816-21.

Mali P, Yang L, Esvelt KM, Aach J, Guell M, DiCarlo JE, Norville JE, Church GM. (2013) RNA-guided human genome engineering via Cas9. *Science*. Feb 15;339(6121):823-6.

Technical Support

Intact Genomics (IG[®]) is dedicated to customer satisfaction regarding the use of our products for your research needs. We test our products thoroughly to ensure they conform to the highest quality standards and provide excellent results when following the protocol's specifications. Please follow the protocol information provided in this manual carefully and contact our customer and technical support team with any questions or comments you may have regarding this or our other products.

IG Technical Support

Email: sales@intactgenomics.com

Phone: (314) 942-3655 | Toll free: 855-835-7172

Intact Genomics is committed to supporting the worldwide scientific research community by supplying the highest quality life science products to help research scientists worldwide discover solutions to critical challenges in human health, agriculture and the environment. Each new lot of our products is tested to ensure they meet the quality standards and specifications designated for the product. Please follow the instructions carefully and contact us if additional assistance is needed.

We appreciate your business and your feedback regarding the performance of our products in your applications. Our hours are Monday - Friday, 8AM to 5PM, U.S. central standard time.

<http://www.intactgenomics.com>

Legal Notices

Intact Genomics owns the following registered trademarks granted by the United States Patent and Trademark Office (USPTO): Intact Genomics®, IG®, ig®, igTherapeutics®, FastAmp®, i7®, DirectPlate® .

All technology protocols discussed within this manual are assumed proprietary to IG's products, and if registered, are protected by US patent and trademark laws.

Notice to Buyer - Limited License

Buyer agrees by purchasing the product that it's intended use is for research purposes only, not for human or diagnostic use. Purchase of this material conveys to buyer the non-transferable right to use the material purchased in research conducted by buyer, whether for teaching, non-commercial or commercial research purposes. Buyer may not sell or otherwise transfer these materials, its components, or unmodified descendants to a third party. Use of these products maybe covered by IG's patents or trademarks.

Products are understood by buyer to be experimental in nature and may have hazardous properties. Unless prohibited by law, buyer assumes all liability for claims or damages against it by third parties that relate to or arise from the use, storage, or disposal of the purchased materials.

Buyer agrees to use the purchased materials in full compliance with applicable law and regulations. Intact Genomics will not be held liable for activities outside of itself and the products intended use.

For more information or to express legal concerns relating to patent and trademarks, please contact us at: info@intactgenomics.com or sales@intactgenomics.com

Thank you for your business!

India Contact:
Life Technologies (India) Pvt. Ltd.
Ph: +91-11-42208000, 42208111, 42208222
Mobile: +91-9810521400
Email: customerservice@lifetechindia.com
www.lifetechindia.com