

Glyphosate

• Intended Use

For the detection and quantitation of glyphosate in water (groundwater, surface water, well water). For soil, crop, and food use contact the company for application bulletins and/or specific matrix validation guidelines.

• Principle

The Abraxis Glyphosate Kit applies the principles of enzyme linked immunosorbent assay (ELISA) to the determination of glyphosate. The sample to be tested is derivatized and then added, along with paramagnetic particles attached with antibodies specific to glyphosate and incubated for 30 minutes. The glyphosate enzyme conjugate is then added, at this point a competitive reaction occurs between the glyphosate which may be in the sample and the enzyme labeled glyphosate analog for the antibody binding sites on the magnetic particles. The reaction is allowed to continue for thirty (30) minutes. At the end of the incubation period, a magnetic field is applied to hold in the test tube the para-magnetic particles (with glyphosate and labeled glyphosate bound to the antibodies on the particles, in proportion to their original concentration), and allow the unbound reagents to be decanted. After decanting, the particles are washed with Washing Solution.

The presence of glyphosate is detected by adding the "Color Solution", which contains the enzyme substrate (hydrogen peroxide) and the chromogen (3,3',5,5'-tetramethylbenzidine). The enzyme-labeled glyphosate bound to the glyphosate antibody catalyzes the conversion of the substrate/ chromogen mixture to a colored product. After an incubation period, the reaction is stopped and stabilized by the addition of a diluted acid (Stopping Solution). Since the labeled glyphosate (conjugate) was in competition with the unlabeled glyphosate (sample) for the antibody sites, **the color developed is inversely proportional to the concentration of glyphosate in the sample.**

• Reagents

The Abraxis Glyphosate Kit contains the following items:

1. Glyphosate Antibody Coupled Paramagnetic Particles

Glyphosate antibody (rabbit anti-glyphosate) covalently bound to paramagnetic particles suspended in a buffered solution with preservative and stabilizers.

120 test kit: one 65 mL vial

2. Glyphosate Enzyme Conjugate

Horseshoe peroxidase (HRP) labeled glyphosate analog diluted in a buffered solution with preservative and stabilizers.

120 test kit: one 35 mL vial

3. Glyphosate Standards

Four concentrations (75, 200, 750, and 4000 parts per trillion) of glyphosate standards in distilled water with preservative and stabilizers. Each vial contains 2.0 mL.

4. Control

A concentration (approximately 500 ppt of glyphosate in distilled water with preservative and stabilizers. A 2.0 mL volume is supplied in one vial.

5. Diluent/Zero Standard

Distilled water with preservative and stabilizers without any detectable glyphosate.

120 test kit: one 65 mL vial

6. Color Solution

A solution of hydrogen peroxide and 3,3',5,5'-tetramethylbenzidine in an organic base.

120 test kit: one 65 mL vial

7. Stopping Solution

A solution of diluted acid.

120 test kit: one 65 mL vial

8. Washing Solution

Preserved deionized water.

120 test kit: one 250 mL vial

9. Assay Buffer

Dissolved buffer salts.

120 test kit: one 125 mL vial

10. Derivatization Reagent

120 test kit: three 80 uL vials

11. Derivatization Reagent Diluent

Dimethyl Sulfoxide (DMSO): three 4 mL vials

12. Test Tubes

Glass (36) are packaged in boxes.

120 test kit: 4 X 36 tube boxes

• Reagent Storage and Stability

Store all reagents at 2-8°C. Do not freeze.

Reagents may be used until the expiration date on the box, except for derivatization reagent (use the same day as diluted). *The test tubes and Washing Solution require no special storage condition and may be stored separately from the reagents to conserve refrigerator space.*

The Derivatization Reagent Diluent may freeze if stored cool, thaw reagent by placing on a 37 C incubator.

Consult state, local and federal regulations for proper disposal of all reagents.

• Materials Required but Not Provided

In addition to the reagents provided, the following items are essential for the performance of the test:

Pipets* Precision pipets capable of delivering 100, 250, 500, 750 uL and a 1.0 mL repeating pipet. Disposable 5 mL pipette.

Vortex Mixer* Thermolyne Maxi Mix, Scientific Industries Vortex Genie, or equivalent

Magnetic Separation System

Photometer* capable of readings at 450 nm

*Please contact Abraxis for supplier information.

• Sample Information

This procedure is recommended for use with water samples. Other samples may require modifications to the procedure and should be

Samples containing gross particulate matter should be filtered (e.g. 0.2 um Anotop™ 25 Plus, Whatman, Inc.) to remove particles.

Samples which have been preserved with monochloroacetic acid or other acids, should be neutralized with strong base e.g. 6N NaOH, prior to assay.

If the glyphosate concentration of a sample exceeds 4 ppb, the sample is subject to repeat testing using a diluted sample. A ten-fold or greater dilution of the sample is recommended with an appropriate amount of Diluent/Zero Standard or Sample Diluent. For example, in a separate test tube make a ten-fold dilution by adding 100 uL of the sample to 900 uL of Diluent/Zero Standard. Mix thoroughly before assaying. Perform the assay according to the Assay Procedure and obtain final results by multiplying the value obtain by the dilution factor e.g. 10.

The presence of the following substances up to 20,000 ppm were found to have no significant effect on the Glyphosate Assay results: calcium, magnesium, nitrate, sodium fluoride, copper, carbonate. Sulfate and potassium up to 2,000 ppm. Phosphate up to 100 ppm. Humic acid up to 20 ppm. Sodium chloride up to 1.0 M; and HCl up to 0.25 N.

Solvents usually used to extract pesticides from soil or plant matrices such as methanol and acetone were found to be acceptable for use in the Glyphosate immunoassay up to 50%.

• Reagent Preparation

All reagents must be allowed to come to room temperature. The antibody coupled paramagnetic particles should be mixed thoroughly before use.

Derivatization of Standards, Control, and Samples

1. Dilute Derivatization Reagent with 3.5 mL of Derivatization Reagent Diluent (Diluted Reagent needs to be used within the same day). Mixed thoroughly.

2. Label single test tubes for standards, control, and samples.

3. Pipette 250 uL of standard, control, sample(s) into separate disposable tubes.

4. Add 1.0 mL of Assay buffer, vortex to mix.

5. Add 100 uL of the diluted derivatization reagent, **vortex each tube immediately after addition of reagent.** We recommend vortexing until no swirl lines are seen in the tube.

6. Incubate at room temperature for 10 minutes.

7. Perform the ELISA as in Assay Procedure, start with step 1 of Assay procedure.

Alternative Derivatization Procedure

Performing the alternative derivatization procedure allows the user to use the same derivatization tubes in the performance of the assay, therefore eliminating the use of additional assay tubes.

1. Dilute Derivatization Reagent with 3.5 mL of Derivatization Reagent Diluent (Diluted Reagent needs to be used within the same day). Mixed thoroughly.

2. Label test tubes in duplicate for standards, control, and samples.

3. Pipette in duplicate, 50 uL of standard, control, sample(s) into disposable assay tubes.

4. Add 200 uL of Assay buffer, vortex to mix.

5. Add 20 uL of the diluted derivatization reagent, **vortex each tube immediately after addition of reagent.** We recommend vortexing until no swirl lines are seen in the tube.

6. Incubate at room temperature for 10 minutes.

7. Perform the ELISA as in Assay Procedure, starting with step 3 of Assay Procedure.

• Procedural Notes and Precautions

As with all immunoassays, a consistent technique is the key to optimal performance. To obtain the greatest precision, be sure to treat each tube in an identical manner.

Add reagents directly to the bottom of the tube while **avoiding contact between the reagents and the pipet tip.** This will help assure consistent quantities of reagent in the test mixture.

Avoid cross-contaminations and carryover of reagents by using clean pipets for each sample addition and by avoiding contact between reagent droplets on the tubes and pipet tips.

Avoid foam formation during vortexing.

The magnetic separation system consists of two parts: an upper rack which will securely hold the test tubes and a lower separator which contains the magnets used to attract the antibody coupled paramagnetic particles. During incubations the upper rack is removed from the lower separator so that the paramagnetic particles remain suspended during the incubation. **For separation steps, the rack and the separator**

are combined to pull the paramagnetic particles to the sides of the tubes.

To obtain optimum assay precision, it is important to perform the separation steps carefully and consistently. Decant the rack by slowly inverting away from the operator using a smooth turning action so the liquid flows consistently along only one side of the test tube. While still inverted, place the rack on an absorbent pad and allow to drain. Lifting the rack and replacing gently onto the pad several times will ensure complete removal of the liquid from the rim of the tube. Do not bang the rack.

Mix the antibody coupled paramagnetic particles just prior to pipetting.

Do not use any reagents beyond their stated shelf life.

Do not use the diluted derivatization reagent after 24 hours from dilution.

Avoid contact of Stopping Solution (diluted sulfuric acid) with skin and mucous membranes. If this reagent comes in contact with skin, wash with water.

• Limitations

The Abraxis Glyphosate Assay will detect glyphosate. Refer to specificity table for data on several of related compounds. The Abraxis Glyphosate Assay kit provides screening results. As with any analytical technique (GC, HPLC, etc...) positive results requiring some action should be confirmed by an alternative method.

The total time required for pipetting the magnetic particles should be kept to **two (2) minutes or less**, therefore the total number of tubes that can be assayed in a run should be adjusted accordingly.

• Quality Control

A control solution at approximately 500 ppt of Glyphosate is provided with the Abraxis Glyphosate Assay kit. It is recommended that it be included in every run and treated in the same manner as unknown samples. Acceptable limits should be established by each laboratory.

• Assay Procedure

Read Reagent Preparation, Procedural Notes and Precautions before proceeding.

1. Label test tubes for standards, control, and samples.

Tube Number	Contents of Tube
1,2	Diluent/Zero Standard, 0 ppt
3,4	Standard 1, 75 ppt
5,6	Standard 2, 200 ppt
7,8	Standard 3, 750 ppt
9,10	Standard 4, 4000 ppt
11,12	Control
13,14	Sample 1
15,16	Sample 2
17,18	Sample 3

2. Add 300 uL of the appropriate derivatized standard, control, or sample.
3. Mix the Glyphosate Antibody Coupled Paramagnetic Particles thoroughly and add 500 uL to each tube.
4. Vortex for 1 to 2 seconds minimizing foaming.
5. Incubate for 30 minutes at room temperature.
6. Add 250 uL of Glyphosate Enzyme Conjugate to each tube.
7. Vortex for 1 to 2 seconds minimizing foaming.
8. Incubate for 30 minutes at room temperature.
9. Separate in the Magnetic Separation System for **two (2) minutes**.
10. Decant and **gently** blot all tubes briefly in a consistent manner.
11. Add 1 mL of Washing Solution to each tube and allow them to remain in the magnetic separation unit for **two (2) minutes**.
12. Decant and **gently** blot all tubes briefly in a consistent manner.

13. Repeat Steps 11 and 12 two (2) additional times.
14. Remove the rack from the separator and add 500 uL of Color Solution to each tube.
15. Vortex for 1 to 2 seconds minimizing foaming.
16. Incubate for 20 minutes at room temperature.
17. Add 500 uL of Stopping Solution to each tube.
18. Add 1 mL Washing Solution to a clean test tube. Use as blank in Step 19.
19. Read results at 450 nm within 15 minutes after adding the Stopping Solution.

• Results

Manual Calculations

1. Calculate the mean absorbance value for each of the standards.
2. Calculate the %B/Bo for each standard by dividing the mean absorbance value for the standard by the mean absorbance value for the Diluent/Zero Standard.
3. Construct a standard curve by plotting the %B/Bo for each standard on vertical linear (Y) axis versus the corresponding glyphosate concentration on horizontal log (X) axis on the graph paper provided.
4. %B/Bo for controls and samples will then yield levels in ppb of glyphosate by interpolation using the standard curve.

(Contact Abraxis for detailed application information on specific photometers.)

Photometric Analyzer

Some instrument manufacturers make available photometers allowing for calibration curves to be automatically calculated and stored. Refer to instrument operating manual for detailed instructions. To obtain results for the Abraxis Glyphosate HS Assay on instruments allowing data transformation the following parameter settings are recommended:

Data Reduct : Lin. Regression
Xformation : Ln/Ln
Read Mode : Absorbance
Wavelength : 450 nm
Units : PPT
Rgt Blk : 0

Calibrators:
of Cals : 5
of Reps : 2

Concentrations:
#1: 0.00 PPT
#2: 75 PPT
#3: 200 PPT
#4: 750 PPT
#5: 4000 PPT

Range : 75 – 4000
Correlation : 0.990
Rep. %CV : 15%

NOTE: Any results obtained with a calculated glyphosate concentration of less than 50 ppt on the print out should be assumed to be below the detection limit of the assay .

• Expected Results

In a study with water samples from various locations, the Abraxis Glyphosate HS Assay was shown to correlate well with another analytical technique.

• Performance Data

Precision

The following results were obtained:

Control	1	2	3
Replicates	5	5	5
Days	5	5	5
n	25	25	25
Mean (ppb)	0.98	2.82	5.80
% CV (within assay)	6.0	3.5	6.9
% CV (between assay)	15.5	11.6	9.5

Sensitivity

The Abraxis Glyphosate HS Assay has an estimated minimum detectable concentration based on a 90% B/Bo of 50 parts per trillion (ppt).

Recovery

Five (5) groundwater samples, were spiked with various levels of glyphosate and then assayed using the Abraxis Glyphosate HS Assay. The following results were obtained:

Amount of Glyphosate Added (ppb)	Mean (ppb)	Recovery (ppb)	S.D. (ppb)	%
0.50	0.47	0.09	95	
1.0	1.04	0.13	104	
2.5	2.70	0.41	108	
Average				102

Specificity

The cross-reactivity of the Abraxis Glyphosate HS Assay for various related analogues can be expressed as the least detectable dose (LDD) which is estimated at 90% B/Bo, or as the dose required for 50% absorbance inhibition (50% B/Bo).

B/Bo Compound	LDD (ppb)	50% (ppb)
Glyphosate	0.05	2.40
Glyphosine	50	3,000
Glufosinate	2000	70,000
AMPA	35,000	>1,000,000
Glycine	>10,000	>1,000,000

The following compounds demonstrated no reactivity in the Abraxis Glyphosate Assay at concentrations up to 1000 ppb: aldicarb, aldicarb sulfoxide, aldicarb sulfone, acetochlor, alachlor, atrazine, ametryn, benomyl, butylate, captan, carbaryl, carbendazim, carbofuran, cyanazine, 2,4-D, 1,3-dichloropropene, dinoseb, MCPA, metolachlor, metribuzin, pentachlorophenol, picloram, propazine, simazine, terbufos, thiabendazole, and thiophanate-methyl.

• Ordering information

Abraxis Glyphosate Assay Kit,120T	PN 500081
Sample Diluent	PN 500082
Standard Set	PN 500083
High Sensitivity Reagent Set	PN 500084
High Sensitivity Standard Set	PN 500085

• Assistance

For ordering or technical assistance contact:

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Road No. 44, Pitampura, Delhi – 110034, India
Mobile: +91-98105-21400, Tel: +91-11-42208000, 8111, 8222, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com, www.atlabs.com ; www.lifetechindia.com

Phone: (215) 357-3911

Fax: (215) 357-5232

Email: info@abraxiskits.com












WEB: www.abraxiskits.com

• General Limited Warranty

Abraxis LLC warrants the products manufactured by the Company, against defects and workmanship when used in accordance with the applicable instructions for a period not to extend beyond the product's printed expiration date.

Abraxis makes no other warranty, expressed or implied. There is no warranty of merchantability or fitness for a particular purpose.

GLYPHOSATE DETAILED FLOWCHART

<p>1.</p>  <p>Remove upper rack from magnetic base. Label test tubes for Standards, Control, and Samples.</p> <table border="1"> <thead> <tr> <th>Tube #</th> <th>Content</th> </tr> </thead> <tbody> <tr> <td>1, 2</td> <td>Diluent/Zero Standard 0 ppb</td> </tr> <tr> <td>3, 4</td> <td>Standard 1, 0.15 ppb</td> </tr> <tr> <td>5, 6</td> <td>Standard 2, 1.0 ppb</td> </tr> <tr> <td>7, 8</td> <td>Standard 3, 5.0 ppb</td> </tr> <tr> <td>9,10</td> <td>Control</td> </tr> <tr> <td>11,12</td> <td>Sample 1</td> </tr> <tr> <td>13,14</td> <td>Sample 2</td> </tr> <tr> <td>15,16</td> <td>Sample 3</td> </tr> </tbody> </table> <p>Add 300 μL of either Derivatized Standards, Control or Samples to the bottom of each test tube by inserting the pipette tip all the way into the bottom of the tube without touching the sides of the tube.</p>	Tube #	Content	1, 2	Diluent/Zero Standard 0 ppb	3, 4	Standard 1, 0.15 ppb	5, 6	Standard 2, 1.0 ppb	7, 8	Standard 3, 5.0 ppb	9,10	Control	11,12	Sample 1	13,14	Sample 2	15,16	Sample 3	<p>7.</p>  <p>Do not separate upper rack from lower base. Using a smooth motion, invert the combined rack assembly over a sink and pour out the tube contents; keep inverted and gently blot the test tube rims on several layers of paper toweling.</p>
Tube #	Content																		
1, 2	Diluent/Zero Standard 0 ppb																		
3, 4	Standard 1, 0.15 ppb																		
5, 6	Standard 2, 1.0 ppb																		
7, 8	Standard 3, 5.0 ppb																		
9,10	Control																		
11,12	Sample 1																		
13,14	Sample 2																		
15,16	Sample 3																		
<p>2.</p>  <p>Add 500 μL of thoroughly mixed Glyphosate Antibody Coupled Magnetic Particles down the inside wall of each tube by using the technique described in Box 2. Vortex for 1 to 2 seconds (at low speed to minimize foaming).</p>	<p>8.</p>  <p>Add 1 mL of Washing Solution down the inside wall of each tube by using the technique described in Box 2. Wait 2 minutes. Using a smooth motion, invert the combined rack assembly over a sink and pour out the tube contents: keep inverted and gently blot the test tube rims on several layers of paper toweling. Repeat this step two times.</p>																		
<p>3.</p>  <p>React 30 minutes at room temperature (15 °- 30°C).</p>	<p>9.</p>  <p>Lift the upper rack (with its tubes) off the magnetic base; add 500 μL of Color Reagent down the inside wall of each tube by using the technique described in Box 2. Vortex for 1 to 2 seconds (at low speed to minimize foaming).</p>																		
<p>4.</p>  <p>Add 250 μL of Glyphosate Enzyme Conjugate down the inside wall of each tube by aiming the pipet tip 1/4" to 1/2" below the tube rim without touching the rim or tube wall with the pipet tip; deliver liquid gently.</p>	<p>10.</p>  <p>React for 20 minutes at room temperature (15°- 30° C). During this period, add 1 mL of Washing Solution into a clean tube for use as an instrument blank in Step 10.</p>																		
<p>5.</p>  <p>React 30 minutes at room temperature (15 °- 30°C).</p>	<p>11.</p>  <p>Add 500 μL of Stopping Solution down the inside wall of each tube by using the technique previously described. Read results at 450 nm within 15 minutes after adding the Stopping Solution. Multiply results of samples by the appropriate dilution factor (if any).</p> <p>[Safety Caution: Stopping Solution contains diluted sulfuric acid.]</p>																		
<p>6.</p>  <p>Combine the upper rack with the magnetic base; press all tubes into base; allow 2 minutes for the particles to separate.</p>																			

India Contact:
















Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Road No. 44, Pitampura, Delhi – 110034, India
Mobile: +91-98105-21400, Tel: +91-11-42208000, 8111, 8222, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com, www.atzlabs.com ; www.lifetechindia.com



Glyphosate Magnetic Particle Kit Part # 500080, 120 Test

GLYPHOSATE CONCISE FLOWCHART

<p>1.</p>  <p>Separate the rack.</p> <p>Add 300 μL of either Derivatized Standards, Control or Samples to the bottom of each test tube.</p> 	<p>7.</p>  <p>Invert the combined rack. Blot gently.</p>
<p>2.</p>  <p>Add 500 μL of mixed Magnetic Particles to each test tube. Vortex.</p>	<p>8.</p>  <p>Add 1 mL of Washing Solution. Wait 2 minutes. Invert the combined rack. Blot gently. Repeat this step two times</p> 
<p>3.</p>  <p>Incubate for 30 minutes.</p>	<p>9.</p>  <p>Separate the rack. Add 500 μL of Color Reagent to each test tube. Vortex.</p> 
<p>4.</p>  <p>Add 250 μL of Enzyme Conjugate to each test tube. Vortex.</p>	<p>10.</p>  <p>Incubate for 20 minutes. Prepare blank.</p>
<p>5.</p>  <p>Incubate for 30 minutes.</p>	<p>11.</p>  <p>Add 500 μL of Stopping Solution to each test tube. Read OD 450</p> 
<p>6.</p>  <p>Combine the rack and magnetic base. Seat all tubes. Wait 2 minutes.</p>	

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Road No. 44, Pitampura, Delhi – 110034, India
Mobile: +91-98105-21400, Tel: +91-11-42208000, 8111, 8222, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com, www.atzlabs.com ; www.lifetechindia.com

Glyphosate Magnetic Particle Kit Part # 500080, 120 Test



Safety Data Sheet

Section 1: Product and Company Identification

1.1 Product Identifiers:

Product Name: Glyphosate Magnetic Particle Kit, Glyphosate Plate Kit

Product Code: 500081, 500086

1.2 Identified Use: Determination of Glyphosate in samples. **Restrictions on Use:** For research use only.

1.3 Company: Abraxis, Inc., 124 Railroad Drive, Warminster, PA 18974 USA, info@abraxiskits.com +1(215) 357-3911, FAX +1(215) 357-5232

1.4 Emergency Telephone Number: +1(215) 357-3911

Section 2: Hazard(s) Identification

2.1 Classification of the mixture:

GHS Classification in accordance with 29 CFR 1910 (OSHA HCS)

Flammable liquids (Category 4), H227 Combustible liquid

HMIS Rating: Dimethyl sulfoxide, CAS No. 67-68-5: Health hazard: 0, Chronic Health Hazard: *, Flammability: 2, Physical Hazard 0

NFPA Rating: Dimethyl sulfoxide, CAS No. 67-68-5: Health hazard: 0, Fire Hazard: 2, Reactivity Hazard: 0

2.2 GHS Label elements, including precautionary statements:

Pictogram(s): none

Signal word(s): Warning

Hazard statement(s):

H227 Combustible liquid.

Precautionary statement(s):

P210 Keep away from heat/sparks/open flames/hot surfaces. No smoking.

P280 Wear protective gloves/eye protection/face protection.

P370 + P378 In case of fire: Use dry sand, dry chemical or alcohol-resistant foam to extinguish.

P403 + P235 Store in a well-ventilated place. Keep cool.

P501 Dispose of contents/ container to an approved waste disposal plant.

2.3 Hazards not otherwise classified (HNOC) or not covered by GHS: Readily absorbed through skin

2.4 Unknown acute toxicity: None known.

Section 3: Composition / Information on Ingredients

3.2 Mixtures: Mixture(s) of the hazardous substance(s) listed below, with nonhazardous additions.

Hazardous component(s):

Name and Synonym(s): DMSO, Dimethyl sulfoxide, Methyl sulfoxide Formula: C_2H_6OS Molecular weight: 78.13 g/mol

CAS No.: 67-68-5 EC-No.: 200-664-3

Classification: Flammable Liquid 4; H227

Percentage in Mixture: 1.91-3.81 %

For full text of H-Statements mentioned in this Section, see Section 2.

Section 4: First Aid Measures

4.1 Description of first aid measures: Consult a physician. Show this safety data sheet to the doctor in attendance. Move out of dangerous area.

If inhaled: If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult a physician.

In case of skin contact: Wash off with soap and plenty of water. Consult a physician.

In case of eye contact: Rinse thoroughly with plenty of water as a precaution.

If swallowed: Do NOT induce vomiting. Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a physician.

4.2 Most important symptoms and effects, both acute and delayed: The most important known symptoms and effects are described in the labelling (see section 2.2) and/or in section 11

4.3 Indication of any immediate medical attention and special treatment needed: No data available. Treat symptomatically.

Section 5: Fire-fighting Measures

5.1 Suitable extinguishing media: Water spray, alcohol-resistant foam, dry chemical or carbon dioxide

5.2 Special hazards arising from the substance or mixture: Carbon oxides, Sulfur oxides

5.3 Advice for firefighters: Wear self-contained breathing apparatus for fire-fighting if necessary.

5.4 Further information: Use water spray to cool unopened containers.

Section 6: Accidental Release Measures

6.1 Personal precautions, protective equipment and emergency procedures: Use personal protective equipment (see section 8). Avoid dust formation. Avoid breathing vapors, mist, dust, or gas. Ensure adequate ventilation. Beware of vapors accumulating to form explosive concentrations. Vapors can accumulate in low areas. Remove all sources of ignition. Evacuate personnel to safe areas.

6.2 Environmental precautions: Prevent further leakage or spillage if safe to do so. Do not let product enter drains. Discharge into the environment must be avoided.

6.3 Methods and materials for containment and cleaning up: Contain spillage, and then collect with an electrically protected vacuum cleaner or by wet-brushing and place in container for disposal according to local regulations. Keep in suitable, closed containers for disposal.

6.4 Reference to other sections: For information on safe handling see section 7. For information on personal protection see section 8. For information on disposal see section 13.

Section 7: Handling and Storage

7.1 Precautions for safe handling: See section 2. Avoid inhalation of vapor or mist. Keep away from sources of ignition. No smoking. Take measures to prevent the buildup of electrostatic charge.

7.2 Precautions for safe storage: Keep container(s) tightly closed in a dry, well-ventilated place. Store under inert gas. Hygroscopic. See label or product insert for appropriate storage temperature and additional specific information. Storage class (TRGS 510): Combustible liquids.

7.3 Specific end use(s): Other than use(s) specified in section 1, no other uses are stipulated.

Section 8: Exposure Controls / Personal Protection

8.1 Control parameters:

Component(s) with workplace control parameters

Dimethyl sulfoxide, CAS No. 67-68-5

Value	Control parameters	Basis
TWA	250.000000 ppm	USA. Workplace Environmental Exposure Levels (WEEL)

8.2 Exposure controls:

Appropriate engineering controls: Provide adequate ventilation. Handle in accordance with good industrial hygiene and safety practice. Wash hands before breaks and at the end of workday. Keep away from food and beverages.

Personal protective equipment

Eye protection: Use equipment for eye protection with side shields tested and approved under appropriate government standards such as NIOSH (US) or EN 166 (EU).

Skin protection: Handle with chemical resistant gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands.

Respiratory protection: Do not breathe vapors. Use a chemical fume hood or approved respiratory protection equipment.

Body protection: Lightweight, protective clothing to prevent skin exposure.

Section 9: Physical and Chemical Properties

9.1 Information on basic physical and chemical properties of mixture

Appearance: Multiple

Odor: Characteristic

Odor Threshold: No data available

pH: Multiple

Melting point/freezing point: No data available

Initial boiling point and boiling range: No data available

Flash point: No data available

Evaporation rate: No data available

Flammability (solid, gas): No data available

Upper/lower flammability or explosive limits: No data available

Vapor pressure: No data available

Vapor density: No data available

Relative density: No data available

Water solubility: Various

Partition coefficient: n-octanol/water: No data available

Auto-ignition temperature: Not applicable

Decomposition temperature: No data available

Viscosity: No data available

Explosive properties: No data available

Oxidizing properties: No data available

9.2 Other information: No data available

Section 10: Stability and Reactivity

10.1 Reactivity: No data available

10.2 Chemical stability: Stable under recommended storage conditions.

10.3 Possibility of hazardous reactions: No data available

10.4 Conditions to avoid: Keep away from open flame, hot surfaces, heat sources, and sources of ignition.

10.5 Incompatible materials: Acid chlorides, phosphorus halides, strong acids, strong oxidizing agents, strong reducing agents

10.6 Hazardous decomposition products: No data available. In the event of fire: see section 5.

Section 11: Toxicological Information

11.1 Information on toxicological effects

To the best of our knowledge, the chemical, physical, and toxicological properties of this product have not been thoroughly investigated.

Acute toxicity (*Dimethyl sulfoxide*, CAS No. 67-68-5):

Inhalation LC50 Inhalation - Rat - 4 h – 40250 ppm

Ingestion LD50 Oral - Rat – 14,500 mg/kg

Skin contact LD50 Dermal - Rabbit - > 5,000 mg/kg; mild skin irritation; **Eye contact** No data available

Respiratory or skin sensitization No data available; **Aspiration hazard** No data available

Mutagenicity (*Dimethyl sulfoxide*, CAS No. 67-68-5): Cytogenetic analysis (Mouse lymphocyte, Rat)--Result: mutation, DNA damage in mammalian somatic cells

Carcinogenicity:

(*Dimethyl sulfoxide*, CAS No. 67-68-5): Rat (oral)—Tumorigenic: Equivocal tumorigenic agent by RTECS criteria. Skin and Appendages: Other: Tumors. Mouse (oral)—Tumorigenic: Equivocal tumorigenic agent by RTECS criteria. Leukemia Skin and Appendages: Other: Tumors.

IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

ACGIH: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

NTP: No component of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.

OSHA: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.

Teratogenicity: No data available

Reproductive/fertility toxicity:

(*Dimethyl sulfoxide*, CAS No. 67-68-5): Reproductive toxicity (Rat, intraperitoneal and subcutaneous)--Effects on Fertility: Abortion; post-implantation mortality (e.g., dead and/or resorbed implants per total number of implants); litter size (e.g.; # fetuses per litter; measured before birth).

Reproductive toxicity (Mouse, oral)--Effects on Fertility: Pre-implantation mortality (e.g., reduction in number of implants per female; total number of implants per corpora lutea). Effects on Embryo or Fetus: Fetotoxicity (except death, e.g., stunted fetus). Specific Developmental Abnormalities: Musculoskeletal system.

Developmental Toxicity (Mouse, intraperitoneal)--Effects on Embryo or Fetus: Fetotoxicity (except death, e.g., stunted fetus). Specific Developmental Abnormalities: Musculoskeletal system.

Specific target organ toxicity, single exposure (*Dimethyl sulfoxide*, CAS No. 67-68-5): No data available

Specific target organ toxicity, repeated exposure: The substance or mixture is not classified as specific target organ toxicant, repeated exposure.

Additional information (*Dimethyl sulfoxide*, CAS No. 67-68-5): RTECS: PV6210000 Exposure to large amounts can cause redness of skin, itching, burning, sedation, headache, nausea, dizziness. Eyes - Eye disease - Based on Human Evidence

Section 12: Ecological Information

12.1 Toxicity: (*Dimethyl sulfoxide*, CAS No. 67-68-5) Toxicity to fish LC50 - Pimephales promelas (fathead minnow) - 34,000 mg/l - 96 h; LC50 - Oncorhynchus mykiss (rainbow trout) - 35,000 mg/l - 96 h. Toxicity to daphnia and other aquatic invertebrates EC50 - Daphnia magna (Water flea) - 24,600 mg/l - 48 h (OECD Test Guideline 202). Toxicity to algae EC50 - Pseudokirchneriella subcapitata (green algae) - 17,000 mg/l - 72 h (OECD Test Guideline 201)

12.2 Persistence and degradability: Dimethyl sulfoxide result: 31 %, not readily biodegradable (OECD Test Guideline 301D)

12.3 Bioaccumulative potential: No data available

12.4 Mobility in soil: No data available

12.5 Results of PBT and vPvB assessment: No data available

12.6 Other adverse effects: An environmental hazard cannot be excluded in the event of unprofessional handling or disposal.

Section 13: Disposal Considerations

13.1 Waste treatment methods

Product: All waste must be handled and disposed according to local, state, and federal regulations. Avoid disposing large volumes in sewer.

Contaminated packaging: All waste must be handled and disposed according to local, state, and federal regulations.

Refer to sections 7 and 8 for safe handling guidance.

Section 14: Transport Information

DOT, Land Transport ADR/RID (cross-border), Maritime Transport IMDG, Air Transport ICAO-TI and IATA-DGR

UN Number: 3316

UN Proper shipping name: Chemical Kit, (contains DMSO)

Transport hazard class(es): 9

Packing group: III

Environmental hazard: See section 12

Bulk transport: Excepted/Limited quantity

Special considerations: See section 7 for handling

Section 15: Regulatory Information

EU Regulations, Hazard Symbol(s): Dimethyl sulfoxide: Xi (Irritant)

Safety Phrases: Dimethyl sulfoxide: S 24/25 Avoid contact with skin and eyes

SARA Title III, Section 302 Components: No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

SARA Title III, Section 313 Components: This material does not contain any chemical components with known CAS numbers that exceed the threshold (DeMinimis) reporting levels established by SARA Title III, Section 313.

SARA 311/312 Hazards: Dimethyl sulfoxide, CAS No. 67-68-5: Fire Hazard, Chronic Health Hazard

State Right-to-Know

Massachusetts: No components are subject to the Massachusetts Right to Know Act.

Pennsylvania: Dimethyl sulfoxide, CAS No. 67-68-5

New Jersey: Dimethyl sulfoxide, CAS No. 67-68-5

California Prop. 65 Components: This product does not contain any chemicals known to State of California to cause cancer, birth defects, or any other reproductive harm.

Section 16: Other information

This information is based on our present knowledge. While Abraxis , INC. believes that the data contained herein are factual and the opinions expressed represent a best effort to present accurate information, the data are not to be taken as a warranty or representation for which Abraxis , INC. assumes legal responsibility. The information shall not be taken as being all-inclusive and is to be used only as a guide. The data are offered solely for the user's consideration, investigation, and verification. These suggestions should not be confused with either state, municipal, or insurance requirements, or with national safety codes and constitute no warranty. Any use of these data and information must be determined by the user to be in accordance with applicable federal, state, and local regulations.

All materials and mixtures may present unknown hazards and should be used with caution. Since Abraxis , INC. cannot control the methods, volumes, or conditions of use of this product, Abraxis , INC. shall not be held liable for any damages or losses resulting from the handling or from contact with the product as described herein. An individual technically qualified to handle potentially hazardous chemicals must supervise the use of this material. This product is sold for research use only. It is not for any human or animal therapeutic or clinical diagnostic use.

Date this SDS was prepared: 5/20/2016

Version: 3

Changes from previous version: Abraxis, LLC changed to Abraxis, Inc.