

Product Data Sheet

Wash Buffer Concentrate (100X) for ELISA

Cat. No. WB-150 **Size:** 50 ml **Form:** 100X Solution (makes 5 Liter solution)
Storage: Store at 4⁰C

Introduction

Antigen or antibody coating on ELISA plates is the first step in ELISA. This is followed by washing off the unbound antigens. This step is also necessary after the incubation of plates with primary antibody and enzyme conjugates. This wash buffer can be used to remove unbound or non-specific binding of antigen/antibody/conjugates in ELISA. It is also suitable for use in removing non-specific binding in Western blotting procedures.

This buffer formulation is the same as used in ADI autoimmune ELISA kits with cat#WB-100. It is supplied as a separate item for users who need this buffer in bulk.

Form

The Wash buffer concentrate is provided as 100 X solution preserved in 0.02% merthiolate as preservative. Therefore, care should be taken when disposing this solution. Some salts may crystallize upon storage in cold but redissolves after dilution.

Recommended Use

Dilute wash Buffer concentrate 1:100 with distilled water before use. For most ELISA applications, plates should be washed 3-5 times with 300 ul wash buffer/per well. Briefly shake the plate, discard solution, and tap over paper towels to remove traces of liquid. Repeat this step 3-5 times. Alternatively, ELISA Washers can be programmed to wash the plates with 300 ul buffer, with or without soak time (time between washings).

For Western blotting, filters/membranes should be completely immersed in wash buffer (50 ml buffer/filter is sufficient) and shaken on an orbital shaker for 3-5 min. Discard buffer and remove traces of buffer by draining the container. This should be repeated 3-5 times. If background persists, volume of buffer and soak time may be increased up to 15 min.

This product is for in vitro research use only.

MSDS: A Material Safety Data Sheet is not required for this product. The product does not contain any hazardous components above 1% or any carcinogens above 0.1% as defined in 29 CFR 1910.1200, the OSHA Hazard Communication Standard. The product contains natural purified proteins in a buffer containing 0.1% or less sodium azide. MSDS can be consulted for these items. Please note that although the product is defined as not hazardous, it is still advisable to follow prudent laboratory practices when handling laboratory reagents.

This material is sold for research purposes only and is not required to appear on the TSCA inventory. It is not intended for food, diagnostic, drug, household, agricultural or cosmetic use. Its use must be supervised by a technically qualified individual experienced in handling potentially hazardous chemicals.

Related items available from ADI

Catalog#	ProdDescription
80012	High binding ELISA Strips plates (8 wellsx12 strips)
80050	ELISA Plate Coating buffer concentrate (10X)
80060	ELISA Plate Blocking Buffer concentrate (10X) BSA-based
80062	ELISA Plate Blocking Buffer concentrate (10X) milk-based
80070	Antibody and Conjugate Diluent for ELISA concentrate (10X)
80080	Wash buffer concentrate (20X) for ELISA
80081	Phosphate Buffered Saline (PBS, pH 7.4) concentrate (20X)
80091	TMB substrate (1-component) for ELISA
80100	Stop solution for TMB substrate (ELISA) concentrate (10X)
80110	Single solution TMB substrate (membrane) solution for western blot
80160	Rabbit Serum Antibody detection ELISA kit, Qualitative (mouse, goat, sheep, monkey, human, chicken kits are also available)
WB-150	110329A

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com