

Human Endothelial Nitric Oxide Synthase (eNOS) cDNA Probe

Cat#	eNOS51-D	Size: 2.0 ug DNA	Store at -20°C
		Form:	Lyophilized Lyophilized

The cDNA probe is provided as a double stranded cDNA fragment as lyophilized powder. It can be re-suspended in a desired buffer (50-100ul buffer: 10mM Tris-HCl, pH 7.5, 0.1 mM EDTA). Store in suitable size aliquots at -20°C for 6-12 months. The probe has 5' overhangs and >98 % purity. The length of the probe is derived from a segment of human endothelial constitutive nitric oxide synthase (eNOS) gene, which has a length of approximately 1.23 kbp. The probe may be labeled by radioactive or non-radioactive compounds using either random primer labeling or nick translation.

Northern blot analysis reveals an eNOS mRNA with an approximate size of 4.4 kbp.

Form and usage:

For best results, reconstitute the probe in just prior to use in either water or TE buffer. If the reconstituted probe cannot be used within two weeks, it should be aliquoted and into smaller vials and stored at -20°C. It is stable for 1 yr at -20°C. For Northern and Southern blot techniques, we recommend following procedures described in molecular biology manuals such as Current Protocols in Molecular Biology or Molecular Cloning: A Laboratory manual.

Northern blot analyses reveal that eNOS mRNA from bovine and rat endothelial cells has a size of approx. 4.4 kbp (1). In addition, a size of approx. 4.7 kbp has been reported for eNOS mRNA from human bronchiolar epithelial cells, ovine fetal pulmonary artery endothelial cells, and rat sheep, and pig lung. Bovine eNOS is 94% homologous to human eNOS at the amino acid levels. It is only 50-60% related with iNOS and bNOS from various species.

References:

1. Sessa et al. Molecular cloning and expression of a cDNA encoding endothelial cell nitric oxide synthase. *J. Biol. Chem.* **267**, 15274-15276 (1992)
2. Feinburg and Vogelstein. A technique for radiolabeling DNA restriction endonuclease fragments to high specific activity. *Anal. Biochem.* **132**, 6-10 (1983)
3. Koch et al. AN improved method for labeling of DNA probes by nick translation. *Nucleic Acid Res.* **14**, 7132 (1986).
4. Sebi et al. Non Radioactive labeling and detection of nucleic acids III. Application of the digoxigenin system. *Biol. Chem. Hoppe-Seyler* **371**, 939-951 (1990).
5. Kessler et al. Non radioactive labeling and detection of nucleic acids. A novel DNA labeling and detection system based on digoxigenin: Anti-digoxigenin ELISA principle (Digoxigenin system). *Biol. Chem. Hoppe-Seyler* **371**, 916-927 (1990).
6. Shaul et al. Endothelial nitric oxide synthase is expressed in cultured human bronchiolar epithelium. *J. Clin. Invest.* **94**, 2231-2236 (1994).

These materials should be used for In Vitro Research Use and Manufacturing Only.

Related material available from ADI

Anti-bNOS (NOSI); Anti-iNOS (NOSII), Anti-eNOS (NOSIII) Antibodies and Control Peptide
Purified/Recombinant NOSs as Positive Controls.
Anti-Rabbit IgG-HRP conjugates Single solution, ready to use TMB substrates for Blotting

Chemiluminescence Substrates and Western blot kits

eNOS51-D-2 71212A

Alpha Diagnostic Intl Inc., 6203 Woodlake Center Dr, S an Antonio, T X 7 8 24 4 , U S A;

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400 Fax: +91-11-42208444 Email: customerservice@lifetechindia.com Website: www.lifetechindia.com