

Product Data Sheet

**Anti-Ferret IgA (alpha) Antibodies and conjugates**

<b>Cat#</b> 70519	Goat Anti-Ferret IgA, unlabeled	<b>Size:</b> 0.5 mg	<b>Form:</b> Liquid;	Powder
<b>Cat#</b> 70520	Goat Anti-Ferret IgA-HRP conjugate	<b>Size:</b> 0.5 ml	<b>Form:</b> Liquid;	Powder

**Description**

Goats were immunized with antigen grade **Ferret IgA**. Antibodies have been isolated using ammonium sulfate, ion-exchange, and affinity chromatography to remove non-specific antibodies. Specificity of antibody has been tested using IEP, immunodiffusion, and ELISA. The product reacts with Ferret IgA with no significant reactivity with other Ig's or serum proteins. The antibody may recognize other species IgA. Purified anti-Ferret IgA antibody is supplied as unlabeled or conjugated to Biotin, FITC, AP or HRP.

**Form and Storage**

**Cat# 70519, Anti-Ferret IgA, unlabeled**

The antibody is supplied in PBS, pH 7.4, and 0.05% azide in either **lyophilized** (0.5 mg) or **liquid** form (1 mg/ml). Reconstitute powder in PBS in 1.0 ml to prepare 1 mg/ml solution. Store at -20°C in suitable aliquots. Stability is ~6-12 months. Do not freeze and thaw.

**Cat# 70520, Anti-Ferret IgA-HRP-conjugate**

Purified antibody was coupled to HRP (RZ>3.0) using periodate method. The molar enzyme to protein (E/P) ratio = 4.0. The antibody is supplied in stabilizing buffer, 0.1% prolcin-300 as preservative in either **lyophilized** (0.5 ml) or **liquid** form (0.5 ml; 0.5-1 mg/ml). Reconstitute powder in PBS in 1 ml. Store at 4°C in suitable aliquots. Stability is ~6-12 months. Do not freeze and thaw.

Suggested conjugate dilutions are 1:1,000-1:10,000 ELISA, 1:1K-1:5K for western, and 1:200-1:1000 (IHC).

**Recommended Working Dilution for ELISA**

Working dilution for the specific application should be determined by the investigator to obtain the best conditions. Working solution should be prepared immediately before use and diluted solution should be discarded.

**Related Material available for ADI**

**Ferret IgG, IgM ELISA kits**

Anti-Ferret IgG (H+L), IgG, IgG (H+L), IgG (H+L) antibodies and conjugates

ELISA kits for the detection of Ferret Antibodies

**Western blot recycling kit** (Use the same blot/strip to probe with multiple antibodies) **Strip Antibodies in 5 minutes** (no heating or mercaptoethanol)

70519, 70520

70724A