


AmplifyRP[®] Amplicon Detection Chamber



Part Number: ADC 98800/0001

Overview

Intended Use:

Amplicon Detection Chambers are used in conjunction with AmplifyRP[®] Acceler8[®] Discovery and pathogen specific kits.

A lateral flow strip contained inside the chamber detects dual labeled (FAM & Biotin) DNA / RNA amplicons. Amplicon Detection Chambers are specifically designed to provide reliable amplicon detection while preventing their escape and contamination of future assays.

Includes:

- One, single use, amplicon detection chamber

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com