

## ImmunoStrip® for Plantago asiatica mosaic virus (PIAMV)



**Test Label:** NA  
**Test Format:** ImmunoStrip®  
**Capture Reagent:** Polyclonal  
**Detection Reagent:** Polyclonal

### Intended Use:

The PIAMV ImmunoStrip® is an on-site tool to quickly identify PIAMV-infected plants with or without external disease symptoms.

This assay has been confirmed to detect a diverse panel of PIAMV isolates from Asia, Europe, and the United States.

The PIAMV ImmunoStrip requires no special equipment or expertise, making it ideal for use in the field.

### Kit (ISK) Includes:

- ImmunoStrips®
- Sample bags containing SEB1 buffer
- [User Guide](#)

### ImmunoStrips® Only (STX) Includes:

- ImmunoStrips®
- SEB1 extraction buffer purchased separately

#### Option

ImmunoStrip for PIAMV, 5 strips and buffer filled bags  
ISK 91500/0005

ImmunoStrip for PIAMV, 25 strips and buffer filled bags  
ISK 91500/0025

ImmunoStrip for PIAMV, 25 ct \*\*\*BUFFER FILLED BAGS NOT INCLUDED\*\*\*  
STX 91500/0025

#### India Contact:

#### Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)