

ELISA Reagent Set for Potato virus V (PVV)

Test Label: Alkaline Phosphatase

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: Polyclonal Detection Reagent: Polyclonal

Intended Use:

This test is used for detection of Potato virus V (PVV) in foliage and tuber sprouts of potato and

other solanaceous crops, particularly pepino (*Solanum muricatum*) in which the virus is common.

A small amount of cross-reaction can occur with samples that contain *Tobacco etch virus* (TEV).

Includes:

- Capture antibody
- · Alkaline phosphatase enzyme conjugate
- 96-well microtiter plates, solid* (*strip in single plate set)
- User Guide



Reagent Set for Potato virus V (PVV), 96 testwells, SRA 10100/0096

Reagent Set for Potato virus V (PVV), 500 testwells, SRA 10100/0500

Reagent Set for Potato virus V (PVV), 1000 testwells, SRA 10100/1000

Reagent Set for Potato virus V (PVV), 5000 testwells, SRA 10100/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444 Email: customerservice@lifetechindia.com Website: www.lifetechindia.com