

ELISA Reagent Set for Pepino mosaic virus (PepMV)

Test Label: Alkaline Phosphatase

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Polyclonal

Intended Use:

Agdia's PepMV ELISA is used for detection of *Pepino mosaic virus* (PepMV) in several hosts.

Includes:

- Capture antibody
- Alkaline phosphatase enzyme conjugate
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)



Reagent Set for PepMV, 96 testwells,SRA 13001/0096

Reagent Set for PepMV, 500 testwells,SRA 13001/0500

Reagent Set for PepMV, 1000 testwells,SRA 13001/1000

Reagent Set for PepMV, 5000 testwells,SRA 13001/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com