

ELISA Reagent Set for *Acidovorax avenae* subsp. *citrulli* (Aac)

Test Label: Alkaline Phosphatase

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: Monoclonal

Detection Reagent: Monoclonal

Intended Use:

The ELISA and ImmunoStrip® are reliable tools for *Acidovorax avenae* subsp. *citrulli*(Aac), the causal agent of bacterial fruit blotch in watermelon, squash, pumpkins and cantaloupe.

The test can be used with leaves, fruit, bacterial cultures, and seedlings exhibiting symptoms of Aac. This test was designed to produce a low background signal with necrotic seedlings.

Includes:

- Capture antibody
- Alkaline phosphatase enzyme conjugate
- RUB3 enzyme conjugate diluent
- General extract buffer 2, GEB2
- 96-well microtiter plates, solid* (*strip in single plate kit)
- [User Guide](#)



Reagent Set for Aac, 96 testwells,SRA 14800/0096

Reagent Set for Aac, 500 testwells,SRA 14800/0500

Reagent Set for Aac, 1000 testwells,SRA 14800/1000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com



Reagent Set for Aac, 5000 testwells, SRA 14800/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com