

ELISA Reagent Set for Potato virus Y necrotic strain (PVY-n)

Test Label: Alkaline Phosphatase

Test Format: Triple Antibody Sandwich (TAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Monoclonal

Intended Use:

This test is used for detection of isolates of necrotic strain group of PVY. This test also detects the ntn strains that cause the tuber necrotic ringspot disease (PTNRD) of potatoes.

Includes:

- Capture antibody
- Detection antibody, bottle A
- Alkaline phosphatase enzyme conjugate, bottle B
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)
- [Tech Note](#)



Reagent Set for Potato virus Y-necrotic strain (PVY -n), 96 testwells,SRA 26001/0096

Reagent Set for Potato virus Y-necrotic strain (PVY -n), 500 testwells,SRA 26001/0500

Reagent Set for Potato virus Y-necrotic strain (PVY -n), 1000 testwells,SRA 26001/1000

Reagent Set for Potato virus Y-necrotic strain (PVY -n), 5000 testwells,SRA 26001/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com