

ELISA Reagent Set for *Xanthomonas hortorum* pv. *pelargonii* (Xhp)

Test Label: Alkaline Phosphatase

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Monoclonal

Synonym: *Xanthomonas campestris* pv. *pelargonii* (Xcp)

Intended Use:

Using this test system, you can test geranium petioles, stems, leaves and roots directly. No bacterial isolation step is required.

We recommend testing symptomatic tissue if available, or petioles if no symptoms are present. This test system can also be used to test bacterial cultures.

This test system reacts positively with more than 24 pathogenic isolates of *Xanthomonas hortorum* pv. *pelargonii*. The sensitivity of this test system is extraordinarily high. No cross-reaction has been observed with any of 17 other bacterial isolates tested, including *Xanthomonas campestris* pvs. *begoniae*, *campestris*, *dieffenbachiae*, and *phaseoli*.

Includes:

- Capture antibody
- Alkaline phosphatase enzyme conjugate
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)



India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com

Reagent Set for Xanthomonas hortorum pv. pelargonii (Xhp), 96 testwells, SRA 32503/0096

Reagent Set for Xanthomonas hortorum pv. pelargonii (Xhp), 500 testwells, SRA 32503/0500

Reagent Set for Xanthomonas hortorum pv. pelargonii (Xhp), 1000 testwells, SRA 32503/1000

Reagent Set for Xanthomonas hortorum pv. pelargonii (Xhp), 5000 testwells, SRA 32503/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com