



## ELISA Reagent Set for BBWV-1,2

Test Label: Alkaline Phosphatase

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: Polyclonal Detection Reagent: Polyclonal

## Intended Use:

This ELISA is a qualitative serological assay for detection of both *Broad bean wilt virus 1* (BBWV 1) and *Broad bean wilt virus 2* (BBWV 2) in plant foliage.

The test is proven to detect over 29 BBWV-1 and BBWV-2 isolates, but cannot distinguish between them.

The test has been evaluated in over 165 different plant hosts including ornamental and vegetable. Hydrangea, Heliotropium and Salvia were the only species observed to exhibit host effect.

## **Includes:**

- Capture antibody
- Alkaline phosphatase enzyme conjugate
- General extract buffer 3, GEB3
- Tween 20
- 96-well microtiter plates, solid\* (\*strip in single plate set)
- User Guide







Reagent Set for Broad Bean Wilt Virus 1 and 2 (BBWV 1,2), 96 testwells, SRA 46202/0096

Reagent Set for Broad Bean Wilt Virus 1 and 2 (BBWV 1,2), 500 testwells, SRA 46202/0500

Reagent Set for Broad Bean Wilt Virus 1 and 2 (BBWV 1,2), 1000 testwells, SRA 46202/1000

Reagent Set for Broad Bean Wilt Virus 1 and 2 (BBWV 1,2), 5000 testwells, SRA 46202/5000