

AmplifyRP® XRT for *Dickeya* spp.



Part Number: XCS 55500/0048

Test Label: Biotin labeled primer / FAM probe

Test Format: AmplifyRP® XRT

AmplifyRP XRT for *Dickeya* spp. is an isothermal DNA amplification and detection system that rapidly amplifies small portions of DNA conserved among all *Dickeya* spp.

AmplifyRP XRT for *Dickeya* spp. is a real-time isothermal DNA amplification and detection system that offers unrivaled detection capabilities in an easy-to-use testing format. It offers comparable sensitivity and specificity to published PCR methods while eliminating laborious and costly nucleic acid extractions.

The test can be performed virtually anywhere using the battery operated Amplifire fluorometer. Assay parameters are loaded via barcode and results are automatically displayed as (+) or (-).

Prior molecular diagnostic experience is not required to perform AmplifyRP XRT tests.

NOTE: This assay requires a fluorometer to work properly.

Includes:

- XRT reaction pellets for *Dickeya* spp. (48)
- Pre-filled 100 µL PD1 Pellet diluent tubes (48)
- Amp1 extraction buffer (55 mL)
- [User Guide](#)

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com