

AmplifyRP® XRT for HSVd (Hop stunt viroid)



Part Number: XCS 64200/0048

Test Label: FAM probe Test Format: AmplifyRP[®] XRT

AmplifyRP® XRT for HSVd is a real-time isothermal amplification and detection platform designed for field-based or laboratory testing of hops, citrus, prunus, grape and tomato for *Hop stunt viroid*.

The test can be performed virtually anywhere using the battery operated AmpliFire fluorometer. Assay parameters are loaded via barcode and results are automatically displayed as (+) or (-). Prior molecular diagnostic experience is not required to perform AmplifyRP XRT tests.

NOTE: This assay requires a fluorometer to work properly.

Includes:

- XRT reaction pellets for HSVd (48)
- Pre-filled 100 µL PD1 Pellet diluent tubes (48)
- GEB-filled mesh extraction bags (50)
- User Guide