

1kb Sharp DNA Marker



Cat. No. RD001

Characteristics: 5 μ l x 100 reactions

1 kb Sharp DNA Marker: 500 μ l

6X DNA Loading Buffer: 1 ml

Ideal for Size Determination of Digested DNA and PCR Product

Description

Real Biotech Corporation's combined size and mass marker provides superior resolution on agarose gels. The ready-to-load format is highly convenient. Marker is pre-mixed with loading dye and comes with additional 6X loading buffer.

1 kb Sharp DNA Marker contains 9 fragments: 500 bp, 1000 bp, 2000 bp, 3000 bp, 4000 bp, 5000 bp, 6000 bp, 8000 bp, 10000 bp.

Applications

Sharp and clear banding patterns.

Publication quality.

Mass marker bands for easy DNA quantitation.


Accurate DNA migration on agarose gels.

Reference bands for more exact molecular weight determination.

Storage Conditions

1 kb Sharp DNA Marker should be stored immediately upon receipt at -20°C in a constant temperature freezer. 1 kb Sharp DNA Marker can be stored for up to 12 months without showing any deduction in performance and quality with proper storage.

For best results, load 5-10 μ l of marker per well.



India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
 Email: customerservice@lifetechindia.com Website: www.lifetechindia.com