

**HUMAN / ANIMAL ELISA KITS**

Category	Product description	Catalog#	Species Reactivity	Standards/format etc	Suggested applications	
ELISA Kits/Bacterial/Viral/Vaccines	Hepatitis B Surface Antigen (HBsAg) ELISA Kit, Quantitative	4100	Human	Sds=0-10 ng/ml ; Sample=50 ul; Sensitivity= 0.3 ng/ml; 75 min assay; Quantitative, Sandwich, 96 tests	For the measurement of recombinant HBeAg in manufacturing process or vaccine formulations	
ELISA Kits/Human	Human Free T4 (fT4) ELISA Kit, 96 tests, Quantitative	1110	Human	Sds=0-100 pg/ml ; Sample=25 ul; Sensitivity= 1 pg/ml; 75 min assay; Quantitative, Competitive assay; 96 tests		
ELISA Kits/Human	Human T3 (total) ELISA Kit, 96 tests, Quantitative	1700	Human	Sds=0-1000 ng/dL ; Sample=50 ul; Sensitivity= 20 ng/dL; 90 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Free triiodothyronine (fT3) ELISA Kit, 96 tests, Quantitative	1650	Human	Sds=0-40 pg/ml ; Sample=25 ul; Sensitivity= 0.3 pg/ml; 75 min assay; Quantitative, Competitive assay; 96 tests		
ELISA Kits/Human	Human Thyroid Stimulating Hormone (TSH) ELISA Kit, 96 tests, Quantitative	0600	Human	Sds=0-40 mIU/L; Sample=50 ul; Sensitivity= 0.5 mIU/L; 90 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Total Thyroxine (T4) ELISA Kit, 96 tests, Quantitative	1100	Human	Sds=0-25 ug/dL; Sample=25 ul; Sensitivity= 0.5 ug/dL; 90 min assay; Quantitative, Competitive assay		
ELISA Kits/Human	Human Cortisol ELISA Kit, 96 tests, Quantitative	1850	Human	Sds=0-60 ug/dL ; Sample=10 ul; Sensitivity= 0.36 ug/dL; 90 min assay; Quantitative, Competitive assay		
ELISA Kits/Human	Human Progesterone ELISA Kit, 96 tests, Quantitative	1860	Human	Sds=0-30 ng/ml ; Sample=10 ul; Sensitivity= 0.2 ng/ml; 90 min assay; Quantitative, Competitive assay		
ELISA Kits/Human	Human Testosterone ELISA Kit, 96 tests, Quantitative	1880	Human	Sds=0-20 ng/ml ; Sample=10 ul; Sensitivity= 0.125 ng/ml; 90 min assay; Quantitative, Competitive assay		
ELISA Kits/Human	Human Androstenedione ELISA Kit, 96 tests, Quantitative	1910	Human	Sds=0-10 ng/ml ; Sample=25 ul; Sensitivity= 0.05 ng/ml; 75 min assay; Quantitative, Competitive assay		
ELISA Kits/Human	Human Dihydrotestosterone (DHT) ELISA Kit, 96 tests, Quantitative	1940	Human	Sds=0-2500 pg/ml ; Sample=50 ul; Sensitivity= 6 pg/ml; 75 min assay; Quantitative, Competitive assay		
ELISA Kits/Human	Human Pregnenolone ELISA Kit, 96 tests, Quantitative	1865	Human	Sds=0-25 ng/ml; Sample=50 ul; Sensitivity= 0.05 ng/ml; 75 min assay; Quantitative, Competitive assay; 96 tests		
ELISA Kits/Human	Human Estradiol ELISA Kit, 96 tests, Quantitative	1920	Human	Sds=0-3200 pg/ml ; Sample=50 ul; Sensitivity= 10 pg/ml; 70 min assay; Quantitative, Competitive assay		
ELISA Kits/Human	Human Free Testosterone ELISA Kit, 96 tests, Quantitative	1885	Human	Sds=0-100 pg/ml; Sample=25 ul; Sensitivity= 0.17 pg/ml; 75 min assay; Quantitative, Competitive assay; 96 tests		
ELISA Kits/Human	Human Aldosterone ELISA Kit, 96 tests, Quantitative	1875	Human	Sds=0-1000 pg/ml ; Sample=50 ul; Sensitivity= 10 pg/ml; 60 min assay; Quantitative, Competitive assay		
ELISA Kits/Human	Human Androstenedione (salva) ELISA Kit, 96 tests, Quantitative	1915	Human	Sds=0-1000 pg/ml ; Sample=100 ul; Sensitivity= 1 pg/ml; 75 min assay; Quantitative, Competitive assay		
ELISA Kits/Human	Human Progesterone (salva) ELISA Kit, 96 tests, Quantitative	1855	Human	Sds=0-5000 pg/ml ; Sample=50 ul; Sensitivity= 20 pg/mL; 75 min assay; Quantitative, Sandwich assay		
ELISA Kits/Human	Human DHEA-sulphate (DHEA-S) ELISA Kit, 96 tests, Quantitative	1950	Human	Sds=0-10000 ng/ml ; Sample=25 ul; Sensitivity= 5 ng/ml; 60 min assay; Quantitative, Competitive assay; 96 tests		
ELISA Kits/Human	Human Estrone ELISA Kit, 96 tests, Quantitative	1925	Human	Sds=0-2000 pg/ml; Sample=50 ul; Sensitivity= 10 pg/ml; 75 min assay; Quantitative, Competitive assay; 96 tests		
ELISA Kits/Human	Human Cortisol (salva) ELISA Kit, 96 tests, Quantitative	3500	Human	Sds=0-100 ng/ml ; Sample=50 ul; Sensitivity= 1 ng/ml; 65 min assay; Quantitative, Competitive assay; 96 tests		
ELISA Kits/Human	Human 5a-Androstane-3a, 17b-diol Glucuronide (3a-Diol G) ELISA Kit, 96 tests, Quantitative	1905	Human	Sds=0-50 ng/ml ; Sample=50 ul; Sensitivity= 0.1 ng/ml; 45 min assay; Quantitative, Competitive assay		
ELISA Kits/Human	Human Prolactin (PRL) ELISA Kit, 96 tests, Quantitative	0300	Human	Sds=0-100 ng/ml ; Sample=25 ul; Sensitivity= 2 ng/ml; 45 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Chorionic Gonadotropin (hCG) ELISA Kit, 96 tests, Quantitative	0400	Human	Sds=0-200 mIU/ml ; Sample=25 ul; Sensitivity= 3 mIU/ml; 40 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Folic Acid Stimulating Hormone (FSH) ELISA Kit, 96 tests, Quantitative	0200	Human	Sds=0-100 mIU/ml; Sample=50 ul; Sensitivity= 1 mIU/ml; 90 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Luteinizing Hormone (LH) ELISA Kit, 96 tests, Quantitative	0100	Human	Sds=0-200 mIU/ml; Sample=25 ul; Sensitivity= 2 mIU/ml; 40 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Chorionic Gonadotropin (hCG-free beta) ELISA Kit, 96 tests, Quantitative	0410	Human	Sds=0-100 mIU/ml ; Sample=25 ul; Sensitivity= 0.5 mIU/ml; 40 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Insulin ELISA Kit, 96 tests, Quantitative, 96 tests, Quantitative	0030	Human	Sds=0-200 uU/ml ; Sample=25 ul; Sensitivity= 1.5 uU/ml; 45 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human C-peptide ELISA Kit, 96 tests, Quantitative, 96 tests, Quantitative	0040	Human	Sds=0-0.01 nmol/L; Sample=10 ul; Sensitivity= 0.3 ng/ml; 45 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Neopterin ELISA Kit, 96 tests, Quantitative	3410	Human	Sds=0-100 ng/ml; Sample=25 ul; Sensitivity=0.2 ng/ml; 150 min assay; Quantitative; Competitive assay; 96 tests		
ELISA Kits/Human	Human Pancreatic Colorectal cancer (CA-242) ELISA Kit, 96 tests, Quantitative	0060	Human	Sds=0-200 U/ml ; Sample=25 ul; Sensitivity= 1 U/ml; 210 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Alpha Fetoprotein (AFP) ELISA Kit, 96 tests, Quantitative	0500	Human	Sds=0-200 IU/ml ; Sample=25 ul; Sensitivity= 1.0 IU/ml; 75 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Pancreatic & GI Cancer (CA199) ELISA Kit, 96 tests, Quantitative	1840	Human	Sds=0-200 U/ml ; Sample=25 ul; Sensitivity= 3 U/ml; 150 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Ovarian Cancer (CA153) ELISA Kit, 96 tests, Quantitative	1820	Human	Sds=0-200 U/ml ; Sample=25 ul; Sensitivity= 3 U/ml; 150 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Ovarian Cancer (CA125) ELISA Kit, 96 tests, Quantitative	1820	Human	Sds=0-500 U/ml ; Sample=25 ul; Sensitivity= 5 U/ml; 210 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Prostatic Acid Phosphatase (PAP) ELISA Kit, 96 tests, Quantitative	1400	Human	Sds=0-30 ng/L ; Sample=25 ul; Sensitivity= 0.3 ng/L; 40 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Pancreatic Lipase ELISA Kit, 96 tests, Quantitative	1310	Human	Sds=0-400 U/L ; Sample=25 ul; Sensitivity= 1 ng/ml; 40 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Prostate Specific Antigen (PSA) ELISA Kit, 96 tests, Quantitative	1500	Human	Sds=0-60 ng/ml ; Sample=25 ul; Sensitivity= 0.3 ng/ml; 40 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Free Prostate Specific Antigen (PSA-free) ELISA Kit, 96 tests, Quantitative	1510	Human	Sds=0-800 ng/ml ; Sample=10 ul; Sensitivity= 5 ng/ml; 40 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Growth Hormone (GH) ELISA Kit, 96 tests, Quantitative	1600	Human	Sds=0-25 ng/ml ; Sample=25 ul; Sensitivity= 0.2 ng/ml; 90 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Angiogenin ELISA Kit, 96 tests, Quantitative	100-160-ANH	Human	Sds=125-2000 pg/ml ; Sample=100 ul (diluted); Sensitivity= 75 pg/ml; 195 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Neuron Specific Enolase (NSE) ELISA Kit, 96 tests, Quantitative	0050	Human	Sds=0-200 ug/L ; Sample=25 ul; Sensitivity= 1 ug/L; 150 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Sex Hormone Binding Globulin (SHBG) ELISA Kit, 96 tests, Quantitative	0700	Human	Sds=0-300 nmol/L; Sample=20 ul; Sensitivity= 0.2 nmol/L; 60 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Ferritin ELISA Kit, 96 tests, Quantitative	1810	Human	Sds=0-800 ng/ml ; Sample=10 ul; Sensitivity= 5 ng/ml; 40 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Transferrin (Tf) ELISA Kit, 96 tests, Quantitative	1210	Human	Sds=20-400 ng/ml ; Sample=100 ul (diluted); Sensitivity= 10 ng/ml; 105 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human IGF-Binding Protein 1 (IGFBP1) ELISA Kit, 96 tests, Quantitative	0900	Human	Sds=0-250ug/L; Sample=25 ul; Sensitivity= 0.4 ug/L; 70 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Neutrophil Peptide-1 (NP-1/HNHP-1) ELISA Kit, 96 tests, Quantitative	100-270-NP1	Human	Sds = 10 - 160 ng/ml ; Sample = 100 ul; Sensitivity = 5 ng/ml [165 min assay] 1 kit; Quantitative; Sandwich assay; 96 tests		
ELISA Kits/Human	Human Zeta (z) Globulin (alpha thalassaemia) ELISA Kit, 96 tests, Semi-Quantitative	1970	Human	Negative and +ve Calibrators; Sample=100 ul; 75 min; Semi-Quant.; Sandwich; 96 tests		
ELISA Kits/Human	Human globular Adiponectin (gAcrp30) ELISA Kit, 96 tests, Quantitative	200-120-AGH	Human	Sds=0-100 ng/ml ; Sample=50 ul; Sensitivity= 1 ng/ml; 195 min; Quantitative; Sandwich, 96 tests		
ELISA Kits/Human	Human Leptin ELISA Kit, 96 tests, Quantitative, 96 tests, Quantitative	0010	Human	0-50 ng/ml   50 ul   0.5 ng/ml   120 min; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Adiponectin (Acrp30) ELISA Kit, 96 tests, Quantitative	100-140-ADH	Human	Sds=0-12 ng/ml ; Sample=50 ul; Sensitivity= 0.6 ng/ml; 225 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Resistin FIZZ3 ELISA Kit, 96 tests, Quantitative	100-110-RSH	Human	Sds=0-2000 pg/ml ; Sample=20 ul; Sensitivity= 45 pg/ml; 165 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human C-Reactive Protein (CRP) ELISA Kit, 96 tests, Quantitative	1000	Human	Sds=0-100 ng/ml ; Sample=10 ul; Sensitivity= 0.3 ng/ml; 40 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Serum Albumin ELISA Kit, 96 tests, Quantitative	1190	Human	Sds=10-200 ng/ml ; Sample=100 ul (diluted); Sensitivity= 5 ng/ml; 105 min assay; Quantitative; Sandwich assay; 96 tests		
ELISA Kits/Human	Human Beta-2 microglobulin (B2M) ELISA Kit, 96 tests, Quantitative	0020	Human	Sds=0-150 ng/ml ; Sample=10 ul; Sensitivity= 0.5 ng/ml; 90 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Albumin ELISA Kit, 96 tests, Quantitative	1200	Human	Sds=0-100 ug/ml ; Sample=25 ul; Sensitivity= 1 ug/ml; 40 min assay; Quantitative, Competitive assay		
ELISA Kits/Human	Human beta-Defensin-3 (BD-3) ELISA Kit, 96 tests, Quantitative	100-260-BD3	Human	Sds = 50-800 pg/ml ; Sample = 100 ul; Sensitivity = 20 pg/ml [105 min assay] 1 kit; Quantitative; Sandwich assay; 96 tests		
ELISA Kits/Human	Human beta-Defensin-1 (BD-1) ELISA Kit, 96 tests, Quantitative	100-240-BD1	Human	Sds = 100 - 1600 pg/ml ; Sample = 100 ul; Sensitivity = 25 pg/ml [165 min assay] 1 kit; Quantitative; Sandwich assay; 96 tests		
ELISA Kits/Human	Human beta-Defensin-2 (BD-2) ELISA Kit, 96 tests, Quantitative	100-250-BD2	Human	Sds = 12.5-200 pg/ml ; Sample = 100 ul; Sensitivity = 5 pg/ml [155 min assay] 1 kit; Quantitative; Sandwich assay; 96 tests		
ELISA Kits/Human	Human Angiopoietin-2 (Ang-2) ELISA Kit, 96 tests, Quantitative	100-180-APH	Human	Sds=0-8000 pg/ml ; Sample=100 ul; Sensitivity= 200 pg/ml; 225 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Erythropoietin (EPO) Protein ELISA Kit, 96 tests, Quantitative	0070	Human	Sds=0-350 IU/ml; Sample=200 ul; Sensitivity= 1.5 IU/ml (WHO refs calibrated); overnight assay (16-20 hrs); Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Granulocyte Colony Stimulating Factor (G-CSF) Protein ELISA Kit, 96 tests, Quantitative	0080	Human	Sds=0-350 IU/ml; Sample=200 ul; Sensitivity= 1.5 IU/ml (WHO refs calibrated); overnight assay (16-20 hrs); Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Bone Morphogenic Protein 7 (BMP-7) ELISA Kit, 96 tests, Quantitative	100-190-B7H	Human	Sds=100-1600 pg/ml ; Sample=100 ul; Sensitivity= 20 pg/ml; 165 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human IgE ELISA Kit, 96 tests, Quantitative	1800	Human	Sds=0-600 IU/ml ; Sample=10 ul; Sensitivity= 2 IU/ml; 70 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human IgM ELISA Kit, 96 tests, Quantitative	1760	Human	Sds=20-200 ng/ml ; Sample=100 ul; Sensitivity= 5 ng/ml; 105 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human IgG (total) ELISA Kit, 96 tests, Quantitative	1750	Human	Sds=10-200 ng/ml ; Sample=100 ul; Sensitivity= 5 ng/ml; 105 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/Human	Human Rheumatoid Factor IgM (RF) ELISA Kit, 96 tests, Semi-Quantitative	3000	Human	Sds=0, 25, 75 IU/ml calibrators; Sample=100 ul; 70 min; Semi-Quant.; Direct Antigen assay		
ELISA Kits/Human	Human Anti-Nuclear Antibodies (ANA) ELISA Kit, 96 tests, Semi-Quantitative	3200	Human	Negative and +ve Calibrators; Sample=100 ul; 75 min; Semi-Quant.; Direct Antigen assay		
ELISA Kits/Human	Human anti-dsDNA ELISA Kit, 96 tests, Quantitative	3100	Human	Sds=0-300 IU/ml; Sample=100 ul; 75 min; Semi-Quant.; Direct Antigen assay		
ELISA Kits/Human	Human Anti-Mouse Antibodies (HAMA) ELISA Kit, 96 tests, Quantitative	1770	Human	Sds=0-1000 ng/ml ; Sample=25 ul (1:100); Sensitivity= 40 ng/ml; 45 min assay; Quantitative; Sandwich assay; 96 tests		
Cardiac Biomarkers	Human Myoglobin ELISA Kit	600-610-HMY	Human	(5 -1000 ng/ml)		Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Cardiac Biomarkers	Human Cardiac Troponin I (Tn-I) ELISA Kit	600-410-CTN	Human	(2 -75 ng/ml)		Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
ELISA Kits/animals/mouse	Mouse serum amyloid P (SAP) ELISA Kit, 96 tests, Quantitative	6260	Mouse	Sds=12.5-200 ng/ml; Sample=100 ul (diluted); Sensitivity= 10 ng/ml; 165 min assay; Quantitative; Sandwich assay; 96 tests		Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Cardiac Biomarkers	Mouse Myoglobin ELISA Kit	600-630-MMY	Mouse	for serum, plasma & urine samples (3.9 - 250 ng/ml)		Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
ELISA Kits/animals/mouse	Mouse Haptoglobin ELISA Kit, 96 tests, Quantitative	6250	Mouse	Sds=0-125 ng/ml; Sample=100 ul; Sensitivity= 2 ng/ml; 110 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/animals/mouse	Mouse Adiponectin (Acrp30) ELISA Kit, 96 tests, Quantitative	100-150-ADM	Mouse	Sds=0-1500 pg/ml ; Sample=50 ul; Sensitivity= 50 pg/ml; 225 min assay; Quantitative; Sandwich assay; 96 tests		
ELISA Kits/animals/mouse	Mouse Resistin/FIZZ3 ELISA Kit, 96 tests, Quantitative	100-120-RSM	Mouse	Sds=125-2000 pg/ml; Sample=100 ul; Sensitivity= 50 pg/ml; 195 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/animals/mouse	Mouse Leptin ELISA Kit, 96 tests, Quantitative	100-130-LEM	Mouse	Sds=250-4000 pg/ml ; Sample=100 ul; Sensitivity= 100 pg/ml; 225 min assay; Quantitative; Sandwich assay; 96 tests		
ELISA Kits/animals/mouse	Mouse Transferrin (Tf) ELISA Kit, 96 tests, Quantitative	6380	Mouse	Sds=10-200 ng/ml; Sample=100 ul (diluted); Sensitivity= 10 ng/ml; 105 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/animals/mouse	Mouse TNF-alpha ELISA Kit, 96 tests, Quantitative	100-210-TNM	Mouse	Sds=50-1000 pg/ml; Sample=100 ul; Sensitivity= 25 pg/ml; 195 min assay; Quantitative, Sandwich, 96 tests		
Liver Biomarkers	Mouse Liver Fatty acid binding protein (FABP-L) ELISA kit	600-320-LFP	Mouse	(15 - 1000 ng/ml)		Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
ELISA Kits/animals/mouse	Mouse Albumin ELISA Kit, 96 tests, Quantitative	6300	Mouse	Sds=20-400 ng/ml; Sample=100 ul (diluted); Sensitivity= 10 ng/ml; 105 min assay; Quantitative, Sandwich, 96 tests		
ELISA Kits/animals/mouse	Mouse C-Reactive Protein (CRP) reference serum (~1 ug/ml)	1040-RS	mouse	control serum		
ELISA Kits/animals/mouse	Mouse IgG (total) ELISA Kit, 96 tests, Quantitative	6320	Mouse	Sds=10-200 ng/ml; Sample=100 ul (diluted); Sensitivity= 10 ng/ml; 105 min assay; Quantitative, Sandwich, 96 tests		
Acute phase proteins, Liver Biomarkers	Mouse Hemoexin ELISA kit	600-700-HPX	Mouse	(0.156-2.5 ng/ml)		Determination of protein in biological samples (serum, plasma, culture media or purified proteins)

ELISA Kits/animals/mouse	Mouse C-Reactive Protein (CRP) ELISA Kit, 96 tests, Quantitative	1040	mouse	Stds=0-250 ng/ml ; Sample=100 ul (1:100) Sensitivity=2 ng/ml, 125 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/animals/mouse	Mouse Serum Amyloid A (SAA) ELISA Kit, 96 tests, Quantitative	6240	Human	Stds = 7.8 - 500 ng/ml ; Sample = 100 ul (-1:100 diluted); Sensitivity = 2.5 ng/ml 90 min assay 1 Quantitative, Sandwich assay; 96 tests	
Acute phase proteins, Liver Biomarkers	Mouse alpha 2-Macroglobulin ELISA kit	600-720-A2M	Mouse	(0.195 - 6.25 mg/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
ELISA Kits/animals/mouse	Mouse C3 ELISA Kit, 96 tests, Quantitative	6270	Mouse	Stds=12-200 ng/ml; Sample=100 ul (diluted); Sensitivity=5 ng/ml, 105 min assay; Quantitative; Sandwich assay; 96 tests	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Acute phase proteins, Liver Biomarkers	Mouse Clusterin ELISA kit	610-100-MCL	Mouse	(0.39 - 25 ng/ml)	
ELISA Kits/animals/mouse	Mouse VEGF ELISA Kit, 96 tests, Quantitative	100-230-VEM	Mouse	Stds=31-500 pg/ml ; Sample=100 ul; Sensitivity=10 pg/ml; 225 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/animals/mouse	Mouse IgA ELISA Kit, 96 tests, Quantitative	6310	Mouse	Stds=20-400 ng/ml; Sample=100 ul (diluted); Sensitivity=10 ng/ml; 105 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/animals/mouse	Mouse IgM ELISA Kit, 96 tests, Quantitative	6320	Mouse	Stds=10-100 ng/ml; Sample=100 ul (diluted); Sensitivity=5 ng/ml; 105 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/animals/mouse	Mouse IgG2b ELISA Kit, 96 tests, Quantitative	6350	Mouse	Stds=5-100 ng/ml; Sample=100 ul; Sensitivity=2 ng/ml; 105 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/animals/mouse	Mouse IgG3 ELISA Kit, 96 tests, Quantitative	6360	Mouse	Stds=50-1000 ng/ml; Sample=100 ul (diluted); Sensitivity=25 ng/ml; 105 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/animals/mouse	Mouse IgE ELISA Kit, 96 tests, Quantitative	6370	Mouse	Stds=10-200 ng/ml; Sample=100 ul (diluted); Sensitivity=5 ng/ml; 105 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/animals/mouse	Mouse IgG2a ELISA Kit, 96 tests, Quantitative	6340	Mouse	Stds=20-400 ng/ml; Sample=100 ul (diluted); Sensitivity=10 ng/ml; 105 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/animals/mouse	Mouse IgM ELISA Kit, 96 tests, Quantitative	6300	Mouse	Stds=10-200 ng/ml; Sample=100 ul (diluted); Sensitivity=10 ng/ml; 105 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/animals/mouse	Mouse Serum albumin control (level 1 - 1ug/ml) for ELISA	6300-C1	Mouse	pure protein	
Cardiac Biomarkers	Rat & Mouse Cardiac MLC-1 ELISA Kit	610-400-MLC	Mouse	(0.078 - 2.5 ng/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
ELISA Kits/animals/Cat	Cat IgG (total) ELISA Kit, 96 tests, Quantitative	300-100-CGG	Cat	Stds=0-500 ng/ml; Sample=20 ul; Sensitivity=20 ng/ml; 105 min assay; Quantitative; Sandwich assay; 96 tests	
ELISA Kits/animals/mouse	Horse Haptoglobin ELISA Kit, 96 tests, Quantitative	6250-20	Mouse	Stds=0-100 ng/ml; Sample=100 ul; Sensitivity=1 ng/ml; 95 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/animals/Cat	Cat IgM ELISA Kit, 96 tests, Quantitative	300-170-CGM	Cat	Stds=0-1000 ng/ml; Sample=20 ul; Sensitivity=35 ng/ml; 105 min assay; Quantitative; Sandwich assay; 96 tests	
ELISA Kits/animals/mouse	Dog Haptoglobin ELISA Kit, 96 tests, Quantitative	6250-10	Mouse	Stds=0-125 ng/ml; Sample=100 ul; Sensitivity=1 ng/ml; 95 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/Animals/Automune	Rabbit Anti-Myelin Basic Protein (MBP) Ig's ELISA Kit, 96 tests, Quantitative	630-120-MBP	Rabbit	630-110-MBP	
ELISA Kits/Mouse Automune	Rat Anti-Myelin Basic Protein (MBP) Ig's ELISA Kit, 96 tests, Quantitative	630-110-MBP	Rat		
ELISA Kits/animals/Bovine	Bovine IgG (total) ELISA Kit, 96 tests, Quantitative	8010	Bovine	Stds=10-200 ng/ml ; Sample=100 ul; Sensitivity=5 ng/ml; 105 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/animals/Bovine	Bovine IgM ELISA Kit, 96 tests, Quantitative	8080	Bovine	Stds=0-750 ng/ml ; Sample=20 ul; Sensitivity=15 ng/ml; 105 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/animals/Bovine	Bovine Transferrin ELISA Kit, 96 tests, Quantitative	8070	Bovine	Stds=0-4000 ng/ml ; Sample=20 ul; Sensitivity=50 ng/ml; 105 min assay; Quantitative; Sandwich assay; 96 tests	
ELISA Kits/animals/Bovine	Bovine Lactoferrin ELISA Kit, 96 tests, Quantitative	8090	Bovine	Stds=10-160 ng/ml ; Sample=100 ul ; Sensitivity=5 ng/ml ; 105 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/animals/Sheep	Sheep IgG (total) ELISA Kit, 96 tests, Quantitative	7620	Sheep	Stds=0-250 ng/ml ; Sample=20 ul ; Sensitivity=70 ng/ml ; 75 min assay; Quantitative, Competitive assay	
ELISA Kits/animals/Pig	Pig Albumin ELISA Kit, 96 tests, Quantitative	9000	Pig	Stds=0-100 ng/ml ; Sample=100 ul ; Sensitivity=2 ng/ml; 105 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/animals/Pig	Pig IgG (total) ELISA Kit, 96 tests, Quantitative	9020	Pig	Stds=0-1000 ng/ml ; Sample=20 ul ; Sensitivity=20 ng/ml; 105 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/animals/Rat	Rat IgM ELISA Kit, 96 tests, Quantitative	6480	Rat	Stds=0-10,000 ng/ml ; Sample=20 ul; Sensitivity=500 ng/ml; 105 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/animals/Rat	Rat IgG (total) ELISA Kit, 96 tests, Quantitative	6420	Rat	Stds=10-200 ng/ml ; Sample=100 ul ; Sensitivity=5 ng/ml ; 105 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/animals/mouse	Rat C3 ELISA Kit, 96 tests, Quantitative	6230	Mouse	Stds=0-80 ng/ml; Sample=100 ul (diluted); Sensitivity=5 ng/ml; 75 min assay; Quantitative; Sandwich assay; 96 tests	
ELISA Kits/animals/Rat	Rat Albumin ELISA Kit, 96 tests, Quantitative	6400	Rat	Stds=10-200 ng/ml ; Sample=100 ul ; Sensitivity=5 ng/ml; 105 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/animals/Rat	Rat Alpha-1 Glycoprotein (A1-AGP) ELISA Kit, 96 tests, Quantitative	6490	Rat	Stds=0-100 ng/ml ; Sample=100 ul; Sensitivity=15 ng/ml; 120 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/animals/Rat	Rat C-Reactive Protein (CRP) ELISA Kit, 96 tests, Quantitative	1010	Rat	Stds=2-5-60 ng/ml ; Sample=100 ul (diluted); Sensitivity=2 ng/ml; 105 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/animals/Chicken	Chicken IgY ELISA Kit, 96 tests, Quantitative	6030	Chicken	Stds=0-200 ng/ml ; Sample=100 ul ; Sensitivity=5 ng/ml ; 105 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/animals/Turkey	Turkey IgG (total) ELISA Kit, 96 tests, Quantitative	8020	Turkey	Stds=0-100 ug/ml ; Sample=20 ul ; Sensitivity=70 ng/ml ; 75 min assay; Quantitative, Competitive assay	
ELISA Kits/animals/Chicken	Chicken IgM ELISA Kit, 96 tests, Quantitative	6040	Chicken	Stds=0-1000 ng/ml ; Sample=20 ul ; Sensitivity=12 ng/ml ; 105 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/animals/Chicken	Chicken Egg Ovalbumin ELISA Kit, 96 tests, Quantitative	6050	Chicken	Stds=0-4000 ng/ml ; Sample=100 ul ; Sensitivity=100 ng/ml ; 105 min assay; Quantitative; sandwich assay; 96 tests	
ELISA Kits/animals/Chicken	Chicken IgG ELISA Kit, 96 tests, Quantitative	6020	Chicken	Stds=0-1000 ng/ml ; Sample=20 ul ; Sensitivity=15 ng/ml ; 105 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/animals/Chicken	Chicken Egg Ovalbumin ELISA Kit, 96 tests, Quantitative	6010	Chicken	Stds=0-5000 ng/ml ; Sample=20 ul ; Sensitivity=1 ng/ml ; 75 min assay; Quantitative, Competitive assay	
ELISA Kits/animals/Rabbit	Rabbit IgG (total) ELISA Kit, 96 tests, Quantitative	6520	Rabbit	Stds=10-160 ng/ml ; Sample=100 ul ; Sensitivity=5 ng/ml ; 105 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/animals/Rabbit	Rabbit C-Reactive Protein (CRP) ELISA Kit, 96 tests, Quantitative	1030	Rabbit	Stds=0-250 ng/ml ; Sample=100 ul ; Sensitivity=4 ng/ml ; 120 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/animals/Dog	Dog IgM ELISA Kit, 96 tests, Quantitative	400-170-DGM	Dog	Stds=0-800 ng/ml ; Sample=20 ul ; Sensitivity=30 ng/ml ; 105 min assay; Quantitative; Sandwich assay; 96 tests	
ELISA Kits/animals/Dog	Dog C-Reactive Protein (CRP) ELISA Kit, 96 tests, Quantitative	1020	Dog	Stds=0-2000 ng/ml ; Sample=20 ul ; Sensitivity=35 ng/ml ; 105 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/animals/Dog	Dog IgG (total) ELISA Kit, 96 tests, Quantitative	400-100-DDG	Dog	Stds=0-1000 ng/ml ; Sample=20 ul ; Sensitivity=15 ng/ml ; 105 min assay; Quantitative; Sandwich assay; 96 tests	
ELISA Kits/animals/Monkey	Chimp IgE ELISA Kit, 96 tests, Quantitative	7075	Monkey	Stds=0-20 U/ml ; Sample=100 ul ; Sensitivity<1 U/ml; 165 min assay; Quantitative; Sandwich assay; 96 tests	
ELISA Kits/animals/Pig	Pig IgM ELISA Kit, 96 tests, Quantitative	9080	Pig	Stds=0-3000 ng/ml ; Sample=20 ul ; Sensitivity=150 ng/ml; 105 min assay; Quantitative, Sandwich, 96 tests	
Acute phase proteins, Liver Biomarkers	Rat a-1-Acid Glycoprotein ELISA kit	610-410-MLC	Rat	(0.078 - 5 mg/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Cardiac Biomarkers	Monkey Cardiac Troponin 1 (Tn-I) ELISA Kit	600-420-CTN	Monkey	(For use with serum, 0.156 - 10 ng/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Acute phase proteins, Liver Biomarkers	Rat Haptoglobin ELISA kit	6250-30	Rat	(0.098 - 6.25 mg/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Cardiac Biomarkers	Monkey Skeletal Muscle Troponin 1 (Tn-I) ELISA Kit	600-430-MTN	Monkey	(3.1 - 50 ng/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Cardiac Biomarkers	Mouse Cardiac Tn-I ELISA kit for plasma samples	600-440-CTN	Mouse/Rat	(For use with plasma, 0.078 - 5 ng/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Cardiac Biomarkers	Mouse Cardiac Tn-I ELISA kit for serum samples	600-450-CTN	Mouse	(For use with serum, 0.156 - 10 ng/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Acute phase proteins, Liver Biomarkers	Bovine Haptoglobin ELISA kit	6250-60	Bovine	(15.6 - 500 mg/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Cardiac Biomarkers	Pig Cardiac Troponin 1 (Tn-I) ELISA Kit	600-470-CTN	Pig	(For use with serum, 0.039 - 2.5 ng/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Cardiac Biomarkers	Rabbit Cardiac Troponin 1 (Tn-I) ELISA Kit	600-480-CTN	Rabbit	(For use with serum, 0.156 - 10 ng/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Cardiac Biomarkers	Rat Cardiac Troponin 1 (Tn-I) ELISA Kit for serum sample	610-500-CTN	Rat	(For use with serum, 0.156 - 10 ng/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Acute phase proteins, Liver Biomarkers	Pig Haptoglobin ELISA kit	6250-40	Pig	(0.15 - 2.4 mg/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Acute phase proteins, Liver Biomarkers	Rat alpha 2-Macroglobulin ELISA kit	610-420-A2M	Rat	(0.008 - 25 mg/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Cardiac Biomarkers	Rat Skeletal Muscle Troponin 1 (Tn-I) ELISA Kit	600-510-MTN	Rat	(3.1 - 50 ng/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Acute phase proteins, Liver Biomarkers	Rat Hemopexin ELISA kit	610-710-HPX	Rat	(0.098-6.25 mg/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Acute phase proteins, Liver Biomarkers	Rat Clusterin ELISA kit	610-110-RCL	Rat	(0.39 - 25 mg/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Liver Biomarkers	Rat Thioastin ELISA kit	600-730-RTS	Rat	(0.012-0.75 mg/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Cardiac Biomarkers	Rat Myoglobin ELISA Kit	600-660-RMY	Rat	for serum, plasma & urine samples, 1.56 - 100 ng/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Cardiac Biomarkers	Rabbit Myoglobin ELISA Kit	600-650-RMY	Rabbit	for serum, plasma & urine samples, 2 - 125 ng/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Cardiac Biomarkers	Pig Myoglobin ELISA Kit	600-640-PMY	Pig	for serum, plasma & urine samples, 0.5 - 31.25 ng/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Cardiac Biomarkers	Dog Myoglobin ELISA Kit	600-600-DMY	Dog	for serum, plasma & urine samples, 3.9 - 250 ng/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Cardiac Biomarkers	Monkey Myoglobin ELISA Kit	600-620-MMY	Monkey	for serum, plasma & urine samples, 2 - 125 ng/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Cardiac Biomarkers	Rat Cardiac Troponin 1 (Tn-I) ELISA Kit for plasma samples	610-490-CTN	Rat	(For use with plasma, 0.078 - 5 ng/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Cardiac Biomarkers	Pig Cardiac Fatty acid binding protein (FABP) ELISA kit	600-360-CFP	Pig	(0.78 - 25 ng/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Cardiac Biomarkers	Mouse Skeletal Muscle Troponin 1 (Tn-I) ELISA Kit	600-460-MTN	Mouse	(3.1 - 50 ng/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Acute phase proteins, Liver Biomarkers	Cat Haptoglobin ELISA kit	6250-50	Cat	(0.030 - 2.5 mg/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Cardiac Biomarkers	Dog Cardiac Troponin 1 (Tn-I) ELISA Kit	600-400-CTN	Dog	(For use with serum, 0.156 - 10 ng/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Cardiac Biomarkers	Rat Cardiac Fatty acid binding protein (FABP) ELISA kit	600-330-CFP	Rat	(0.078 - 5 ng/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
ELISA Kits/animals/mouse	Rat serum amyloid P (SAP) ELISA Kit, 96 tests, Quantitative	6410	Rat	Stds=0-160 ng/ml; Sample=100 ul (diluted); Sensitivity=1.25 ng/ml; 75 min assay; Quantitative; Sandwich assay; 96 tests	
Cardiac Biomarkers	Dog Cardiac Fatty acid binding protein (FABP) ELISA kit	600-380-CFP	Dog	(0.31 - 20 ng/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Liver Biomarkers	Rat Liver Fatty acid binding protein (FABP-L) ELISA kit	600-340-LFP	Rat	(15 - 1000 ng/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)

Cardiac Biomarkers	Rabbit Cardiac Fatty acid binding protein(FABP) ELISA kit	600-350-CFP	Rabbit	(1.56 - 100 ng/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
Cardiac Biomarkers	Mouse Cardiac Fatty acid binding protein (FABP) ELISA kit	600-310-CFP	Mouse	(0.31 - 20 ng/ml)	Determination of protein in biological samples (serum, plasma, culture media or purified proteins)
ELISA Kits/Animals/Monkey	Monkey IgG ELISA Kit, 96 tests, Quantitative	7010	Monkey	Stds=200-3000 ng/ml ; Sample=100 ul (diluted); Sensitivity: 100 ng/ml; 105 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/Animals/Monkey	Monkey Haptoglobin ELISA Kit, 96 tests, Quantitative	1060	Monkey	Stds=0-156.25 ng/ml ; Sample=100 ul (diluted -1:10,000); Sensitivity: 0.75 ng/ml; 110 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/Animals/Monkey	Monkey C-Reactive Protein (CRP) ELISA Kit, 96 tests, Quantitative	1050	Monkey	Stds=0-150 ng/ml ; Sample=100 ul; Sensitivity= 3 ng/ml; 135 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/Animals/Monkey	Monkey Alpha-1 acid glycoprotein ELISA Kit, 96 tests, Quantitative	1070	Monkey	Stds=3-100 ng/ml; Sample=100 ul (diluted -1:20,000); Sensitivity= 1.5 ng/ml; 110 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/Animals/Monkey	Monkey IgG (total) ELISA Kit, 96 tests, Quantitative	7050	Monkey	Stds=10-200 ng/ml ; Sample=100 ul; Sensitivity= 5 ng/ml; 105 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/Animals/Monkey	Monkey IgE ELISA Kit, 96 tests, Quantitative	7070	Monkey	Stds=0.4 U/ml ; Sample=100 ul; Sensitivity=0.5 U/ml; 165 min assay; Quantitative, Sandwich assay, 96 tests	
ELISA Kits/Animals/Monkey	Monkey IgM ELISA Kit, 96 tests, Quantitative	7060	Monkey	Stds=50-1000 ng/ml ; Sample=100 ul (diluted); Sensitivity=25 ng/ml; 105 min assay; Quantitative, Sandwich, 96 tests	
ELISA Kits/Mouse Autoimmune	Mouse RF Ig's (total (A+G+M) ELISA kit, 96 Tests, Qualitative	6200	Mouse	Stds=ve and +ve controls; Sample=100 ul (1:100); 75 min assay; Qualitative; Ag-based Ab assay, 96 tests	
ELISA Kits/Mouse Autoimmune	Monkey Anti-Histones Ig's ELISA Kit, 96 tests, Quantitative	670-120-HIM	Monkey		
ELISA Kits/Mouse Autoimmune	Mouse Anti-SSB Ig's (total (A+G+M) ELISA Kit, 96 tests, Quantitative	5810	Mouse	Stds=0-75 ng/ml; Sample: 100 ul (1:100); Sensitivity: 2 ng/ml; 75 min assay; Quantitative; Ag-based Ab assay, 96 tests	
ELISA Kits/Mouse Autoimmune	Mouse Anti-RNP Ig's (total (A+G+M) ELISA Kit, 96 tests, Qualitative	6410	Mouse	Stds=0-4000 ng/ml; Sample: 100 ul (1:100); Sensitivity: 150 ng/ml; 75 min assay; Quantitative; Ag-based Ab assay, 96 tests	
ELISA Kits/Mouse Autoimmune	Mouse Anti-SSA/Ro Ig's (total (A+G+M) ELISA Kit, 96 tests, Quantitative	5710	Mouse	Stds=0-4000 ng/ml; Sample: 100 ul (1:100); Sensitivity: 125 ng/ml; 75 min assay; Quantitative; Ag-based Ab assay, 96 tests	
ELISA Kits/Mouse Autoimmune	Mouse Anti-esDNA IgM-specific ELISA Kit, 96 tests, Quantitative	5330	Mouse	Stds=0-125 ng/ml; Sample: 100 ul (1:100); Sensitivity: 2 ng/ml; 75 min assay; Quantitative; Ag-based Ab assay, 96 tests	
ELISA Kits/Mouse Autoimmune	Mouse Circulating Immune Complexes (CIC) Ig's (total (A+G+M) ELISA kit, 96 Tests, Qualitative	5900	Mouse	Stds=ve and +ve controls; Sample=100 ul (1:100); 75 min assay; Qualitative; Ag-based Ab assay, 96 tests	
ELISA Kits/Mouse Autoimmune	Mouse anti-dsDNA IgM-specific ELISA Kit, 96 tests, Quantitative	5130	Mouse	Stds=0-2500 ng/ml; Sample: 100 ul (1:100); Sensitivity: 75 ng/ml; 75 min assay; Quantitative; Ag-based Ab assay, 96 tests	
ELISA Kits/Mouse Autoimmune	Mouse Anti-dsDNA Ig's (total (A+G+M) ELISA kit, 96 tests, Quantitative	5310	Mouse	Stds=0-5000 ng/ml; Sample: 100 ul (1:100); Sensitivity: 150 ng/ml; 75 min assay; Quantitative; Ag-based Ab assay, 96 tests	
ELISA Kits/Mouse Autoimmune	Mouse Anti-Histones Ig's (total (A+G+M) ELISA Kit, 96 tests, Quantitative	5610	Mouse	Stds=0-4000 ng/ml; Sample: 100 ul (1:100); Sensitivity: 125 ng/ml; 75 min assay; Quantitative; Ag-based Ab assay, 96 tests	
ELISA Kits/Mouse Autoimmune	Mouse anti-Scl70 Ig's (total (A+G+M) ELISA Kit, 96 tests, Quantitative	6110	Mouse	Stds=0-250 ng/ml; Sample: 100 ul (1:100); Sensitivity: 10 ng/ml; 75 min assay; Quantitative; Ag-based Ab assay, 96 tests	
ELISA Kits/Mouse Autoimmune	Mouse anti-dsDNA IgG-specific ELISA Kit, 96 tests, Quantitative	5120	Mouse	Stds=0-5000 ng/ml; Sample: 100 ul (1:100); Sensitivity: 75 ng/ml; 75 min assay; Quantitative; Ag-based Ab assay, 96 tests	
ELISA Kits/Mouse Autoimmune	Mouse Anti-dsDNA Ig's (Total A+G+M) ELISA Kit, 96 tests, Quantitative	5110	Mouse	Stds=0-4650 ng/ml; Sample: 100 ul (1:100); Sensitivity: 400 ng/ml; 75 min assay; Quantitative; Stds=0-5000 ng/ml; Sample: 100 ul (1:100); Sensitivity: 75 ng/ml; 75 min assay, 96 tests	
ELISA Kits/Mouse Autoimmune	Monkey Anti-Ovalbumin Ig's ELISA Kit, 96 tests, Quantitative	670-130-OVM	Monkey		
ELISA Kits/Mouse Autoimmune	Mouse Anti-Nuclear Antibodies (ANA) Ig's (total (A+G+M) ) ELISA Kit, 96 tests, Quantitative	5210	Mouse	Stds=0-10 ug/ml; Sample: 100 ul (1:100); Sensitivity: 0.2 ng/ml; 75 min assay; Quantitative; Ag-based Ab assay, 96 tests	
ELISA Kits/Mouse Autoimmune	Mouse Anti-dsDNA IgG-specific ELISA Kit, 96 tests, Quantitative	5320	Mouse	Stds=0-10,000 ng/ml; Sample: 100 ul (1:100); Sensitivity: 300 ng/ml; 75 min assay; Quantitative; Ag-based Ab assay, 96 tests	
ELISA Kits/Animals/Rat	Rat Anti-Myelin Oligodendrocyte protein (MOG) Ig's ELISA Kit, 96 tests, Quantitative	600-240-MOG	Rat		
ELISA Kits/Animals/Mouse	Mouse Anti-Myelin Oligodendrocyte protein (MOG35-55) Ig's ELISA Kit, 96 tests, Quantitative	600-200-MOG	Mouse	Stds=0-1000 ng/ml ; Sample=100 ul (1:100); Sensitivity= 50 ng/ml; 75 min assay; Quantitative, Ag-based Ab Assay, 96 tests	
ELISA Kits/Mouse Autoimmune	Mouse ANA IgG +ve control serum; standard F (10 ug/ml)	5217	Mouse	Reference serum	
ELISA Kits/Mouse Autoimmune	Mouse Anti-Sm Ig's (total (A+G+M) ELISA Kit, 96 tests, Quantitative	5405	Mouse	Stds=0-1000 ng/ml; Sample=100 ul; Sensitivity=25 ng/ml; 75 min assay; Quantitative; Ag-based Ab Assay, 96 tests	
ELISA Kits/Mouse Autoimmune	Mouse Anti-Jo-1 Ig's (total) ELISA Kit, 96 tests, Quantitative	6005	Mouse	Stds=0-500 ng/ml; Sample=100 ul (1:100); Sensitivity: 15 ng/ml; 75 min assay; Quantitative; Ag-based Ab Assay, 96 tests	
ELISA Kits/Mouse Autoimmune	Mouse Anti-C1q Ig's ELISA kit, 96 Tests, Qualitative	5020	Mouse	Mouse Anti-C1q Ig's ELISA Kit; Qualitative; Ag-based Ab assay, 96 tests	
ELISA Kits/Mouse Autoimmune	Mouse Anti-PCNA Ig's ELISA kit, 96 Tests, Quantitative	5010	Mouse	Stds=0-50 ng/ml; Sample: 100 ul (1:100); Sensitivity: 1 ng/ml; 75 min assay; Qualitative; Ag-based Ab assay, 96 tests	
ELISA Kits/Mouse Autoimmune	Monkey Anti-dsDNA Ig's ELISA Kit, 96 tests, Quantitative	670-100-DNM	Monkey		
ELISA Kits/Animals/Mouse	Mouse Anti-Myelin Oligodendrocyte protein (MOG) Ig's ELISA Kit, 96 tests, Quantitative	600-230-MOG	Mouse		
ELISA Kits/Mouse Autoimmune	Mouse Anti-beta-2 Glycoprotein I Apopt Ig's ELISA Kit, 96 Tests, Quantitative	5000	Mouse		
ELISA Kits/Mouse Autoimmune	Mouse Anti-Myelin Basic Protein (MBP) Ig's ELISA Kit, 96 tests, Quantitative	630-100-MBP	Mouse	Stds=0-500 ng/ml; Sample: 100 ul (1:100); Sensitivity: 15 ng/ml; 75 min assay; Qualitative; Ag-based Ab assay; 96 tests	
ELISA Kits/Mouse Autoimmune	Mouse Anti-Myelin Associated Glycoprotein (MAG) Ig's ELISA Kit, 96 tests, Quantitative	640-100-MAG	Mouse		
ELISA Kits/Mouse Autoimmune	Mouse Anti-Myelin Oligodendrocyte protein (MOG35-55) Ig's ELISA Kit, 96 Tests, Qualitative	4000	Mouse	Stds=0-500 ng/ml; Sample: 100 ul (1:100); Sensitivity: 15 ng/ml; 75 min assay; Qualitative; Ag-based Ab assay; 96 tests	
ELISA Kits/Mouse Autoimmune	Mouse Anti-Cardiolipin Ig's -ve control for ELISA	5502	Mouse		
ELISA Kits/Animals/Rabbit	Rabbit Anti-Myelin Oligodendrocyte protein (MOG35-55) Ig's ELISA Kit, 96 tests, Quantitative	600-210-MOG	Mouse	Stds=0-1000 ng/ml ; Sample=100 ul (1:100); Sensitivity= 50 ng/ml; 75 min assay; Quantitative, Ag-based Ab Assay, 96 tests	
ELISA Kits/Mouse Autoimmune	Mouse Anti-Cardiolipin Ig's -ve control for ELISA	5503	Mouse		
ELISA Kits/Mouse Autoimmune	Monkey Anti-ANA Ig's ELISA Kit, 96 tests, Quantitative	670-110-ANM	Monkey		
ELISA Kits/Mouse Autoimmune	Mouse Anti-Glutamic Acid Decarboxylase 65 (GAD67) Ig's ELISA Kit, 96 tests, Quantitative	660-110-G67	Mouse		
ELISA Kits/Mouse Autoimmune	Mouse Anti-Glutamic Acid Decarboxylase 65 (GAD65) Ig's ELISA Kit, 96 tests, Quantitative	660-100-G65	Mouse		
ELISA Kits/Animals/Mouse	Mouse Anti-Ovalbumin IgG1 ELISA Kit, 96 tests, Quantitative	600-110-OG1	Mouse	Stds=0-50 ng/ml ; Sample=100 ul (1:100); Sensitivity= 1 ng/ml; 105 min assay; Quantitative, Ag-based Ab Assay, 96 tests	
ELISA Kits/Animals/Mouse	Mouse Anti-Ovalbumin IgG2a ELISA Kit, 96 tests, Quantitative	600-120-O2A	Mouse	Stds=0-300 ng/ml ; Sample=100 ul (1:100); Sensitivity= 5 ng/ml; 105 min assay; Quantitative, Ag-based Ab Assay, 96 tests	
ELISA Kits/Animals/Mouse	Mouse Anti-Ovalbumin Ig's (total (A+G+M) ELISA Kit, 96 tests, Quantitative	600-100-OGG	Mouse	Stds=0-400 ng/ml ; Sample=100 ul (1:100); Sensitivity= 5 ng/ml; 105 min assay; Quantitative, Ag-based Ab Assay, 96 tests	
ELISA Kits/Animals/Mouse	Mouse Anti-Ovalbumin IgG2b ELISA Kit, 96 tests, Quantitative	600-130-O2B	Mouse	Stds=0-200 ng/ml ; Sample=100 ul (1:100); Sensitivity= 10 ng/ml; 105 min assay; Quantitative, Ag-based Ab Assay, 96 tests	
ELISA Kits/Animals/Mouse	Mouse anti-Dinitrophenol (DNP) IgG ELISA Kit, 96 tests, Quantitative	600-120-DGO	Mouse	Stds=0-60 ng/ml ; Sample=50 ul; Sensitivity= 20 ng/ml; 105 min assay; Quantitative, Ag-based Ab Assay, 96 tests	
ELISA Kits/Animals/Mouse	Mouse Anti-Ovalbumin IgM ELISA Kit, 96 tests, Quantitative	600-170-OGM	Mouse	Stds=0-300 ng/ml ; Sample=100 ul (1:100); Sensitivity= 10 ng/ml; 105 min assay; Quantitative, Ag-based Ab Assay, 96 tests	
ELISA Kits/Animals/Mouse	Mouse Anti-Ovalbumin IgG (IgG specific) ELISA Kit, 96 tests, Quantitative	600-105-OGG	Mouse	Stds=0-300 ng/ml ; Sample=100 ul (1:100); Sensitivity= 10 ng/ml; 105 min assay; Quantitative, Ag-based Ab Assay, 96 tests	
ELISA Kits/Animals/Mouse	Mouse Anti-Dinitrophenol (DNP) IgE ELISA Kit, 96 tests, Quantitative	640-200-DGE	Mouse	Stds=0-500 ng/ml ; Sample=50 ul; Sensitivity= 20 ng/ml; 105 min assay; Quantitative, Ag-based Ab Assay, 96 tests	
ELISA Kits/Animals/Mouse	Mouse Anti-Ovalbumin IgE ELISA Kit, 96 tests, Quantitative	600-165-OGE	Mouse	Stds=0-250 ng/ml ; Sample=100 ul (1:100); Sensitivity= 5 ng/ml; 135 min assay; Quantitative; Ag-based Ab Assay; 96 tests	
ELISA Kits/Animals/Rat	Rat Anti-Dinitrophenol (DNP) IgG (total) ELISA Kit, 96 tests, Quantitative	600-110-DGO	Rat	Rat Anti-Dinitrophenol (DNP) IgG (total) ELISA Kit; Quantitative; Ag-based Ab Assay, 96 tests	
ELISA Kits/Animals/Rat	Rat Anti-Ovalbumin IgG (total Ig's A+G+M) ELISA Kit, 96 tests, Quantitative	610-100-OGG	Rat	Stds=0-5000 ng/ml ; Sample=20 ul; Sensitivity= 60 ng/ml; 105 min assay; Quantitative, Ag-based Ab Assay, 96 tests	
ELISA Kits/Animals/Rat	ELISA Kit/Animals/Rat	600-100-OGG	Rat	Quantitative, Ag-based Ab Assay, 96 tests	
ELISA Kits/Animals/Mouse	Mouse Anti-Ovalbumin IgG3 ELISA Kit, 96 tests, Quantitative	600-140-OG3	Mouse	Stds=0-300 ng/ml ; Sample=100 ul (1:100); Sensitivity= 105 ng/ml; 105 min assay; Quantitative, Ag-based Ab Assay, 96 tests	
ELISA Kits/Animals/Rabbit	Rabbit Anti-Ovalbumin IgG (total Ig's A+G+M) ELISA Kit, 96 tests, Quantitative	610-100-OGG	Rabbit		
ELISA Kits/Animals/Goat	Goat IgG (total) ELISA Kit, 96 tests, Quantitative	7520	Goat	Stds=0-400 ug/ml ; Sample=20 ul; Sensitivity= 250 pg/ml; 75 min assay; Quantitative, Competitive assay	
ELISA Kits/Animals/Chicken, ELISA Kits Bioprocess contamination	Chicken IgY Antibody detection in whole egg yolks ELISA kit, Qualitative	80177	Chicken	ELISA Kit for the detection of Chicken Antibodies in serum, 75 min; Qualitative, Direct Antigen based assay	
ELISA Kits/Animals/Chicken, ELISA Kits Bioprocess contamination	Chicken Serum Antibody detection ELISA kit, Qualitative	80176	Chicken	ELISA Kit for the detection of Chicken Antibodies in serum, 75 min; Qualitative, Direct Antigen based assay	
ELISA Kits/Animals/Rabbit, ELISA Kits Bioprocess contamination	Rabbit Serum Antibody detection ELISA kit, Qualitative	80160	Rabbit	ELISA Kit for the detection of Rabbit Antibodies in serum, 75 min; Qualitative, Direct Antigen based assay	
ELISA Kits/Animals/Bovine, ELISA Kits Bioprocess contamination	Bovine Serum Antibody detection ELISA kit, Qualitative	80165	Bovine	ELISA Kit for the detection of Bovine Antibodies in serum, 75 min; Qualitative, Direct Antigen based assay	
ELISA Kits/Animals/Rabbit, ELISA Kits Bioprocess contamination	Rabbit Anti-BSA IgG (total) ELISA Kit, 96 tests, Quantitative	710-110-BSR	Rabbit	Stds=0-20 ug/ml ; Sample=50 ul (1:100); Sensitivity= 25 ng/ml; 75 min; Quantitative, Ag-based Ab Assay, 96 tests	
ELISA Kits/Animals/Goat, ELISA Kits Bioprocess contamination	Human Anti-BSA IgG (total) ELISA Kit, 96 tests, Quantitative	710-140-BSM	Human		
ELISA Kits Bioprocess contamination	Mouse anti-Bovine fetuin Ig's ELISA kit	800-180-BFM	Mouse		
ELISA Kits Bioprocess contamination	Mouse Brain proteins (total) ELISA kit, 96 tests	800-190-MBN	Mouse		
ELISA Kits Bioprocess contamination	Rabbit anti-Bovine fetuin Ig's ELISA kit	800-170-BFR	Rabbit		
ELISA Kits/Animals/Rabbit, ELISA Kits Bioprocess contamination	Rabbit Anti-GST IgG (total) ELISA Kit, 96 tests, Quantitative	720-110-GSR	Rabbit	Rabbit Anti-GST IgG (total) ELISA Kit; Quantitative, Ag-based Ab Assay, 96 tests	
ELISA Kits/Animals/Mouse, ELISA Kits Bioprocess contamination	Mouse Anti-BSA IgG (total) ELISA Kit, 96 tests, Quantitative	710-130-BSM	Mouse		
ELISA Kits Bioprocess contamination	Protein-A ELISA Kit, 96 tests, Quantitative	800-110-PRA	Chicken	Stds=125-2000 pg/ml ; Sample=100 ul; Sensitivity= 50 pg/ml; 105 min assay; Quantitative, sandwich, 96 tests	
ELISA Kits Bioprocess contamination	Protein-G ELISA Kit, 96 tests, Quantitative	800-120-PRG	Chicken	Stds=0-40 ng/ml ; Sample=100 ul; Sensitivity= 2 ng/ml; 105 min assay; Quantitative, Ag-based Ab Assay, 96 tests	
ELISA Kits/Animals/Rabbit, ELISA Kits Bioprocess contamination	Keyhole Limpet Hemocyanin (KLH) coated ELISA plate, 96 wells/5 plates/pk	700-101	ELISA plates		
ELISA Kits/Animals/Rabbit, ELISA Kits Bioprocess contamination	Rabbit Anti-KLH IgG (total) ELISA Kit, 96 tests, Quantitative	700-110-KLR	Rabbit	Stds=0-1000 ng/ml ; Sample=100 ul (1:100); Sensitivity= 50 ng/ml; 75 min; Quantitative, Ag-based Ab Assay, 96 tests	
ELISA Kits/Animals/Mouse, ELISA Kits Bioprocess contamination	Mouse Anti-Vacuum/Immucolite (KLH) IgG (total) ELISA Kit, 96 tests	700-150-VAM	Mouse		
ELISA Kits/Animals/Bovine	Bovine Albumin (BSA) ELISA Kit, 96 tests, Quantitative	8000	Bovine	Stds=0-15.40 ug/ml ; Sample=100 ul; Sensitivity= 100 ng/ml; 75 min assay; Quantitative, Competitive assay	
ELISA Kits/Animals/Bovine, ELISA Kits Bioprocess contamination	Human Anti-Vacuum/Immucolite (KLH) IgG (total) ELISA Kit, 96 tests	700-160-VAH	Mouse		
ELISA Kits/Animals/Chicken, ELISA Kits Bioprocess contamination	Chicken Anti-BSA IgG (total) ELISA Kit, 96 tests, Quantitative	710-120-BSC	Chicken	Stds=0-10000 ng/ml ; Sample=100 ul (1:100); Sensitivity= 50 ng/ml; 75 min; Quantitative, Ag-based Ab Assay, 96 tests	
ELISA Kits/Animals/Goat, ELISA Kits Bioprocess contamination	Monkey Anti-KLH IgG ELISA Kit, 96 tests, Quantitative	700-170-KLG	Monkey	6.25-100 units/ml	
ELISA Kits/Animals/Goat, ELISA Kits Bioprocess contamination	Goat Anti-KLH IgG (total) ELISA Kit, 96 tests, Quantitative	700-100-KLG	Goat	Stds=0-10000 ng/ml ; Sample=100 ul (1:100); Sensitivity= 25 ng/ml; 75 min; Quantitative, Ag-based Ab Assay, 96 tests	
ELISA Kits/Animals/Mouse, ELISA Kits Bioprocess contamination	Mouse Anti-KLH IgG (total) ELISA Kit, 96 tests, Quantitative	700-130-KLM	Mouse		
ELISA Kits/Animals/Bovine	Bovine Fetuin (Alpha-2 HS-Glycoprotein, AHSG, A2HS), BSE-TSE free (Austrian Origin), Low endotoxin	FETB16-N-1	Bovine	BSE-TSE free, Low endotoxins	
ELISA Kits/Animals/Chicken, ELISA Kits Bioprocess contamination	Monkey Anti-KLH IgM ELISA Kit, 96 tests, Quantitative	700-180-KLM	Monkey	6.25-100 units/ml	
ELISA Kits/Animals/Goat, ELISA Kits Bioprocess contamination	Goat Anti-BSA IgG (total) ELISA Kit, 96 tests, Quantitative	710-100-BSG	Goat	Stds=0-1000 ng/ml ; Sample=50 ul (1:100); Sensitivity= 40 ng/ml; 105 min; Quantitative, Ag-based Ab Assay, 96 tests	
ELISA Kits/Animals/Bovine	Bovine Fetuin (Alpha-2 HS-Glycoprotein, AHSG, A2HS), BSE-TSE free (New Zealand Origin), Low endotoxin	FETB15-N-1	Bovine	BSE-TSE free, Low endotoxins; 3CFR bacteria and virus tested	
ELISA Kits/Animals/Chicken, ELISA Kits Bioprocess contamination	Chicken Anti-KLH IgG (total) ELISA Kit, 96 tests, Quantitative	700-120-LCK	Chicken	Stds=0-1000 ng/ml ; Sample=100 ul (1:100); Sensitivity= 100 ng/ml; 75 min; Quantitative, Ag-based Ab Assay, 96 tests	
ELISA Kits/Animals/Chicken, ELISA Kits Bioprocess contamination	Chicken Anti-GST IgG (total) ELISA Kit, 96 tests, Quantitative	720-120-GSC	Chicken	Chicken Anti-GST IgG (total) ELISA Kit; Quantitative, Ag-based Ab Assay, 96 tests	
ELISA Kits/Animals/Bovine	Bovine Albumin (BSA) ELISA Kit, 96 tests, Quantitative (for measuring residual BSA in protein/antibody samples not recommended for measuring BSA in serum)	8100	Bovine	Stds=0-80 ng/ml ; Sample=100 ul; Sensitivity= 0.5 ng/ml; 135 min assay; Quantitative, sandwich assay	
ELISA Kits/Animals/Human, ELISA Kits Bioprocess contamination	Human Anti-KLH IgG (total) ELISA Kit, 96 tests, Quantitative	700-140-KLM	Human		
ELISA Kits/Animals/Goat, ELISA Kits Bioprocess contamination	Goat Anti-GST IgG (total) ELISA Kit, 96 tests, Quantitative	720-100-GSG	Goat	Goat Anti-GST IgG (total) ELISA Kit; Quantitative, Ag-based Ab Assay, 96 tests	
ELISA Kits/Animals/Hamster, ELISA Kits for the detection of Abs	Mouse Anti-Hepatitis B Surface Antigen (anti-HBsAg) ELISA kit, Semi-Quantitative	4210	Mouse	Anti-Hepatitis B Surface Antigen (anti-HBsAg) ELISA kit Human; Semi-Quantitative, Direct Antigen based assay	
ELISA Kits/Animals/Hamster, ELISA Kits for the detection of Abs	Ferret Serum Antibody detection ELISA kit, Qualitative	80188	Ferret	Ferret Serum Antibody detection ELISA kit; Qualitative, Direct Antigen based assay	
ELISA Kits/Animals/Rat, ELISA Kits for the detection of Abs	Rat Serum Antibody detection ELISA kit, Qualitative	80155	Rat	ELISA Kit for the detection of Rat Antibodies in serum, 75 min; Qualitative, Direct Antigen based assay	
ELISA Kits/Animals/Turkey, ELISA Kits for the detection of Abs	Turkey Serum Antibody detection ELISA kit, Qualitative	80187	Turkey	ELISA Kit for the detection of Turkey Antibodies in serum, 75 min; Qualitative, Direct Antigen based assay	
ELISA Kits/Animals/Sheep, ELISA Kits for the detection of Abs	Sheep Serum Antibody detection ELISA kit, Qualitative	80166	Sheep	ELISA Kit for the detection of Sheep Antibodies in serum, 75 min; Qualitative, Direct Antigen based assay	
ELISA Kits/Animals/Goat, ELISA Kits for the detection of Abs	Goat Serum Antibody detection ELISA kit, Qualitative	80165	Goat	ELISA Kit for the detection of Goat Antibodies in serum, 75 min; Qualitative, Direct Antigen based assay	
ELISA Kits/Animals/G, pig, ELISA Kits for the detection of Abs	G, Pig Serum Antibody detection ELISA kit, Qualitative	80161	G, Pig	ELISA Kit for the detection of G, Pig Antibodies in serum, 75 min; Qualitative, Direct Antigen based assay	
ELISA Kits/Animals/Hamster, ELISA Kits for the detection of Abs	Hamster Serum Antibody detection ELISA kit, Qualitative	80180	Hamster	ELISA Kit for the detection of Hamster Antibodies in serum, 75 min; Qualitative, Direct Antigen based assay	
ELISA Kits for the detection of Abs, ELISA Kits/Human	Human Serum Antibody detection ELISA kit, Qualitative	80170	Human	ELISA Kit for the detection of Human Antibodies in serum, 75 min; Qualitative, Direct Antigen based assay	
ELISA Kits/Animals/Monkey, ELISA Kits for the detection of Abs	Monkey Serum Antibody detection ELISA kit, Qualitative	80175	Monkey	ELISA Kit for the detection of Monkey Antibodies in serum, 75 min; Qualitative, Direct Antigen based assay	

ELISA Kits/Animals/Mouse, ELISA Kits for the detection of Abs	Mouse Serum Antibody detection ELISA kit, Qualitative	80150	Mouse	ELISA Kit for the detection of Mouse Antibodies in serum, 75 min, Qualitative, Direct Antigen based assay	
ELISA Kits/Animals/Pig, ELISA Kits for the detection of Abs	Pig Serum Antibody detection ELISA kit, Qualitative	80186	Pig	ELISA Kit for the detection of Pig Antibodies in serum, 75 min, Qualitative, Direct Antigen based assay	
Fusion proteins (GST, His, Myc, Beta-Galactosidase, MBP, VSV)	Mouse Anti-E. coli proteins Ig's ELISA kit, 96 tests	500-100-ECP	Mouse	Kit	
Fusion proteins (GST, His, Myc, Beta-Galactosidase, MBP, VSV)	E. coli proteins (5 strains) ELISA kit, 96 tests	800-130-ECP	E. coli	Kit	
Fusion proteins (GST, His, Myc, MBP, VSV)	SP20 Mouse Myeloma Cells (SP20) Proteins (culture medium) ELISA kit, 96 tests	800-160-SP2	Ver0 Cells	Kit	
Fusion proteins (GST, His, Myc, Beta-Galactosidase, MBP, VSV)	SP20 Mouse Myeloma Cells (SP20) Proteins (culture medium) ELISA kit, 96 tests	800-160-SP2	SP20 Cells	Kit	
Fusion proteins (GST, His, Myc, Beta-Galactosidase, MBP, VSV)	Chinese Hamster Ovary Cell (CHO) Proteins (culture medium) ELISA kit, 96 tests	800-140-CHO	CHO Cells	Kit	
Fusion proteins (GST, His, Myc, Beta-Galactosidase, MBP, VSV)	Rabbit Anti-E. coli proteins Ig's ELISA kit, 96 tests	500-120-ECP	Rabbit	Kit	
Anthrax Toxin/Receptors (ATR), Lethal factor/Protective Antigen/Spore	Anthrax Protective Antigen 83 (PA83) Protein ELISA kit	800-100-P83	Anthrax bacillus	Quantitative, Sandwich   Sds = 6.25 – 100 ng/ml   Sample=100 ul; 105 min assay   Sensitivity ~3ng/ml	
Anthrax Toxin/Receptors (ATR), Lethal factor/Protective Antigen/Spore	Mouse Anti-Anthrax Protective Antigen 83 (PA83) Ig's ELISA kit	900-100-83T	Anthrax bacillus	Indirect, Quantitative  Sds = 12.5 – 200 ng/ml   Sample=100 ul; 105 min assay	
Anthrax Toxin/Receptors (ATR), Lethal factor/Protective Antigen/Spore	Goat Anti-Anthrax Protective Antigen 83 (PA83) Ig's ELISA kit	900-105-83G	Anthrax bacillus	Indirect, Quantitative  Sds = 12.5 – 200 ng/ml   Sample=100 ul; 105 min assay	
Anthrax Toxin/Receptors (ATR), Lethal factor/Protective Antigen/Spore	Goat Anti-Anthrax Protective Antigen 83 (PA83) Ig's ELISA kit	900-105-83G	Anthrax bacillus	Indirect, Quantitative  Sds = 12.5 – 200 ng/ml   Sample=100 ul; 105 min assay	
Anthrax Toxin/Receptors (ATR), Lethal factor/Protective Antigen/Spore	Mouse Anti-Anthrax Protective Antigen 83 (PA83) Ig's WB kit, 12 tests	900-170-83WM	Mouse	PA83 recombinant proteins separated on SDS-PAGE and transferred on gel; Kit contains 12 PVDF strips--ve control and positive control rabbit sera	
Anthrax Toxin/Receptors (ATR), Lethal factor/Protective Antigen/Spore	Rabbit Anti-Anthrax Protective Antigen 83 (PA83) Ig's WB kit, 12 tests	900-180-83WR	Rabbit	PA83 recombinant proteins separated on SDS-PAGE and transferred on gel; Kit contains 12 PVDF strips--ve control and positive control rabbit sera	
Anthrax Toxin/Receptors (ATR), Lethal factor/Protective Antigen/Spore	Mouse Anti-lethal Factor (LF) Ig's ELISA kit	900-200-LFM	Mouse		
Anthrax Toxin/Receptors (ATR), Lethal factor/Protective Antigen/Spore	Rabbit Anti-Edema Factor (EF) Ig's ELISA kit	900-320-EFR	Rabbit		
Anthrax Toxin/Receptors (ATR), Lethal factor/Protective Antigen/Spore	Rabbit Anti-Edema Factor (EF) Ig's ELISA kit	900-320-EFR	Rabbit		
Anthrax Toxin/Receptors (ATR), Lethal factor/Protective Antigen/Spore	Mouse Anti-Edema Factor (EF) Ig's ELISA kit	900-300-EFR	Mouse		
Anthrax Toxin/Receptors (ATR), Lethal factor/Protective Antigen/Spore	Rabbit Anti-Anthrax Protective Antigen 83 (PA83) Ig's ELISA kit	900-120-83T	Anthrax bacillus	Indirect, Quantitative  Sds = 12.5 – 200 ng/ml   Sample=100 ul; 105 min assay	
ELISA Kits/Bacterial/Viral/Vaccines	Swine/Pig Anti-Influenza A virus IgM ELISA kit	920-020-PAM	Swine/Pig		
ELISA Kits/Bacterial/Viral/Vaccines	Chicken Anti-Avian Influenza A virus (AIV) IgM ELISA kit	920-105-AIM	Chicken		
ELISA Kits/Bacterial/Viral/Vaccines	Swine/Pig Anti-Influenza A virus IgA ELISA kit	920-030-PAA	Swine/Pig		
ELISA Kits/Bacterial/Viral/Vaccines	Human Anti-Influenza A virus IgG ELISA kit	920-040-HAG	Swine/Pig		
ELISA Kits/Bacterial/Viral/Vaccines	Human Anti-Influenza A virus IgM ELISA kit	920-050-HAM	Swine/Pig		
ELISA Kits/Bacterial/Viral/Vaccines	Human Anti-Influenza A virus IgG ELISA kit	920-060-HAA	Swine/Pig		
ELISA Kits/Bacterial/Viral/Vaccines	Swine/Pig Anti-Influenza A virus IgG ELISA kit	920-010-PAG	Swine/Pig		
ELISA Kits/Bacterial/Viral/Vaccines	Chicken serum (SPP) tested –ve for various chicken viruses, Salmonella and Mycoplasma (suitable for chicken virus ELISA kits)	CSNC11-S	Chicken		
ELISA Kits/Bacterial/Viral/Vaccines	Chicken Allantoic fluid (SPP eggs) tested –ve for various chicken viruses, Salmonella and Mycoplasma (suitable for chicken virus ELISA kits)	CAF11-S	Chicken		
ELISA Kits/Bacterial/Viral/Vaccines	Chicken Anti-Marek's Disease Virus (MDV) IgG ELISA kit	920-130-MDV	Chicken		
ELISA Kits/Bacterial/Viral/Vaccines	Chicken Anti-Infectious Bronchitis Virus (IBV) IgG ELISA kit	920-120-IBV	Chicken		
ELISA Kits/Bacterial/Viral/Vaccines	Chicken Anti-Anemia Virus (AV) IgG ELISA kit	920-110-AV	Chicken		
ELISA Kits/Bacterial/Viral/Vaccines	Chicken Anti-Avian Influenza A virus (AIV) IgG ELISA kit	920-300-H5I	Chicken		
ELISA Kits/Bacterial/Viral/Vaccines	Chicken Anti-Avian Influenza virus (H5N1) IgG ELISA kit (1x96 wells)	920-300-H5I	Chicken		
ELISA Kits/Bacterial/Viral/Vaccines	Chicken Anti-Avian Influenza virus (H5N1) IgG ELISA kit (2x96 wells)	920-300-H52	Chicken		
ELISA Kits/Bacterial/Viral/Vaccines	Chicken Anti-Newcastle Disease Virus (NDV) IgGs ELISA kit	920-140-NDV	Chicken		
ELISA Kits/Bacterial/Viral/Vaccines	Mouse Anti-Tetanus Toxoid IgG ELISA kit	950-110-TTM	Mouse		
ELISA Kits/Bacterial/Viral/Vaccines	Mouse Anti-Adenovirus IgM ELISA kit	950-140-AMM	Mouse		
ELISA Kits/Bacterial/Viral/Vaccines	Mouse Anti-Adenovirus IgG ELISA kit	950-130-AMG	Mouse		
ELISA Kits/Bacterial/Viral/Vaccines	Human Anti-Adenovirus IgM ELISA kit	950-120-AHM	Human		
ELISA Kits/Bacterial/Viral/Vaccines	Human Anti-Adenovirus IgG ELISA kit	950-110-AHG	Human		
ELISA Kits/Bacterial/Viral/Vaccines	Human Anti-Adenovirus IgA ELISA kit	950-100-AHA	Human		
ELISA Kits/Bacterial/Viral/Vaccines	Human Anti-Diphtheria Toxoid IgM ELISA kit	940-110-DHM	Human		
ELISA Kits/Bacterial/Viral/Vaccines	Human Anti-B. pertussis IgA ELISA kit	960-100-PHA	Human		
ELISA Kits/Bacterial/Viral/Vaccines	Human Anti-Diphtheria Toxoid IgG ELISA kit	940-100-DHG	Human		
ELISA Kits/Bacterial/Viral/Vaccines	Mouse Anti-Japanese encephalitis virus (JEV) Ig's ELISA kit	910-100-JEM	Mouse		
ELISA Kits/Bacterial/Viral/Vaccines	Mouse Anti-Japanese encephalitis virus (JEV) Ig's WB kit, 12 tests	910-110-JWM	Mouse		
ELISA Kits/Bacterial/Viral/Vaccines	Human Anti-Diphtheria Toxoid IgM ELISA kit	940-120-DHM	Human		
ELISA Kits/Bacterial/Viral/Vaccines	Human Anti-B. pertussis IgM ELISA kit	960-120-PHM	Human		
ELISA Kits/Bacterial/Viral/Vaccines	Mouse Anti-B. pertussis IgG ELISA kit	960-130-PHG	Mouse		
ELISA Kits/Bacterial/Viral/Vaccines	Human Anti-B. pertussis IgG ELISA kit	960-110-PHG	Human		
ELISA Kits/Bacterial/Viral/Vaccines	Mouse Anti-B. pertussis IgM ELISA kit	960-140-PHM	Mouse		
ELISA Kits/Bacterial/Viral/Vaccines	Human Anti-Tetanus Toxoid IgG ELISA kit	930-100-TTH	Human		
ELISA kits for Animal Diseases (bacterial/Viral)	Porcine Circovirus ELISA kit, 192 tests Control Disease (Serum)	AE-200150-2	Pig/Porcine/Swine	100 min Serum; ; 100 min	
ELISA kits for Animal Diseases (bacterial/Viral)	Bird Flu ELISA kit, 192 tests, Detect H5N1 Antibody (Serum)	AE-200190-2	Bird	100 min Serum; ; 100 min	
ELISA kits for Animal Diseases (bacterial/Viral)	PRRSV Antibody ELISA kit, 192 tests, Control Disease (Serum)	AE-200180-2	Pig/Porcine/Swine	100 min Serum; ; 100 min	
ELISA kits for Animal Diseases (bacterial/Viral)	Porcine Epidemic Encephalitis B Ab ELISA kit, 192 tests Control Disease (Serum)	AE-200160-2	Pig/Porcine/Swine	100 min Serum; ; 100 min	
ELISA kits for Animal Diseases (bacterial/Viral)	Hog Cholera Virus Antibody ELISA kit, 192 tests Control Disease	AE-200140-2	Pig/Cholera Virus	100 min Serum; ; 100 min	
ELISA kits for Animal Diseases (bacterial/Viral)	Porcine Pseudorabies Virus IgE Antibody Distinguishing kit, 192 tests Detect Pig Immunized with Inactivated Vaccine from Pig	AE-200130-2	Pig/Porcine/Swine	100 min 100 min Serum; ; 100 min	
ELISA kits for Animal Diseases (bacterial/Viral)	Porcine Pseudorabies Antibody ELISA kit, 192 Control Disease (Serum)	AE-200130-2	Pig/Porcine/Swine	100 min Serum; ; 100 min	
ELISA kits for Animal Diseases (bacterial/Viral)	Swine Foot and Mouth Disease ELISA kit, 192 tests Detect Ab (Serum)	AE-200120-2	Pig/Porcine/Swine	100 min Serum; ; 100 min	
ELISA kits for Animal Diseases (bacterial/Viral)	Swine Foot and Mouth Disease IgG Distinguishing kit, 192 tests Detect pig Immunized with Inactivated Vaccine Pig, (Serum)	AE-200120-2	Pig/Porcine/Swine	100 min 100 min Serum; ; 100 min	
ELISA kits for Animal Diseases (bacterial/Viral)	Porcine Toxoplasmosis IgG Antibody kit, 96 tests Detect Ab(Serum)	AE-200110-2	Pig/Porcine/Swine	100 min Serum; ; 100 min	
ELISA kits for Animal Diseases (bacterial/Viral)	Porcine Parvovirus Antibody ELISA kit, 192 tests Monitor Disease (Serum)	AE-200170-2	Pig/Porcine/Swine	100 min Serum; ; 100 min	
ELISA kits-Food Safety & Drug Residue	Nitrofuran (AOZ) ELISA kit, 96 tests (For Fish, Shrimp, Honey, Chicken/Liver)	DE-100080	All Species	Fish, Shrimp: 0.2; 0.16; 80 min Honey: 0.1; 16; 80 min ; ; 120 min Chicken/Liver: 0.1; 16; h;	
ELISA kits-Food Safety & Drug Residue	Clenbuterol ELISA kit, 96 tests (For Urine, Serum, Feed, Meat, Chicken)	DE-100010	ELISA Kits	Pork/Liver: 0.1; 50 min; 80 min; Serum: 0.1; 10 min; 80 min; Urine: 0.1; 10 min; 80 min; Feed: 50; 30 min; 80 min; Meat: 0.2; 20 min; 80 min;	
ELISA kits-Food Safety & Drug Residue	Ractopamine ELISA kit, 96 tests (For Liver, Urine, Feed)	DE-100020	All Species	Liver: 0.4; 20 min; 80 min Urine: 0.2; 10 min; 80 min Feed: 15; 30 min; 50 min	
ELISA kits-Food Safety & Drug Residue	Salbutamol ELISA kit, 96 tests (For Urine, Tissue, Feed, Animal Tissue, Aquatic, Honey, Intestine, Chloramphenicol ELISA kit, 96 tests (For Animal Tissue, Aquatic, Honey, Intestine, Urine, Egg, Milk, Serum)	DE-100030	All Species	Urine: 0.5; 10 min; 80 min Tissue: 0.5; 20 min; 80 min Feed: 50; 30 min; 80 min	
ELISA kits-Food Safety & Drug Residue	Flortencol ELISA kit, 96 tests (For Animal Tissue, Aquatic, Honey)	DE-100040	All Species	Animal tissue: 0.05; 20 min; 80 min Aquatic: 0.05; 20 min; 80 min Honey: 0.15; 20 min; 0 min Intestine: 0.1; 20 min; 80 min Urine: 0.1; 2.5h; 80 min Egg: 0.1; 30 min; 80 min Milk: 0.05; 30 min; 80 min Serum: 0.1; 15 min; 80 min	
ELISA kits-Food Safety & Drug Residue	Tylosin ELISA kit, 96 tests (For Meat, Liver, Honey, Egg)	DE-100160	All Species	Animal tissue: 0.5; 30 min; 80 min Aquatic: 0.5; 30 min; 80 min Honey: 1.5; 30 min; 80 min	
ELISA kits-Food Safety & Drug Residue	Sulfagin-oxaline ELISA kit, 96 tests (For Pork/Liver, Honey/Egg, Serum/Urine, Milk)	DE-100240	All Species	Meat: 3; 20 min; 60 min Liver: 3; 20 min; 60 min Honey/Egg: 1.5; 10 min; 60 min	
ELISA kits-Food Safety & Drug Residue	Olaquinox ELISA kit, 96 tests (Tissue)	DE-100230	All Species	Pork/Liver: 1; 20 min; 90 min Honey/Egg: 1; 10 min; 90 min Serum/Urine: 4; 10 min; 90 min Milk; 20; 10 min; 90 min	
ELISA kits-Food Safety & Drug Residue	Streptomycin ELISA kit, 96 tests (Chicken/Liver, Honey, Milk)	DE-100210	All Species	Tissue: 1; 40 min; 80 min 80 min	
ELISA kits-Food Safety & Drug Residue	Gentamicin ELISA kit, 96 tests (Chicken/Liver)	DE-100200	All Species	Chicken/Liver: 20; 40 min; 60 min Honey: 10; 40 min; 60 min Milk; 20; 5 min; 60 min	
ELISA kits-Food Safety & Drug Residue	Quinolones (QNS) ELISA kit, 96 tests (For Pork/Liver, Chicken/Liver, Shrimp, Fish, Serum, Honey)	DE-100120	All Species	80 min Chicken/Liver: 4; 40 min; 80 min	
ELISA kits-Food Safety & Drug Residue	Diethylstilbestrol (DES) ELISA kit, 96 tests (Fish, Shrimp, Liver, Meat, Feed, Urine)	DE-100190	All Species	Pork/Liver: 3; 50 min; 90 min Chicken/Liver: 3; 50 min; 90 min Shrimp, Fish: 3; 20 min; 90 min Serum: 3; 10 min; 90 min Honey: 6; 10 min; 80 min	
ELISA kits-Food Safety & Drug Residue	Nitrofuran (SEM) ELISA kit, 96 tests (Honey, Fish, Shrimp, Chicken/Liver, Fish/Shrimp)	DE-100070	All Species	Fish, Shrimp: 0.2; 20 min; 80 min Liver, Meat: 2; 20 min; 80 min Feed: 20; 20 min; 80 min Urine: 0.6; 20 min; 80 min	
ELISA kits-Food Safety & Drug Residue	Trenbolone ELISA kit, 96 tests (For Animal Tissue, Aquatic, Urine)	DE-100170	All Species	Honey: 0.2; 16; 80 min Fish, Shrimp: 0.3; 16; 80 min ; ; ; 80 min Chicken/Liver: 0.2; 16; 80 min	
ELISA kits-Food Safety & Drug Residue	Nitrofuran (AMCZ) ELISA kit, 96 tests (For Fish, Shrimp, Honey, Chicken/Liver)	DE-100060	All Species	Animal tissue: 5; 30 min; 80 min Aquatic: 5; 30 min; 80 min Urine: 0.5; 30 min; 80 min	
ELISA kits-Food Safety & Drug Residue	Sulfonamides Residues (SAs) ELISA kit, 96 tests (For Chicken/Liver, Pork/Liver, Honey/Egg, Serum/Urine, Milk)	DE-100090	All Species	Fish, Shrimp: 0.3; 2 h; 80 min Honey: 0.2; 2 h; 80 min Chicken/Liver: 0.2; 2 h; 80 min	
ELISA kits-Food Safety & Drug Residue	Sulfathiazine (SMZ) ELISA kit, 96 tests (For Chicken/Liver, Pork/Liver, Honey/Egg, Serum/Urine, Milk)	DE-100100	All Species	Chicken/Liver: 1; 20 min; 90 min Pork/Liver: 1; 20 min; 90 min Honey/Egg: 1; 10 min; 90 min Serum/Urine: 4; 10 min; 90 min Milk; 20; 10 min; 90 min	
ELISA kits-Food Safety & Drug Residue	Sulfamethoxydiazine (SMD) ELISA kit, 96 tests (For Chicken/Liver, Pork/Liver, Honey/Egg, Serum/Urine)	DE-100110	All Species	Chicken/Liver: 1; 20 min; 90 min Pork/Liver: 1; 20 min; 90 min Honey/Egg: 1; 10min; 90 min Serum/Urine: 4; 10 min; 90 min	
ELISA kits-Food Safety & Drug Residue	Benzyl Penicillin ELISA kit, 96 tests (For Pork/Liver, Chicken, Duck, Shrimp, Fish, Honey, Milk)	DE-100150	All Species	Pork/Liver: 2; 15 min; 80 min Chicken(meat,liver): 2; 15 min; 80 min Duck(meat,liver): 2; 15 min; 80 min Shrimp, Fish: 2; 15 min; 80 min Honey/Milk; 2; 10 min; 80 min	
ELISA kits-Food Safety & Drug Residue	Ampicillin ELISA kit, 96 tests (For Pork/Liver, Chicken, Duck, Shrimp, Fish, Honey, Milk)	DE-100140	All Species	Pork/Liver: 2; 15 min; 80 min 80 min 80 min 80 min	
ELISA kits-Food Safety & Drug Residue	Enrofloxacin ELISA kit, 96 tests (For Pork/Liver, Chicken/Liver, Shrimp, Fish, Serum, Honey)	DE-100130	All Species	Pork/Liver: 1; 50 min; 80 min 80 min Chicken Liver: 1; 50 min; 80 min Shrimp, Fish: 1; 20 min; 80 min Serum: 1; 10 min; 80 min Honey: 2; 10 min; 80 min	
ELISA kits-Food Safety & Drug Residue	Nitrofuran (AHJ) ELISA kit, 96 tests (For Fish, Shrimp, Honey, Chicken/Liver)	DE-100170	All Species	Fish, Shrimp: 0.3; 16; 80 min Honey: 0.2; 16; 80 min Chicken/Liver: 0.2; 16; 80 min	
ELISA kits-Food Safety & Drug Residue	Diazepam ELISA kit, 96 tests (For Tissue, Urine, Feed)	DE-100180	All Species	Tissue: 1; 30 min; 80 min Urine: 1; 30 min; 80 min Feed: 20; 30 min; 80 min	

Contact Us: Life Technologies (India) Pvt. Ltd., 306, Agarwal City Mall, opposite M2K Pitampura Delhi-110034 (India)

Tel # +91-11-4220-8000/4220-8111/4220-8222 Fax# +91-11-4220-8444/Mobile# +91-98105-21400, Email - customerservice@lifetechnia.com or www.lifetechnia.com