



LIFE TECHNOLOGIES (INDIA) PVT. LTD.

your molecular & cell technology partner

Bringing zest to your research



BIO BASIC INC

Plastic Ware

CODE	DESCRIPTION	QTY
1	BK044-2 GRADUATED PLASTIC BEAKERS, 50ML	1
2	BK044-3 GRADUATED PLASTIC BEAKERS, 100ML	1
3	BK044-4 GRADUATED PLASTIC BEAKERS, 250ML	1
4	BK044-5 GRADUATED PLASTIC BEAKERS, 500ML	1
5	BK044-1K GRADUATED PLASTIC BEAKERS, 1000ML	1
6	BK044-2K GRADUATED PLASTIC BEAKERS, 2000ML	1
7	B003N-15 PLASTIC BOTTLE, WIDE NECK, 15ML, NATURAL	10
8	B003N-S15 PLASTIC BOTTLE, WIDE NECK, STERILE, 15ML, NATURAL	10
9	B003N-30 PLASTIC BOTTLE, WIDE NECK, 30ML, NATURAL	10
10	B003N-S30 PLASTIC BOTTLE, WIDE NECK, STERILE, 30ML, NATURAL	10
11	B003N-70 PLASTIC BOTTLE, WIDE NECK, 70ML, NATURAL	10
12	B003N-S70 PLASTIC BOTTLE, WIDE NECK, STERILE, 70ML, NATURAL	10
13	B003N-125 PLASTIC BOTTLE, WIDE NECK, 125ML, NATURAL	10
14	B003N-S125 PLASTIC BOTTLE, WIDE NECK, STERILE, 125ML, NATURAL	10
15	B003N-250 PLASTIC BOTTLE, WIDE NECK, 250ML, NATURAL	10
16	B003N-S250 PLASTIC BOTTLE, WIDE NECK, STERILE, 250ML, NATURAL	10
17	B003N-500 PLASTIC BOTTLE, WIDE NECK, 500ML, NATURAL	5
18	B003N-S500 PLASTIC BOTTLE, WIDE NECK, STERILE, 500ML, NATURAL	5
19	B010-250 PLASTIC BOTTLE, SMALL NECK, 250ML, NATURAL	1
20	B011-500 PLASTIC BOTTLE, SMALL NECK, 500ML, NATURAL	1
21	BV040 GLASS BOTTLE, 250ML, BLUE CAP	1
22	BV041 GLASS BOTTLE, 500ML, BLUE CAP	1
23	XP001 WASHING PLASTIC BOTTLE, NATURAL	1
24	XP002 WASHING PLASTIC BOTTLE, RED	1
25	XP003 WASHING PLASTIC BOTTLE, BLUE	1
26	XP004 WASHING PLASTIC BOTTLE, YELLOW	1
27	BR001-10 BURETTES BRUSH, 10X100MM	1
28	BR001-20 BURETTES BRUSH, 20X100MM	1
29	BR001-30 BURETTES BRUSH, 30X100MM	1
30	BR004-10 BURETTES BRUSH, 10ML	1
31	BR004-20 BURETTES BRUSH, 20ML	1
32	BR004-50 BURETTES BRUSH, 50ML	1
33	BR005-100 VOLUMETRIC FLASK-BRUSHES, 100ML	1
34	BR005-250 VOLUMETRIC FLASK-BRUSHES, 250ML	1
35	BR005-500 VOLUMETRIC FLASK-BRUSHES, 500ML	1
36	BR006-100 FLASK-BRUSHES, 100ML	1
37	BR006-250 FLASK-BRUSHES, 250ML	1
38	BR006-500 FLASK-BRUSHES, 500ML	1
39	BR007-250 BEAKER BRUSHES, 250ML	1
40	BR007-500 BEAKER BRUSHES, 500ML	1
41	BR007-1K BEAKER BRUSHES, 1L	1
42	BR008-250 CYLINDER BRUSH, 250ML	1
43	BR008-500 CYLINDER BRUSH, 500ML	1
44	BR008-1K CYLINDER BRUSH, 1000ML	1
45	CH001 10/200mm Communis Chromatography Columns	1
46	CH002 10/300mm Communis Chromatography Columns	1

Contact Us: **Life Technologies (India) Pvt. Ltd .**, 306, Agarwal City Mall, opposite M2K Pitampura ,
Delhi-110034 (India) Tel # +91-11-4220-8000;4220-8111;4220-8222 Fax# +91-11-4220-8444,Mobile#
+91-98105-21400, Email - customerservice@lifetechindia.com Web: www.lifetechindia.com



47	CH003	10/400mm Communis Chromatography Columns	1
48	CH004	10/500mm Communis Chromatography Columns	1
49	CH007	16/200mm Communis Chromatography Columns	1
50	CH008	16/300mm Communis Chromatography Columns	1
51	CH009	16/400mm Communis Chromatography Columns	1
52	CH010	16/500mm Communis Chromatography Columns	1
53	CH013	26/200mm Communis Chromatography Columns	1
54	CH014	26/300mm Communis Chromatography Columns	1
55	CH015	26/400mm Communis Chromatography Columns	1
56	CH016	26/500mm Communis Chromatography Columns	1
57	CH022	35/200mm Communis Chromatography Columns	1
58	CH023	35/300mm Communis Chromatography Columns	1
59	CH024	35/400mm Communis Chromatography Columns	1
60	CH025	35/500mm Communis Chromatography Columns	1
61	CH028	45/200mm Communis Chromatography Columns	1
62	CH029	45/300mm Communis Chromatography Columns	1
63	CH030	45/400mm Communis Chromatography Columns	1
64	CH031	45/500mm Communis Chromatography Columns	1
65	CH034	55/200mm Communis Chromatography Columns	1
66	CH035	55/300mm Communis Chromatography Columns	1
67	CH036	55/400mm Communis Chromatography Columns	1
68	CH037	1055/500mm Communis Chromatography Columns	1
69	CX011	Cellulose Paper, 1mm thick, 60x60cm, 25/pk	1
70	CX012	Cellulose Paper, 3mm thick, 60x60cm, 25/pk	1
71	CT331-50	50G Cotton Degrease Cotton & Adhesive Plaster	1
72	CT331-100	100G Cotton Degrease Cotton & Adhesive Plaster	1
73	CT331-500	500G Cotton Degrease Cotton & Adhesive Plaster	1
74	CD001-1	GRADUATED PLASTIC CYLINDERS, 100ML	1
75	CD001-2	GRADUATED PLASTIC CYLINDERS, 250ML	1
76	CD001-3	GRADUATED PLASTIC CYLINDERS, 500ML	1
77	CD001-4	GRADUATED PLASTIC CYLINDERS, 1000ML	1
78	CD001-5	GRADUATED PLASTIC CYLINDERS, 2000ML	1
79	DS001-40	40x9mm Cell Culture Plastic Dishes	20
80	DS001-60	60X12mm Cell Culture Plastic Dishes	20
81	DS001-90	90X15mm Cell Culture Plastic Dishes	10
82	DS001-100	100X15mm Square Shaped Plastic Dishes for cell culture	10
83	DSG110	DRYING REAGENT, SILICA GEL	1KG
84	DSG111	SILICA GEL IN BAG, 5G, 100/PK	1PK
85	DB031	DRYING BOX, 100 X 40 X 10MM	1
86	PF002	Parafilm, 4x125'	1
87	PF003	Parafilm, 2x250'	1
88	HP131	Cellulose Acetate Membrane,13MM,0.45UM, 100/BOX	1
89	HP132	Cellulose Acetate Membrane,25MM,0.45UM, 100/BOX	1
90	HP133	Cellulose Acetate Membrane, 50MM,0.22UM, 100/BOX	1
91	HP134	Cellulose Acetate Membrane, 50MM,0.45UM, 100/BOX	1
92	HP135	Cellulose Acetate Membrane, 25MM,0.80UM, 100/BOX	1
93	HP136	NYLON MEMBRANCE, 13MM,0.45UM, 100/BOX	1
94	HP137	NYLON MEMBRANCE, 25MM,0.45UM, 100/BOX	1
95	HP138	NYLON MEMBRANCE, 50MM,0.45UM, 100/BOX	1
96	HP139	NYLON MEMBRANCE, 50MM,0.22UM, 100/BOX	1
97	HP141	Syringe Filtration - aqueous solutions, non sterile, 0.22UM, 13MM, 10/PK	1
98	HP142	Syringe Filtration - aqueous solutions, non sterile, 0.22UM, 25MM, 10/PK	1
99	HP143	Syringe Filtration - aqueous solutions, non sterile, 0.45UM, 13MM, 10/PK	1
100	HP144	Syringe Filtration - aqueous solutions, non sterile, 0.45UM, 25MM, 10/PK	1
101	HP145	Syringe Filtration - organic solutions, non sterile, 0.22UM, 13MM, 10/PK	1
102	HP146	Syringe Filtration - organic solutions, non sterile, 0.22UM, 25MM, 10/PK	1
103	HP147	Syringe Filtration - organic solutions, non sterile, 0.45UM, 13MM, 10/PK	1
104	HP148	Syringe Filtration - organic solutions, non sterile, 0.45UM, 25MM, 10/PK	1



105	HP151	Syringe Filtration - aqueous solutions, sterile, 0.22UM, 13MM, STERILE, 10/PK	1
106	HP152	Syringe Filtration - aqueous solutions, sterile, 0.22UM, 25MM, STERILE, 10/PK	1
107	HP153	Syringe Filtration - aqueous solutions, sterile, 0.45UM, 13MM, STERILE, 10/PK	1
108	HP154	Syringe Filtration - aqueous solutions, sterile, 0.45UM, 25MM, STERILE, 10/PK	1
109	HP155	Syringe Filtration - organic solutions, sterile, 0.22UM, 13MM, STERILE, 10/PK	1
110	HP156	Syringe Filtration - organic solutions, sterile, 0.22UM, 25MM, STERILE, 10/PK	1
111	HP157	Syringe Filtration - organic solutions, sterile, 0.45UM, 13MM, STERILE, 10/PK	1
112	HP158	Syringe Filtration - organic solutions, sterile, 0.45UM, 25MM, STERILE, 10/PK	1
113	FC001	FORCEPS, 12.5CM, SHARP END	1
114	FC002	FORCEPS, 16CM, SHARP END	1
115	FC003	FORCEPS, 12.5CM, BLUNT END	1
116	FC004	FORCEPS, 16CM, BLUNT END	1
117	FC005	FORCEPS, 20CM, BLUNT END	1
118	FN001	REGULAR PLASTIC FUNNELS, 60MM	1
119	FN002	REGULAR PLASTIC FUNNELS, 68MM	1
120	FN003	REGULAR PLASTIC FUNNELS, 100MM	1
121	FN004	REGULAR PLASTIC FUNNELS, 150MM	1
122	FN005	REGULAR PLASTIC FUNNELS, 250MM	1
123	GL001S	DISPOSABLE LATEX GLOVES, SMALL, 100/BOX	1PK
124	GL002M	DISPOSABLE LATEX GLOVES, MEDIUM, 100/BOX	1PK
125	GL003L	DISPOSABLE LATEX GLOVES, LARGE, 100/BOX	1PK
126	GL003XL	DISPOSABLE LATEX GLOVES, XLARGE, 100/BOX	1PK
127	GL111BS	NITRILE POWDER-FREE GLOVES, SMALL, BLUE, 100/BOX	1PK
128	GL111BM	NITRILE POWDER-FREE GLOVES, MEDIUM, BLUE, 100/BOX	1PK
129	GL111BL	NITRILE POWDER-FREE GLOVES, LARGE, BLUE, 100/BOX	1PK
130	GL111BXL	NITRILE POWDER-FREE GLOVES, EXT LARGE, BLUE, 100/BOX	1PK
131	GL112GS	NITRILE POWDER-FREE GLOVES, SMALL, GREEN, 100/BOX	1PK
132	GL112GM	NITRILE POWDER-FREE GLOVES, MEDIUM, BLUE, 100/BOX	1PK
133	GL112GL	NITRILE POWDER-FREE GLOVES, LARGE, BLUE, 100/BOX	1PK
134	GL112GXL	NITRILE POWDER-FREE GLOVES, EXT LARGE, BLUE, 100/BOX	1PK
135	GL011S	Disposable Polyethylene Gloves PE, SMALL, 100/BOX	1PK
136	GL011M	Disposable Polyethylene Gloves PE, MEDIUM, 100/BOX	1PK
137	GL011L	Disposable Polyethylene Gloves PE, LARGE, 100/BOX	1PK
138	GL331	CHEM-MASTER GLOVES, 30CM	1PAIRS
139	GL332	CHEM-MASTER GLOVES, 35CM	1PAIRS
140	GL333	CHEM-MASTER GLOVES, 40CM	1PAIRS
141	GL334	CHEM-MASTER GLOVES, 45CM	1PAIRS
142	GL335	CHEM-MASTER GLOVES, 50CM	1PAIRS
143	GL336	CHEM-MASTER GLOVES, 55CM	1PAIRS
144	IB001	ICE BUCKET	1
145	IT771	ICE TRAYS, 1.5L	1
146	IT772	ICE TRAYS, 3.0L	1
147	IL311-1	INOCULATING LOOPS&NEEDLES, 1UL BLUE, STERILE, 25/BAG	1PK
148	IL311-10	INOCULATING LOOPS&NEEDLES, 10UL, YELLOW, STERILE, 25/BAG	1PK
149	IL311-N	INOCULATING LOOPS&NEEDLES, NEEDLE WHITE, STERILE, 25/BAG	1PK
150	IL771	INOCULATION NEEDLE HOLDER, 220MM, 2/PACK	1PK
151	IL772	INOCULATION NEEDLE WIRE, 70MM, 10/PACK	1PK
152	SP114	STAINLESS STEEL SPREADING, 16MM	1
153	SP115	STAINLESS STEEL SPREADING, 25MM	1
154	LB001G	LABELS ON ROLL, GREEN, 30X15MM	1ROLLS
155	LB002R	LABELS ON ROLL, RED, 30X15MM	1ROLLS
156	LB003W	LABELS ON ROLL, WHITE, 30X15MM	1ROLLS
157	LB004Y	LABELS ON ROLL, YELLOW, 30X15MM	1ROLLS
158	LB101G	LABELS ON ROLL, GREEN, 40X20MM	1ROLLS
159	LB102R	LABELS ON ROLL, RED, 40X20MM	1ROLLS
160	LB103W	LABELS ON ROLL, WHITE, 40X20MM	1ROLLS
161	LB104Y	LABELS ON ROLL, YELLOW, 40X20MM	1ROLLS
162	TP001B	LABELING TAPE, BLUE	1ROLLS



LIFE TECHNOLOGIES (INDIA) PVT. LTD.

your molecular & cell technology partner

Bringing zest to your research



BIO BASIC INC

163	TP002R	LABELING TAPE,RED	1ROLLS
164	TP003W	LABELING TAPE,WHITE	1ROLLS
165	TP004Y	LABELING TAPE,YELLOW	1ROLLS
166	LB021B	LABELING COLOR DOTS, BLUE/10SHEETS	1
167	LB022R	LABELING COLOR DOTS, RED/10SHEETS	1
168	LB023G	LABELING COLOR DOTS, GREEN/10SHEETS	1
169	LB024Y	LABELING COLOR DOTS, ASSORTED COLORS/10SHEETS	1
170	GR111	TISSUE GRINDING PESTLES, CONICAL	1
171	GR112	TISSUE GRINDING PESTLES, ROUND	1
172	MGR113	DISPOSABLE PLASTIC PESTLES, CONICAL, 210MM,50/BAG	1PK
173	MGR114	DISPOSABLE PLASTIC PESTLES, ROUND, 205MM,50/BAG	1PK
174	MGR115	DISPOSABLE PLASTIC PESTLES, CONICAL, 70MM,50/BAG	1PK
175	PP3322	Blotting Paper, 10x15cm, 100 SHEETS	1PK
176	PP3323	Blotting Paper, 15x15cm, 100 SHEETS	1PK
177	PP3324	Blotting Paper, 20x20cm, 100 SHEETS	1PK
178	RB333	MULTI-CHANNEL REAGENT BOAT	1
179	PF001-30	PIPETTE FILLER, 30ML	1
180	PF002-60	PIPETTE FILLER, 60ML	1
181	PF003-90	PIPETTE FILLER, 90ML	1
182	PF330	QUANTITATIVE PIPETTES FILLERS	1
183	BT215	0.1-10UL MICRO TIPS, SHORT, 1000/PACK	1PK
184	S2817	0.1-10UL MICRO TIPS, SHORT, 1000/PACK	1PK
185	BT215-RS	0.1-10UL MICRO TIPS, SHORT, 96/RACK, STERILE	1PK
186	BT216	0.1-10UL MICRO TIPS, MEDI, 1000/PACK	1PK
187	BT216-RS	0.1-10UL MICRO TIPS, MEDI, 96/RACK, STERILE	1PK
188	BT225	0.5-10UL MICRO TIPS, LONG, 1000/PACK	1PK
189	S2811	0.5-10UL MICRO TIPS, LONG, 1000/PACK	1PK
190	BT225-RS	0.5-10UL MICRO TIPS, 96/RACK, STERILE	1PK
191	S3332-R	0.5-10UL MICRO TIPS,96/RACK	1PK
192	S3332-RS	0.5-10UL MICRO TIPS,STERILE,96/RACK	1PK
193	BT224-Y	200UL YELLOW TIP, WITH GRADUATION, 1000/PACK	1PK
194	S2819	200UL NATURAL TIP, WITH GRADUATION, 1000/PACK	1PK
195	S2820	200UL YELLOW TIP, WITH GRADUATION, 1000/PACK	1PK
196	BT224-RNS	200UL YELLOW TIP, WITH GRADUATION, STERILE, 96/RACK	1PK
197	S2163-RS	200UL NATURAL TIP, WITH GRADUATION, STERILE, 100/RACK	1PK
198	S2177-RS	200UL YELLOW TIP, WITH GRADUATION, STERILE, 100/RACK	1PK
199	BT224-5RNS	200UL YELLOW TIP, NATURAL-5-LAYER-RACK, WITH GRADUATION, STERILE, 480/RACK	1PK
200	BT227-Y	200UL YELLOW TIP, WITHOUT GRADUATION, 1000/PACK	1PK
201	BT227-RNS	200UL YELLOW TIP, WITHOUT GRADUATION, STERILE, 96/RACK	1PK
202	BT330-N	300UL NATURAL TIPS,1000/PACK	1PK
203	BT330-RNS	300UL NATURAL TIPS, STERILE, 96/PACK	1PK
204	BT222	100-1000UL BLUE TIP, 1000/PACK	1PK
205	S2828	100-1000UL BLUE TIP, 1000/PACK	1PK
206	BT222-RS	100-1000UL BLUE TIP, NATURAL, STERILE, 100/RACK	1PK
207	S2119-R	100-1000UL NATURAL TIP, 100/RACK	1PK
208	S2124-RS	100-1000UL NATURAL TIP, STERILE, 100/RACK	1PK
209	BT999	5.0ML TIP, NATURAL, 250/PACK	1PK
210	BT999-SN	5.0ML TIP, NATURAL, STERILE, 50/PACK	1PK
211	Q214	0.1-10UL LOADING TIP ROUND 0.57 OD 1000/BAG	1PK
212	Q214-R204	0.1-10UL LOADING TIP ROUND 0.57 OD, 204/RACK	1PK
213	FT215	0.1-10UL FILTER TIP, NATURAL, 1000/PACK	1PK
214	FT215-RNS	0.1-10UL FILTER TIP, NATURAL, STERILE, 96/RACK	1PK
215	S2619-RS	0.1-10UL FILTER TIP, NATURAL, STERILE, 96/RACK, 10R/CS	10PK
216	S2532-RS	5-30UL FILTER TIP, NATURAL, STERILE, 96/RACK, 10R/CS	10PK
217	FT216	0.1-10UL FILTER TIP, NATURAL, 1000/PACK	1PK
218	FT216-RNS	0.1-10UL FILTER TIP, NATURAL, STERILE, 96/RACK	1PK
219	FT225	0.1-10UL MICRO TIP, NATURAL, 1000/PACK	1PK

Contact Us: **Life Technologies (India) Pvt. Ltd .**, 306, Agarwal City Mall, opposite M2K Pitampura ,
Delhi-110034 (India) Tel # +91-11-4220-8000;4220-8111;4220-8222 Fax# +91-11-4220-8444,Mobile#
+91-98105-21400, Email - customerservice@lifetechindia.com Web: www.lifetechindia.com



220	FT225-RNS	0.1-10UL MICRO TIP, NATURAL, STERILE, 96/RACK	1PK
221	S2613-RS	0.1-10UL MICRO TIP, NATURAL, STERILE, 96/RACK, 10R/CS	10PK
222	FT224	10-100UL FILTER TIP, NATURAL, 1000/PACK	1PK
223	FT224-RNS	10-100UL FILTER TIP, NATURAL, STERILE, 96/RACK	1PK
224	FT226	5-200UL FLITER TIP, NATURAL, 1000/PACK	1PK
225	FT226-RNS	5-200UL FLITER TIP, NATURAL, STERILE, 96/RACK	1PK
226	S2533-RS	5-200UL FLITER TIP, NATURAL, STERILE, 96/RACK, 10R/CS	10PK
227	FT222	100-1000UL FILTER TIP, NATURAL, 1000/PACK	1PK
228	FT222RS	100-1000UL FILTER TIP, NATURAL, STERILE, 100/RACK	1PK
229	S2531-RS	100-1000UL FILTER TIP, NATURAL, STERILE, 100/RACK	10PK
230	PEL10096	ELISA PLATE, FLAT BOTTOM, 96 WELL 10/BAG	1PK
231	PEL10196	ELISA PLATE, F TYPE, 96 WELL 10/BAG	1PK
232	PEL10148	ELISA PLATE, F TYPE, 48 WELL 10/BAG	1PK
233	PEL10012	ELISA PLATE, 12 STRIP 40/BAG	1PK
234	PEL10008	ELISA PLATE, 8 STRIP 60/BAG	1PK
235	PEL40096	ELISA PLATE FRAME, 12/8 STRIP 10/BAG	1PK
236	PEL40048	ELISA PLATE FRAME, 8 STRIP 10/BAG	1PK
237	BP10	0.5-10UL Pipettor	1
238	BP20	2-20UL Pipettor	1
239	BP200	20-200UL Pipettor	1
240	BP1000	100-1000UL Pipettor	1
241	RA215BR	10UL BLUE RACK	1R
242	RA224BR	200UL BLUE RACK	1R
243	RA225-5LR	5 LAYERS, 200UL BLUE RACK	1R
244	RA226BR	300UL BLUE RACK	1R
245	RA222BR	1ML BLUE RACK	1R
246	R567-N	80-WELL REVERSIBLE, NATURAL	1
247	R567-B	80-WELL REVERSIBLE, BLUE	1
248	R567-G	80-WELL REVERSIBLE, GREEN	1
249	R567-P	80-WELL REVERSIBLE, PINK	1
250	R567-V	80-WELL REVERSIBLE, PURPLE	1
251	R567-Y	80-WELL REVERSIBLE, YELLOW	1
252	R567-C	80-WELL REVERSIBLE, NATURAL COVER ONLY	1
253	R577-N	96-WELL REVERSIBLE, NATURAL	1
254	R577-B	96-WELL REVERSIBLE, BLUE	1
255	R577-G	96-WELL REVERSIBLE, GREEN	1
256	R577-P	96-WELL REVERSIBLE, PINK	1
257	R577-V	96-WELL REVERSIBLE, PURPLE	1
258	R577-Y	96-WELL REVERSIBLE, YELLOW	1
259	R577-C	96-WELL REVERSIBLE, NATURAL COVER ONLY	1
260	R541-N	96-WELL PCR TUBE RACKS, NATURAL	1
261	R541-B	96-WELL PCR TUBE RACKS, BLUE	1
262	R541-Y	96-WELL PCR TUBE RACKS, YELLOW	1
263	R541-R	96-WELL PCR TUBE RACKS, RED	1
264	R005-G	5-WELL CENTRIFUGE TUBE BOX	1
265	R011-N	11-WELL CENTRIFUGE TUBE BOX	1
266	R022-N	50-WELL MICROTUBE RACKS, NATURAL	1
267	R022-R	50-WELL MICROTUBE RACKS, RED	1
268	R022-G	50-WELL MICROTUBE RACKS, GREEN	1
269	R022-Y	50-WELL MICROTUBE RACKS, YELLOW	1
270	RH503-N	100-WELL MICROTUBE RACKS, ROUND, NATURAL	1
271	RH503-B	100-WELL MICROTUBE RACKS, ROUND, BLUE	1
272	RH503-Y	100-WELL MICROTUBE RACKS, ROUND, YELLOW	1
273	RH503-R	100-WELL MICROTUBE RACKS, ROUND, RED	1
274	RH504-N	100-WELL MICROTUBE RACKS, SQUARE, NATURAL	1
275	RH504-B	100-WELL MICROTUBE RACKS, SQUARE, BLUE	1
276	RH504-Y	100-WELL MICROTUBE RACKS, SQUARE, YELLOW	1
277	RH504-R	100-WELL MICROTUBE RACKS, SQUARE, RED	1



278	R023-N	4-WAYS MICROTUBE & CENTRIFUGE TUBES RACKS, NATURAL	1
279	R023-B	4-WAYS MICROTUBE & CENTRIFUGE TUBES RACKS, BLUE	1
280	R023-Y	4-WAYS MICROTUBE & CENTRIFUGE TUBES RACKS, YELLOW	1
281	R023-R	4-WAYS MICROTUBE & CENTRIFUGE TUBES RACKS, RED	1
282	R021-N	FLOATING MICROTUBE RACKS OCTAGON, NATURAL	1
283	R021-R	FLOATING MICROTUBE RACKS OCTAGON, RED	1
284	R021-Y	FLOATING MICROTUBE RACKS OCTAGON, YELLOW	1
285	R021-G	FLOATING MICROTUBE RACKS OCTAGON, GREEN	1
286	RF024	PLASTIC FOAM RACKS, OCTAGONAL, 83X12MM	1
287	RF025	PLASTIC FOAM RACKS, CIRCLE, 83X12MM	1
288	RF026	PLASTIC FOAM RACKS, TRIANGLE, 86X12MM	1
289	IB002A	38 PLACE COLD RACK BOX	1
290	IB003	ICE BAG	1
291	SMK112	MASK, LIGHT BLUE, 100/PK	1
292	SMFEAC4	MASK, 4-LAYERS ACTIVE CHARCOAL, 50/PK	1
293	SMFE3	MASK, 3 LAYERS, EAR LOOP, 50/PK	1
294	SMFEC3	MASK, 3 LAYERS, EAR LOOP CUP STYLE, 50/PK	1
295	SMFE2	MASK, 2 LAYERS, EAR LOOP, 100/PK	1
296	SF001-L	LAB COAT, LARGE	1
297	SF002-M	LAB COAT, MEDIUM	1
298	ZPD001	BENCH TOP LINER, 50CMX15M, 1ROLL	1
299	SP061	SPATULAS, 125MM	5
300	SP062	SPATULAS, 150MM	5
301	SP063	SPATULAS, 160MM	5
302	SP064	SPATULAS, 180MM	5
303	SP065	SPATULAS, 200MM	5
304	SP071	SPATULAS WITH WOODEN HANDLE	1
305	SP072	SCOOP STAINLESS STEEL	1
306	SP073	SPOON WITH A FLAT END	1
307	SP074	SPATULAS ROUND/TAPERED ENDS	1
308	SP075	MICRO SPATULAS FLAT AND ROUNDED ENDS	1
309	SP075S	MICRO SPATULAS ONE FLAT END	1
310	SP076	SCOOPS MATERIAL STAINLESS STEEL, LARGE	1
311	SP077	SCOOPS MATERIAL STAINLESS STEEL, MEDIUM	1
312	S003-G	STANDARD GLASS CELL WITH LID, 45 X 12.5 X 7.5 MM, INSIDE WIDTH: 10MM, 1.7ML	1
313	S004-G	STANDARD GLASS CELL WITH LID, 45 X 12.5 X 12.5 MM, INSIDE WIDTH: 10MM, 3.5ML	1
314	S005-G	STANDARD GLASS CELL WITH LID, 45 X 12.5 X 22.5 MM, INSIDE WIDTH: 10MM, 7.0ML	1
315	SQ003	STANDARD UV QUARTZ CELL WITH LID, 45 X 12.5 X 7.5 MM, INSIDE WIDTH: 10MM, 1.7ML	1
316	SQ004	STANDARD UV QUARTZ CELL WITH LID, 45 X 12.5 X 12.5 MM, INSIDE WIDTH: 10MM, 3.5ML	1
317	SQ005	STANDARD UV QUARTZ CELL WITH LID, 45 X 12.5 X 22.5 MM, INSIDE WIDTH: 10MM, 7.0ML	1
318	SQ244	MICRO UVQUARTZ CELL WITH LID, INSIDE WIDTH: 2MM, 0.7ML	1
319	SQ264	MICRO UVQUARTZ CELL WITH LID, INSIDE WIDTH: 4MM, 1.4ML	1
320	SQ274	MICRO UVQUARTZ CELL WITH LID, INSIDE WIDTH: 5MM, 1.7ML	1
321	S024-G	SEMI-MICRO UV GLASS CELL WITH BLACK WALLS AND WITH LID, INSIDE WIDTH: 4MM, 1.4ML	1
322	S144-G	SEMI-MICRO UV GLASS CELL WITH BLACK WALLS AND WITH LID, INSIDE WIDTH: 5MM, 1.7ML	1
323	SQ024	SEMI-MICRO UV QUARTZ CELL WITH BLACK WALLS AND WITH LID, INSIDE WIDTH: 4MM, 1.4ML	1
324	SQ144	SEMI-MICRO UV QUARTZ CELL WITH BLACK WALLS AND WITH LID, INSIDE WIDTH: 5MM, 1.7ML	1
325	SQ014	STANDARD UV QUARTZ CELL WITH TELFLON STOPPER	1
326	SQ294	SEMI-MICRO UV QUARTZ CELL WITH FROSTED WALLS AND WITH TELFLON STOPPER	1
327	SQ204	STANDARD QUARTZ FLUOROMETER CELL WITH LID	1
328	S204-G	STANDARD GLASS FLUOROMETER CELL WITH LID	1
329	SQ074	MICRO SELF MASKING CONTINUONS FLOWTHROUGH CELL	1
330	S074-G	MICRO SELF MASKING CONTINUONS GLASS CELL	1
331	TP011	SOFT TISSUE PAPER	1
332	ST011	EGG SHAPE, 10.2X6MM	1



333	ST012	EGG SHAPE, 15X7MM	1
334	ST013	EGG SHAPE, 20X8MM	1
335	ST014	EGG SHAPE, 25X9MM	1
336	ST015	EGG SHAPE, 28X9.5MM	1
337	ST016	EGG SHAPE, 30X10MM	1
338	ST017	EGG SHAPE, 40X15MM	1
339	ST018	EGG SHAPE, 50X20MM	1
340	ST019	EGG SHAPE, 70X25MM	1
341	ST021	OCTOGONAL SHAPE, 15X6MM	1
342	ST022	OCTOGONAL SHAPE, 18X7MM	1
343	ST023	OCTOGONAL SHAPE, 22X7MM	1
344	ST024	OCTOGONAL SHAPE, 26.3X7MM	1
345	ST025	OCTOGONAL SHAPE, 30X7MM	1
346	ST026	OCTOGONAL SHAPE, 27.5X9MM	1
347	ST027	OCTOGONAL SHAPE, 33X7MM	1
348	ST028	OCTOGONAL SHAPE, 36X9MM	1
349	ST029	OCTOGONAL SHAPE, 48.5X9MM	1
350	ST030	OCTOGONAL SHAPE, 56.5X9MM	1
351	ST031	OCTOGONAL SHAPE, 66.6X9MM	1
352	ST050	STIRRING-BAR RETRIEVER, RIGID HANDLE	1
353	SR021-1	1ML DISPOSABLE GRADUATED PLASTIC SYRINGES,STERILE	30
354	SR021-2.5	2.5ML DISPOSABLE GRADUATED PLASTIC SYRINGES,STERILE	30
355	SR021-5	5ML DISPOSABLE GRADUATED PLASTIC SYRINGES,STERILE	30
356	SR021-10	10ML DISPOSABLE GRADUATED PLASTIC SYRINGES,STERILE	30
357	SR021-20	20ML DISPOSABLE GRADUATED PLASTIC SYRINGES,STERILE	20
358	SR021-60	60ML DISPOSABLE GRADUATED PLASTIC SYRINGES,STERILE	10
359	SR031	MICRO-SAMPLING SYRINGE, 0.5UL	1
360	SR032	MICRO-SAMPLING SYRINGE, 1UL	1
361	SR033	MICRO-SAMPLING SYRINGE, 5UL	1
362	SR034	MICRO-SAMPLING SYRINGE, 10UL	1
363	SR035	MICRO-SAMPLING SYRINGE, 25UL	1
364	SR036	MICRO-SAMPLING SYRINGE, 50UL	1
365	SR037	MICRO-SAMPLING SYRINGE, 100UL	1
366	AT011	AUTOCLAVE INDICATOR TAPES	1ROLLS
367	TM025	HUMIDITY MONITOR	1
368	TM032	TIMER, ONE CHANNEL	1
369	TM033	TIMER, TWO CHANNEL	1
370	TM034	TIMER, THREE CHANNEL	1
371	TCF70-25	TISSUE CULTURE FLASK, 25ML,10/CS	1PK
372	TCF70-50	TISSUE CULTURE FLASK, 50ML,10/CS	1PK
373	TCF70-250	TISSUE CULTURE FLASK, 250ML,5/CS	1PK
374	TCP20-6	TISSUE CULTURE PLATE, 6 WELL	1
375	TCP20-12	TISSUE CULTURE PLATE, 12 WELL	1
376	TCP20-24	TISSUE CULTURE PLATE, 24 WELL	1
377	TCP20-48	TISSUE CULTURE PLATE, 48WELL	1
378	TCP20-96	TISSUE CULTURE PLATE, 96 WELL	1
379	BT613-N	0.65ML TUBE, NATURAL,1000/BAG	1PK
380	BT613-A	0.65ML TUBE, AMBER,1000/BAG	1PK
381	BT613-B	0.65ML TUBE, BLUE,1000/BAG	1PK
382	BT613-G	0.65ML TUBE, GREEN,1000/BAG	1PK
383	BT613-R	0.65ML TUBE, RED,1000/BAG	1PK
384	BT613-Y	0.65ML TUBE, YELLOW,1000/BAG	1PK
385	BT613-NS	0.65ML TUBE, NATURAL, STERILE,1000TUBES	1PK
386	BT620-N	1.5ML TUBE, NATURAL,1000/BAG	1PK
387	BT620-A	1.5ML TUBE, AMBER,1000/BAG	1PK



388	BT620-B	1.5ML TUBE, BLUE,1000/BAG	1PK
389	BT620-G	1.5ML TUBE, GREEN,1000/BAG	1PK
390	BT620-R	1.5ML TUBE, RED,1000/BAG	1PK
391	BT620-Y	1.5ML TUBE, YELLOW,1000/BAG	1PK
392	BT620-NS	1.5ML TUBE, NATURAL, STERILE,1000/BAG	1PK
393	TU629-B	1.5ML TUBE WITH CODED CAP, BLUE,1000/BAG	1PK
394	TU629-G	1.5ML TUBE WITH CODED CAP, GREEN,1000/BAG	1PK
395	TU629-P	1.5ML TUBE WITH CODED CAP, PURPLE,1000/BAG	1PK
396	TU629-R	1.5ML TUBE WITH CODED CAP, RED,1000/BAG	1PK
397	TU629-W	1.5ML TUBE WITH CODED CAP, WHITE,1000/BAG	1PK
398	TU629-Y	1.5ML TUBE WITH CODED CAP, YELLOW,1000/BAG	1PK
399	BT619-N	2ML TUBE, NATURAL,1000/BAG	1PK
400	BT619-A	2ML TUBE, AMBER,1000/BAG	1PK
401	BT619-B	2ML TUBE, BLUE,1000/BAG	1PK
402	BT619-G	2ML TUBE, GREEN,1000/BAG	1PK
403	BT619-R	2ML TUBE, RED,1000/BAG	1PK
404	BT619-Y	2ML TUBE, YELLOW,1000/BAG	1PK
405	BT619-NS	2ML TUBE, NATURAL, STERILE,1000TUBES	1PK
406	TC131-N	0.5ML CONICAL TUBE WITH SCREW CAP, NATURAL, 500/BAG	1PK
407	TC132-SN	0.5ML CONICAL TUBE WITH SCREW CAP, STERILE, NATURAL, 500/BAG	1PK
408	TC133-SB	0.5ML CONICAL TUBE WITH SCREW CAP, STERILE, BLUE,500/BAG	1PK
409	TC134-SR	0.5ML CONICAL TUBE WITH SCREW CAP, STERILE, RED,500/BAG	1PK
410	TC135-SW	0.5ML CONICAL TUBE WITH SCREW CAP, STERILE, WHITE,500/BAG	1PK
411	TC136-SY	0.5ML CONICAL TUBE WITH SCREW CAP, STERILE, YELLOW,500/BAG	1PK
412	TC137-SG	0.5ML CONICAL TUBE WITH SCREW CAP, STERILE, GREEN,500/BAG	1PK
413	TC137-A	GRADUATED, CONICAL, AMBER TUBE WITH AMBER CAPS,500/BAG	1PK
414	S2664	0.5ML SKIRTED TUBE, 500/BAG	1PK
415	TK341-N	0.5ML CONICAL TUBE WITH SCREW CAP,NATURAL,500/BAG	1PK
416	TK342-SN	0.5ML SKIRTED TUBE WITH SCREW CAP, STERILE, NATURAL,500/BAG	1PK
417	TK343-SB	0.5ML SKIRTED TUBE WITH SCREW CAP, STERILE, BLUE,500/BAG	1PK
418	TK344-SR	0.5ML SKIRTED TUBE WITH SCREW CAP, STERILE, RED,500/BAG	1PK
419	TK345-SW	0.5ML SKIRTED TUBE WITH SCREW CAP, STERILE, WHITE,500/BAG	1PK
420	TK346-SY	0.5ML SKIRTED TUBE WITH SCREW CAP, STERILE, YELLOW,500/BAG	1PK
421	TK347-SG	0.5ML SKIRTED TUBE WITH SCREW CAP, STERILE, GREEN,500/BAG	1PK
422	TK347-A	GRADUATED, SKIRTED, AMBER TUBE WITH AMBER CAPS,500/BAG	1PK
423	TC121-N	1.5ML CONICAL TUBE WITH SCREW CAP, NATURAL, 500/BAG	1PK
424	TC122-SN	1.5ML CONICAL TUBE WITH SCREW CAP, STERILE, NATURAL, 500/BAG	1PK
425	TC123-SB	1.5ML CONICAL TUBE WITH SCREW CAP, STERILE, BLUE,500/BAG	1PK
426	TC124-SR	1.5ML CONICAL TUBE WITH SCREW CAP, STERILE, RED,500/BAG	1PK
427	TC125-SW	1.5ML CONICAL TUBE WITH SCREW CAP, STERILE, WHITE,500/BAG	1PK
428	TC126-SY	1.5ML CONICAL TUBE WITH SCREW CAP, STERILE, YELLOW,500/BAG	1PK
429	TC127-SG	1.5ML CONICAL TUBE WITH SCREW CAP, STERILE, GREEN,500/BAG	1PK
430	S2670	1.5ML CONICAL TUBE , NATURAL,500/BAG	1PK
431	TC151-N	1.65ML CONICAL TUBE WITH SCREW CAP, NATURAL,500/BAG	1PK
432	TC152-SN	1.65ML CONICAL TUBE WITH SCREW CAP, STERILE, NATURAL,500/BAG	1PK
433	TC153-SB	1.65ML CONICAL TUBE WITH SCREW CAP, STERILE, BLUE,500/BAG	1PK
434	TC154-SR	1.65ML CONICAL TUBE WITH SCREW CAP, STERILE, RED,500/BAG	1PK
435	TC155-SW	1.65ML CONICAL TUBE WITH SCREW CAP, STERILE, WHITE,500/BAG	1PK
436	TC156-SY	1.5ML CONICAL TUBE WITH SCREW CAP, STERILE, YELLOW,500/BAG	1PK
437	TC157-SG	1.5ML CONICAL TUBE WITH SCREW CAP, STERILE, GREEN,500/BAG	1PK
438	TC157-A	1.5ML CONICAL TUBE WITH SCREW CAP, AMBER,500/BAG	1PK
439	S2668	1.5ML SKIRTED TUBE , NATURAL,500/BAG	1PK
440	TK351-N	1.5ML SKIRTED TUBE WITH SCREW CAP, STERILE, NATURAL,500/BAG	1PK
441	TK352-SN	1.5ML SKIRTED TUBE WITH SCREW CAP, STERILE, BLUE,500/BAG	1PK
442	TK353-SB	1.5ML SKIRTED TUBE WITH SCREW CAP, STERILE, RED,500/BAG	1PK
443	TK354-SR	1.5ML SKIRTED TUBE WITH SCREW CAP, STERILE, WHITE,500/BAG	1PK
444	TK355-SW	1.5ML SKIRTED TUBE WITH SCREW CAP, STERILE, WHITE,500/BAG	1PK
445	TK356-SY	1.5ML SKIRTED TUBE WITH SCREW CAP, STERILE, YELLOW,500/BAG	1PK



446	TK357-SG	1.5ML SKIRTED TUBE WITH SCREW CAP, STERILE, GREEN,500/BAG	1PK
447	TK357-A	1.5ML SKIRTED TUBE WITH SCREW CAP, AMBER,500/BAG	1PK
448	S2674	2ML CONICAL TUBE , NATURAL,500/BAG	1PK
449	TC161-N	2ML CONICAL TUBE WITH SCREW CAP, NATURAL,500/BAG	1PK
450	TC162-SN	2ML CONICAL TUBE WITH SCREW CAP, STERILE, NATURAL,500/BAG	1PK
451	TC163-SB	2ML CONICAL TUBE WITH SCREW CAP, STERILE, BLUE,500/BAG	1PK
452	TC164-SR	2ML CONICAL TUBE WITH SCREW CAP, STERILE, RED,500/BAG	1PK
453	TC165-SW	2ML CONICAL TUBE WITH SCREW CAP, STERILE, WHITE,500/BAG	1PK
454	TC166-SY	2ML CONICAL TUBE WITH SCREW CAP, STERILE, YELLOW,500/BAG	1PK
455	TC167-SG	2ML CONICAL TUBE WITH SCREW CAP, STERILE, GREEN,500/BAG	1PK
456	TC167-A	2ML CONICAL TUBE WITH SCREW CAP, AMBER,500/BAG	1PK
457	S2672	2ML SKIRTED TUBE,500/BAG	1PK
458	TK361-N	2ML SKIRTED TUBE WITH SCREW CAP, NATURAL,500/BAG	1PK
459	TK362-SN	2ML SKIRTED TUBE WITH SCREW CAP, STERILE, NATURAL,500/BAG	1PK
460	TK363-SB	2ML SKIRTED TUBE WITH SCREW CAP, STERILE, BLUE,500/BAG	1PK
461	TK364-SR	2ML SKIRTED TUBE WITH SCREW CAP, STERILE, RED,500/BAG	1PK
462	TK365-SW	2ML SKIRTED TUBE WITH SCREW CAP, STERILE, WHITE,500/BAG	1PK
463	TK366-SY	2ML SKIRTED TUBE WITH SCREW CAP, STERILE, YELLOW,500/BAG	1PK
464	TK367-SG	2ML SKIRTED TUBE WITH SCREW CAP, STERILE, YELLOW,500/BAG	1PK
465	TK371-N	2ML SKIRTED TUBE WITH SCREW CAP, NATURAL,500/BAG	1PK
466	TK372-SN	2ML SKIRTED TUBE WITH SCREW CAP, STERILE, NATURAL,500/BAG	1PK
467	TK373-SB	2ML SKIRTED TUBE WITH SCREW CAP, STERILE,BLUE,500/BAG	1PK
468	TK374-SR	2ML SKIRTED TUBE WITH SCREW CAP, STERILE,RED,500/BAG	1PK
469	TK375-SW	2ML SKIRTED TUBE WITH SCREW CAP, STERILE,WHITE,500/BAG	1PK
470	TK376-SY	2ML SKIRTED TUBE WITH SCREW CAP, STERILE,YELLOW,500/BAG	1PK
471	TK377-SG	2ML SKIRTED TUBE WITH SCREW CAP, STERILE,GREEN,500/BAG	1PK
472	T666-N	SCREW CAP, DEEP, NATURAL,500/BAG	1PK
473	T666-B	SCREW CAP, DEEP, BLUE,500/BAG	1PK
474	T666-R	SCREW CAP, DEEP, RED,500/BAG	1PK
475	T666-Y	SCREW CAP, DEEP, YELLOW,500/BAG	1PK
476	T666-G	SCREW CAP, DEEP, GREEN,500/BAG	1PK
477	T666-W	SCREW CAP, DEEP, WHITE,500/BAG	1PK
478	T666-A	SCREW CAP, DEEP, Amber,500/BAG	1PK
479	T667-N	SCREW CAP, NATURAL,500/BAG	1PK
480	T667-B	SCREW CAP, BLUE,500/BAG	1PK
481	T667-R	SCREW CAP, RED,500/BAG	1PK
482	T667-Y	SCREW CAP, YELLOW,500/BAG	1PK
483	T667-G	SCREW CAP, GREEN,500/BAG	1PK
484	T667-W	SCREW CAP, WHITE,500/BAG	1PK
485	T667-A	SCREW CAP, DEEP, Amber,500/BAG	1PK
486	S2651	CAPS FOR TUBES, YELLOW, 500/BAG	1PK
487	S2652	CAPS FOR TUBES, RED, 500/BAG	1PK
488	S2657	CAPS FOR TUBES, GREEN, 500/BAG	1PK
489	S2658	CAPS FOR TUBES, BLUE, 500/BAG	1PK
490	S2659	CAPS FOR TUBES, WHITE, 500/BAG	1PK
491	S2660	CAPS FOR TUBES, NATURAL,500/BAG	1PK
492	BP541-N	0.2ML PCR TUBE, FLAT CAP, NATURAL,1000/BAG	1PK
493	BP541-B	0.2ML PCR TUBE, FLAT CAP, BLUE,1000/BAG	1PK
494	BP541-G	0.2ML PCR TUBE, FLAT CAP, GREEN,1000/BAG	1PK
495	BP541-R	0.2ML PCR TUBE, FLAT CAP, RED,1000/BAG	1PK
496	BP541-Y	0.2ML PCR TUBE, FLAT CAP, YELLOW,1000/BAG	1PK
497	BP542-N	0.2ML PCR TUBE, DOMED CAP, NATURAL,1000/BAG	1PK
498	Q541	0.2ML PCR TUBE, DOMED CAP, NATURAL,1000/BAG	1PK
499	BP541-S100	0.2ML PCR TUBE, FLAT CAP, STERILE,1000/BAG	1PK
500	BP542-S100	0.2ML PCR TUBE, DOMED CAP, STERILE,100/BAG	1PK
501	BP602-N	0.5ML PCR TUBE, FLAT CAP, NATURAL,1000/BAG	1PK



502	BP602-B	0.5ML PCR TUBE, FLAT CAP, BLUE,1000/BAG	1PK
503	BP602-G	0.5ML PCR TUBE, FLAT CAP, GREEN,1000/BAG	1PK
504	BP602-R	0.5ML PCR TUBE, FLAT CAP, RED,1000/BAG	1PK
505	BP602-Y	0.5ML PCR TUBE, FLAT CAP, YELLOW,1000/BAG	1PK
506	BP602-S	0.5ML PCR TUBE, FLAT CAP, NATURAL, STERILE,1000/BAG	1PK
507	BP602-S100	0.5ML PCR TUBE, FLAT CAP, NATURAL, STERILE, 100/PK, 10/BAG	1PK
508	BP546-N	8 STRIP 0.2ML PCR TUBE, 120/BAG	1PK
509	BP546-NS	8 STRIP 0.2ML PCR TUBE, STERILE, 120/BAG	1PK
510	BP551	8 STRIP CAP, NATURAL,120/BAG	1PK
511	BP551-NS	8 STRIP CAP, NATURAL, STERILE,120/BAG	1PK
512	BP546-B	8 STRIP 0.2ML PCR TUBE, STERILE, BLUE,120/BAG	1PK
513	BP546-G	8 STRIP 0.2ML PCR TUBE, STERILE, GREEN, 120/BAG	1PK
514	BP546-R	8 STRIP 0.2ML PCR TUBE, STERILE, RED, 120/BAG	1PK
515	BP546-Y	8 STRIP 0.2ML PCR TUBE, STERILE, YELLOW, 120/BAG	1PK
516	PCR-960-CH-S	96 WELL 0.2ML PCR PLATE, CHIMNEY-TOP 10/BAG, STERILE	1PK
517	460-96	96 WELL 200UL PCR, 10/BAG, STERILE	1PK
518	PCR-384-SN	384 WELL PLATE, 10/BAG	1PK
519	BR580-96N	96 WELL 2ML DEEP PLATE, NATURAL,10PLATES/BAG	1PK
520	BR580-96B	96 WELL 2ML DEEP PLATE, BLUE,10PLATES/BAG	1PK
521	BR580-96R	96 WELL 2ML DEEP PLATE, RED,10PLATES/BAG	1PK
522	BR580-96Y	96 WELL 2ML DEEP PLATE, YELLOW,10PLATES/BAG	1PK
523	SD5010	96-WELL 200UL U-SHAPED DEEP-WELL PLATES, STERILE,10PLATE/BAG	1PK
524	BR480-N	48-WELL 4ML U-SHAPED DEEP-WELL PLATES, 10/BAG	1PK
525	TU6416	384 WELL 80UL SMAPLE PLATE,10PLATES/BAG	1PK
526	TB418	SEALING FILM,50SHEETS	1PK
527	TB418-S	SEALING FILM, STERILE,50SHEETS	1PK
528	TSM113	PLATE ROLLER	1
529	TD535	5ML TUBE, DUAL CAP, POLYPROPYLENE,100/BAG	1PK
530	TD535-S	5ML TUBE, DUAL CAP, STERILE, POLYPROPYLENE,50/BAG	1PK
531	TP536-N	5ML TUBE, PLUG CAP, NATURAL, POLYPROPYLENE,100/BAG	1PK
532	TP536-NS	5ML TUBE, PLUG CAP, NATURAL,STERILE, POLYPROPYLENE,50/BAG	1PK
533	TP536-B	5ML TUBE, PLUG CAP, BLUE, POLYPROPYLENE,100/BAG	1PK
534	TP536-SB	5ML TUBE, PLUG CAP, BLUE,STERILE, POLYPROPYLENE,50/BAG	1PK
535	TP536-R	5ML TUBE, PLUG CAP, RED, POLYPROPYLENE,100/BAG	1PK
536	TP536-SR	5ML TUBE, PLUG CAP, RED,STERILE, POLYPROPYLENE,50/BAG	1PK
537	TP536-Y	5ML TUBE, PLUG CAP, YELLOW, POLYPROPYLENE,100/BAG	1PK
538	TP536-SY	5ML TUBE, PLUG CAP, YELLOW,STERILE, POLYPROPYLENE,50/BAG	1PK
539	TD444	12ML TUBE, DUAL CAP, POLYPROPYLENE,100/BAG	1PK
540	TD444-S	12ML TUBE, DUAL CAP, POLYPROPYLENE,50/BAG,STERILE	1PK
541	BT693-N	5ML TUBE & SCREW CAP, NATURAL, POLYPROPYLENE,100/BAG	1PK
542	BT693-SN	5ML TUBE & SCREW CAP, NATURAL, POLYPROPYLENE, STERILE,100/BAG	1PK
543	BT693-SB	5ML TUBE & SCREW CAP, BLUE, POLYPROPYLENE, STERILE,100/BAG	1PK
544	BT693-SR	5ML TUBE & SCREW CAP, RED, POLYPROPYLENE, STERILE,100/BAG	1PK
545	BT693-SY	5ML TUBE & SCREW CAP, YELLOW, POLYPROPYLENE, STERILE,100/BAG	1PK
546	BT694-N	10ML TUBE & SCREW CAP, NATURAL, POLYPROPYLENE,100/BAG	1PK
547	BT694-SN	10ML TUBE & SCREW CAP, NATURAL, POLYPROPYLENE, STERILE,100/BAG	1PK
548	BT694-SB	10ML TUBE & SCREW CAP, BLUE, POLYPROPYLENE, STERILE,100/BAG	1PK
549	BT694-SR	10ML TUBE & SCREW CAP, RED, POLYPROPYLENE, STERILE,100/BAG	1PK
550	BT694-SY	10ML TUBE & SCREW CAP, YELLOW, POLYPROPYLENE, STERILE,100/BAG	1PK
551	BT888-N	15ML CENTRIFUGE TUBE, SCREW CAP, NATURAL,100/BAG	1PK
552	BT888-SN	15ML CENTRIFUGE TUBE, SCREW CAP, NATURAL, STERILE,100/BAG	1PK
553	BT888-SB	15ML CENTRIFUGE TUBE, SCREW CAP, BLUE, STERILE,100/BAG	1PK
554	BT888-SR	15ML CENTRIFUGE TUBE, SCREW CAP, RED, STERILE,100/BAG	1PK
555	BT888-SY	15ML CENTRIFUGE TUBE, SCREW CAP, YELLOW, STERILE,100/BAG	1PK
556	CT788-Y	50ML CONICAL TUBES, YELLOW CAP,25/BAG	1PK
557	CT788-G	50ML CONICAL TUBES, GREEN CAP, 25/BAG	1PK



LIFE TECHNOLOGIES (INDIA) PVT. LTD.

your molecular & cell technology partner

Bringing zest to your research



BIO BASIC INC

558	CT788-YS	50ML CONICAL TUBES, YELLOW CAP,STERILE,25/BAG	1PK
559	CT788-GS	50ML CONICAL TUBES, GREEN CAP, STERILE,25/BAG	1PK
560	TK943-Y	50ML SKIRTED CENTRIGUE TUBE, YELOW CAP,25/BAG	1PK
561	TK943-G	50ML SKIRTED CENTRIGUE TUBE, YELOW CAP,25/BAG	1PK
562	TK943-YS	50ML SKIRTED CENTRIGUE TUBE, YELLOW CAP, STERILE,25/BAG	1PK
563	TK943-GS	50ML SKIRTED CENTRIGUE TUBE, YELLOW CAP, STERILE,25/BAG	1PK
564	CT790-N	50ML HIGH SPEED TUBE, NATURAL CAP,25/BAG	1PK
565	CT790-SN	50ML HIGH SPEED TUBE, NATURAL CAP, STERILE,25/BAG	1PK
566	WT111	WEIGHING PAPER, 3X3,100SHEETS/BAG	1PK
567	WT112	WEIGHING PAPER, 4X4,100SHEETS/BAG	1PK
568	WT113	WEIGHING PAPER, 6X6,100SHEETS/BAG	1PK
569	WT121	WEIGHING PAPER, SQUARE, 13/8X13/8X5/16,200SHEETS/BOX	1PK
570	WT122	WEIGHING PAPER, SQUARE, 7/2X2/2X1,50SHEETS/BOX	1PK
571	WT123	WEIGHING PAPER, SQUARE, 11/2X11/2X7/8,50SHEETS/BOX	1PK