

Petaka G3™

Transforming the art of Cell Culturing

Device Details

- *Dimensions : 124mm× 84mm× 5mm (Handheld device just the size of a microtitre plate)*
- *Cells are grown on both sides*
- *Cell culture surface:150 cm²*
- *20 ml Media Capacity*
- *Cell growth is optimized by providing adequate surface energy (Dynes) on the cell growth surfaces and by maintaining a balanced amount of Oxygen dissolved in the media*
- *The yield depending on specific cell type is approximately 20- 30 × 10⁶ cells per Petaka*
- *PetakaG3™ is a robust and sturdy device if dropped on a hard surface, a certain percentage of cells could detach and after a few minutes of incubation, those cells will generally reattach without incident*



India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura,
Delhi – 110034 (INDIA).

Ph: +91-11-42208000, 42208111, 42208222

Mobile: +91-9810521400

Fax: +91-11-42208444

Email: customerservice@lifetechindia.com

Web: www.lifetechindia.com