

# Why Purchase Two Separate Kits, When a Single Kit Can Serve Your Purpose Most Efficiently.....????

## HiYield™ Gel/PCR DNA Extraction Kit

...for both Gel & PCR DNA Extraction

- Up to 95% Recovery from PCR Product
- Unique dual purpose application. (Gel/PCR DNA Extraction)
- Without phenol extraction or alcohol precipitation
- 20 minutes for Gel Extraction/15 minutes for PCR Cleanup

[Click Here for more details](#)

**1 Kit,  
2 Purposes  
@1/2 Price  
@1/2 Shelf  
Space**



Just Look at the Cross Comparison of kit with Qiagen.....

Company	RBC		QIAGEN			
Product Name	HiYield™ Gel/PCR DNA Extraction Kits		QIAquick Gel Extraction Kits		QIAquick PCR Purification Kits	
Size	Mini		Mini			
Cat No.	YDF100	YDF300	28704	28706	28104	28106
Preparations	100	300	50	250	50	250
DNA fragment	50 bp~10 kb		70 bp~10 kb		100 bp~10 kb	
Sample amount	300 mg agarose gel / 100µl PCR product		400 mg		10µl	
Recovery	Gel:80~90%, PCR:90~95%		70~80%		90~95%	
Operation time	Gel:20 mins, PCR:15 mins					
Price	~ 40% of Qiagen					