

Hepatitis B Virus Vaccines & ELISA Kits

name	Product Type#	Cat #	Product Description
HBsAg Antigen and Antibody kits	Antigen kit	4100	Hepatitis B Surface Antigen (HBsAg) ELISA kit, Quantitative, 96 tests Stds=0-100 ng/ml ; Sample=50 ul; Sensitivity= 3 ng/ml; 75 min assay
	Anti-HBsAg	4200	Human Anti-Hepatitis B Surface Antigen (HBsAg) ELISA kit, 96 tests
	Anti-HBsAg	4210	Mouse Anti-Hepatitis B Surface Antigen (HBsAg) ELISA kit, 96 tests

India Contact:
Life Technologies (India) Pvt. Ltd.

 306, Aggarwal City Mall, Opposite M2K Pitampura,
 Delhi – 110034 (INDIA).

Ph: +91-11-42208000, 42208111, 42208222

Mobile: +91-9810521400

Fax: +91-11-42208444

 Email: customerservice@lifetechindia.com

 Web: www.lifetechindia.com