



Insoluble proteins left you hanging?

Make soluble proteins the easy way!

Espresso™ T7 SUMO CLONING AND EXPRESSION SYSTEM

- Directional cloning of PCR amplicons in seconds
- Ready-to-use vector and competent cells — no ligation step!
- Tightly controlled protein expression
- Cleavable solubility tag

Find out more! www.lucigen.com

LIFETECH

LIFE TECHNOLOGIES (INDIA) PVT. LTD.

your molecular & cell technology partner

Bringing zest to your research

Lucigen®



Expresso™ T7 SUMO Cloning & Expression System

Fastest & easiest SOLUBLE
protein expression solution available

Save days of vector and
cell preparation time!

Unbelievable performance —
inconceivable price!

▶ Find out more, visit
www.lucigen.com



Hotter, Faster Results

Lucigen®

Advanced Products for Molecular Biology



India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura,

Delhi – 110034 (INDIA).

Ph: +91-11-42208000, 42208111, 42208222

Mobile: +91-9810521400

Fax: +91-11-42208444

Email: customerservice@lifetechindia.com

Web: www.lifetechindia.com