



HiYield™ Plasmid Mini Kit for Increased Yield of Plasmid DNA

-High, reproducible recovery of plasmid DNA
-Rapid and simple procedure
-No phenol, chloroform or alcohol

[Click Here For More Details](#)



Just Look at the Cross Comparison of kit with Qiagen.....

Company	RBC		QIAGEN	
Product Name	HiYield™ Plasmid Mini Kit		QIAquick Gel Extraction Kits	
Size	Mini		Mini	
Cat No.	YPD00	YPD300	12123	12125
Preparations	100	300	25	100
Sample class	Bacterial cultures		Bacterial culture in LB medium	
Sample amount	1~4ml		2~5ml	
Yield	low-copy number 3~10µg ,high-copy number 20~30µg		High-copy number 20µg	
Operation time	20 mins		25 mins	
Price	~ 40% of Qiagen			