



Total RNA Isolation from 4-15 Hybridoma cells that were saved in **RNA SAVE** Solution (Cat. No: 01-891-1B) during 4 days at: +4°C, -20°C and -196°C (liquid Nitrogen).



- 1) DNA Ladder
- 2) + 4 ° C
- 3) - 20 ° C
- 4) -196°C (liquid Nitrogen)

=====

Total RNA Isolation from 3D2C2 Hybridoma cells that were saved in **RNA SAVE** Solution (Cat. No: 01-891-1B) during 4 days at: +4°C, -20°C and -196°C (liquid Nitrogen).



- 1) +4°C
- 2) -20°C
- 3) -196°C (liquid Nitrogen)