


### RNA Save

Sample storage solution for RNA stabilization  
Cat. No. 01- 891-1



#### Description

RNA Save is a solution that simplifies the collection of RNA samples from operating rooms and other field locations; it protects and stabilizes the RNA in tissue and cells without need for freezing, thus enabling easy transportation and storage.

#### Features

- Easy to use - collect sample (do not freeze), store in RNA Save at 4°C over night and then freeze
- Scalable
- One solution
- Working with variant samples - cells and tissue
- Working with variant origin - animal, plant and bacteria.

**Storage Conditions of Samples: Keep the samples in RNA Save solution and store them at -20°C or -80°C indefinitely**

#### Related Products

- RNase Exitus Plus™ - for RNase decontamination (cat. 01-895-1B, 01-895-1L )
- EZ-RNA - Total RNA Isolation Kit (cat. 20-400-100)
- EZ-First Strand - cDNA Synthesis Kit for RT-PCR (cat. 20-800-50)