

Cat No. L-201

Product Description

Life Technologies India DNA Ladder L-201 is a ready to use mix of 13 proprietary DNA fragments ranging from 250 to 10,000 base pairs, suitable for use as molecular standards for agarose gel electrophoresis. The mix contains bromophenol blue and Xylene Cyanol FF for tracking. The mix is recommended to be stored at 4°C for use within 12 months and should be stored at -20°C for use within 24 months.

Concentration: 100 µg per ml

Package: 50 µg per 500ul

Recommended Load: 5µl per well

