

FIREScript® KIT



Some applications of this product may require a license which is not provided by the purchase of this product.

For research use only.

Fast and robust reverse transcriptase for wide range of applications

93 /100 1 Citation

Genetically modified MMLV-based reverse transcriptase with increased thermostability and improved performance at elevated temperatures

high specificity and yield

wide reaction temperature from 37°C to 60°C

fast 15 min reaction time

a simple 2-vial kit with enzyme and reaction buffer only

reaction set-up and shipment at without ice

Ordering

Choose Product Size

50 rxn | 50 x 20 µl rxn

200 rxn | 200 x 20 µl rxn

20 rxn | 20 x 20 µl rxn **free sample**

REQUEST FOR BULK SIZE

Reagents Provided

Item	Pcs.	Vial size
FIREScript® Reverse Transcriptase	1	10000 U 50 µl
10x RT Reaction Buffer with DTT	1	125 µl

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura,
Delhi – 110034 (INDIA).Ph: +91-11-42208000, 42208111, 42208222

Mobile: +91-9810521400

Fax: +91-11-42208444

Email: customerservice@lifetechindia.com

Web: www.lifetechindia.com