

Product Data Sheet

**Anti-Human IgM (mu chain specific) Antibodies and conjugates**

<b>Cat#</b> 10219	Goat Anti-human IgM, unlabeled	<b>Size:</b> 0.5 mg	<b>Form:</b> Liquid	Powder
<b>Cat#</b> 10220	Goat Anti-human IgM-HRP conjugate	<b>Size:</b> 0.5 ml	<b>Form:</b> Liquid	Powder
<b>Cat#</b> 10230	Goat Anti-human IgM-FITC conjugate	<b>Size:</b> 0.5 ml	<b>Form:</b> Liquid	Powder
<b>Cat#</b> 10240	Goat Anti-human IgM-Biotin conjugate	<b>Size:</b> 0.5 ml	<b>Form:</b> Liquid	Powder
<b>Cat#</b> 10250	Goat Anti-human IgM-AP conjugate	<b>Size:</b> 0.5 ml	<b>Form:</b> Liquid	Powder

**Description**

Goats were immunized with antigen grade **human IgM**. Antibodies have been isolated using ammonium sulfate, ion-exchange, and affinity chromatography to remove non-specific antibodies. Specificity of antibody has been tested using IEP, immunodiffusion, and ELISA. The product reacts with human IgM with no significant reactivity with human IgG, or other serum proteins. The antibody may recognize other species IgM. Purified anti-human IgM antibody is supplied as unlabeled or conjugated to Biotin, FITC, AP or HRP.

and 0.05% azide in either **lyophilized** (0.5 ml) or **liquid** form (0.5 mg/0.5 ml). Reconstitute powder in PBS in 0.5 ml to prepare 1 mg/ml solution. Store at -20°C in suitable aliquots. Stability is ~6-12 months. Do not freeze and thaw.

Suggested conjugate dilutions are 1:200-1:2000 for immunofluorescence.

**Absorption Wavelength:** 495 nm

**Emission Wavelength:** 528 nm

**Cat# 10240, Anti-human IgM-Biotin-conjugate**

**Form and Storage**

**Cat# 10219, Anti-human IgM, unlabeled**

The antibody is supplied in PBS, pH 7.4, and 0.05% azide in either **lyophilized** (0.5 mg) or **liquid** form (0.5 mg/0.5 ml). Reconstitute powder in PBS in 0.5 ml to prepare 1 mg/ml solution. Store at -20°C in suitable aliquots. Stability is ~6-12 months. Do not freeze and thaw.

Purified anti-human IgM antibody was coupled to Biotin using Biotinamidocaproate N-Hydroxysuccinimide Ester (BAC) at F/P ratio ~10-20:1. The antibody is supplied in PBS, pH 7.4, 0.2% BSA and 0.05% azide in either **lyophilized** (0.5 mg) or **liquid** form (0.5 mg/0.5 ml). Reconstitute powder in PBS in 0.5 ml to prepare 1 mg/ml solution. Store at -20°C in suitable aliquots. Stability is ~6-12 months. Do not freeze and thaw.

Suggested conjugate dilutions are 1:5,000-1:30,000 ELISA, 1:2K-1:10K for western.

**Cat# 10220, Anti-human IgM-HRP-conjugate**

Purified antibody was coupled to HRP (RZ>3.0) using periodate method. The molar enzyme to protein (E/P) ratio = 4.0. The antibody is supplied in stabilizing buffer, 0.1% proclin-300 as preservative in , either **lyophilized** (0.5 ml) or **liquid** form (0.5 mg/0.5 ml). Reconstitute powder in PBS in 0.5 ml. Store at 4°C in suitable aliquots. Stability is ~6-12 months. Do not freeze and thaw.

**Recommended Working Dilution for ELISA**

Working dilution for the specific application should be determined by the investigator to obtain the best conditions. Working solution should be prepared immediately before use and diluted solution should be discarded.

Suggested conjugate dilutions are 1:1,000-1:10,000 ELISA, 1:1K-1:5K for western, and 1:200-1:1000 (IHC).

**Related Material available for ADI**

**Cat# 10250, Anti-human IgM,-AP-conjugate**

The conjugate is provided at ~0.5-1 mg/ml as liquid in a stabilizing buffer (50 mM HEPES pH 7.1, 0.1 M NaCl, 1 mM MgCl<sub>2</sub>, 0.1 mM ZnCl<sub>2</sub>, pH 7.5, containing 0.2% bovine serum albumin, 0.05% sodium azide). The product should be **stored at 4°C** and is stable for a minimum of 1 year. Do not store diluted solutions.

**Human IgM, IgG, IgM, and IgM ELISA kits**

Anti-human IgM, IgG, IgM, IgM antibodies and conjugates

ELISA kits for the detection of human Antibodies

**Western blot recycling kit** (Use the same blot/strip to probe with multiple antibodies) **Strip Antibodies in 5 minutes** (no heating or mercaptoethanol)

Suggested conjugate dilutions are 1:1,000-1:10,000 ELISA, 1:1K-1:5K for western, and 1:200-1:1000 (IHC).

**MSDS:** This product is considered non-hazardous as defined by OSHA (CFR 1910.1200, Nov 25, 1983). It can be disposed of in the drain.

**Cat# 10230, Anti-human IgM-FITC-conjugate**

Purified anti-human IgM antibody was coupled to FITC at F/P ratio ~3:7. The antibody is supplied in PBS, pH 7.4, 0.2% BSA

10219, 10220, 10230, 10240 10250 70219A

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)