

India Contact:

Life Technologies (India) Pvt. Ltd.
Mobile: +91-9810521400, Ph: +91-11-42208000
Email: customerservice@lifetechindia.com
Web: www.lifetechindia.com

Product Information

Mouse Peripheral Blood Mononuclear Cells (mPBMCs)

Catalog Number	10MU-031	Cell Number	5 million cells/vial
Species	<i>Mus Musculus</i>	Storage Temperature	Liquid Nitrogen

Description

Peripheral Blood Mononuclear Cells (PBMCs) is any blood cell with a single nucleus [1]. These blood cells are a critical component in the immune system to fight pathogen. PBMCs consist of lymphocytes (T, B and NK cells), monocytes and dendritic cells [2]. PBMCs are widely used in research fields including immunology, infectious diseases and hematological malignancies etc.

iXCells Biotechnologies provides high quality Mouse Peripheral Blood Mononuclear Cells (mPBMCs), which are isolated from mouse peripheral blood and cryopreserved with 5 million cells in each vial.

Product Details

Tissue	Mouse peripheral blood
Package Size	5 million cells/vial
Shipped	Cryopreserved
Storage	Liquid nitrogen
Growth Properties	Suspension
Media	Blood Cell Culture Medium (Cat # MD-0007)

Protocols

Thawing of Frozen Cells

1. Upon receipt of the frozen cells, it is recommended to thaw the cells and initiate the culture immediately in order to retain the highest cell viability.
2. To thaw the cells, put the vial in 37°C water bath with gentle agitation for ~1 minute. Keep the cap out of water to minimize the risk of contamination.
3. Pipette the cells into a 15 mL conical tube with 5mL fresh **Blood Cell Culture Medium*** (Cat # MD-0007).
4. Centrifuge at 2,000 rpm (~450 g) for 5 minutes under room temperature.
5. Remove the supernatant and cells are ready for use.

**Blood Cell Culture Medium (Cat #MD-0007) doesn't contain active ingredients to active T cells. It is recommended for immune cell recovery and can be used as basal medium for T cell activation.*

Safety Precaution: *it is highly recommended that protective gloves and clothing should be used when handling frozen vials.*

References

[1] Delves, Peter, et al. Roitt's Essential Immunology, 11th Ed.

[2] Miyahira, Andrea (22 Nov 2012). "Types of immune cells present in human PBMC". *sanguinebio.com*. Retrieved 23 Sep 2014.

Disclaimers

This product is intended for laboratory research purposes only. It is not intended for use in humans. While iXCells Biotechnologies uses reasonable efforts to include accurate and up-to-date information on this product sheet, we make no warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. iXCells Biotechnologies does not warrant that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, and use. iXCells Biotechnologies is not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to insure authenticity and reliability of strains on deposit, iXCells Biotechnologies is not liable for damages arising from the misidentification or misrepresentation of cultures.

© iXCells Biotechnologies 2015. All rights reserved.