

Product Description

iXCells Biotechnologies provides high quality Rat Pulmonary Artery Smooth Muscle Cells (RPASMC), which are isolated from adult rat pulmonary artery, and cryopreserved at P1, with >0.5 million cells in each vial. RPASMC are positive for smooth muscle cell specific alpha-actin expression and they are negative for mycoplasma, bacteria, yeast, and fungi. RPASMC can further expand for 10 population doublings in Rat Smooth Muscle Cell Growth Medium (Cat# MD-0068) under the condition suggested by iXCells Biotechnologies.

Product Details

Tissue	Adult rat pulmonary artery
Package Size	0.5x10 ⁶ cells/vial
Passage Number	P2
Shipped	Cryopreserved
Storage	Liquid nitrogen
Growth Properties	Adherent
Media	Rat Smooth Muscle Cell Growth Medium (Cat# MD-0068)

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura,
Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000

Email: customerservice@lifetechindia.com

Web: www.lifetechindia.com