

## Product Description

Peripheral Blood Mononuclear Cells (PBMCs) is any blood cell with a single nucleus [1]. PBMCs consist of lymphocytes (T, B and NK cells), monocytes and dendritic cells [2]. PBMCs are widely used in research fields including immunology, infectious diseases and hematological malignancies etc.

iXCells Biotechnologies provides high quality Rat Peripheral Blood Mononuclear Cells (rPBMC), which are isolated from Sprague Dawley rat peripheral blood and cryopreserved.

## References

[1] Delves, Peter, et al. Roitt's Essential Immunology, 11th Ed.

[2] Miyahira, Andrea (22 Nov 2012). "Types of immune cells present in human PBMC". sanguinebio.com. Retrieved 23 Sep 2014.

India Contact:

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura,  
Delhi – 110034 (INDIA).

Mobile: +91-9810521400, Ph: +91-11-42208000

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com)

Web: [www.lifetechindia.com](http://www.lifetechindia.com)