

Product Data Sheet

Mouse IgG Fab isotype controls and conjugates

<input type="checkbox"/> #20008-1-Fab	Mouse IgG Fab, unlabeled, isotype control	Size: 1 mg
<input type="checkbox"/> # 20008-Fab-B	Mouse IgG Fab-Biotin conjugate, isotype control	Size: 100 ug
<input type="checkbox"/> # 20008-Fab-F	Mouse IgG Fab-FITC conjugate, isotype control	Size: 100 ug
<input type="checkbox"/> # 20008-Fab-HP	Mouse IgG Fab-HRP conjugate, isotype control	Size: 100 ug

Description

A healthy non-immunized population of normal Mouse was used to collect serum. This product was prepared from normal serum by a multi-step process which includes delipidation, salt fractionation and ion exchange chromatography followed by papain digestion and extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Mouse Serum, anti-Mouse IgG and anti-Mouse IgG F(ab')₂. No reaction was observed against anti-Mouse IgG F(c) or anti- Papain.

Highly purified preparation of Mouse IgG Fab is supplied as unlabeled, HRP-, Biotin, FITC, and PE, Cy5 conjugates for use as isotype control or negative control for all mice monoclonal with IgG isotype.

Form and Storage

Cat# 20008-1-Fab, unlabeled IgG

IgG is supplied in PBS, pH 7.4, and 0.05% azide in either **lyophilized** or **liquid** form (1 mg/ml or see lot sp. concn on the vial). Reconstitute powder in PBS as desired. Store at -20oC in suitable aliquots. Stability is ~6-12 months. Do not freeze and thaw. Controls are typically tested at 0.1-1 ug/ml depending upon the application.

Cat# 20008-Fab-B, Biotin-conjugate

Purified Mouse IgG Fab was coupled to Biotin using ADI Biotinylation kit # 80300 at F/P Mouseio ~10-20:1. The antibody is supplied in PBS, pH 7.4, 1 % BSA in either **lyophilized** or **liquid** form (see lot sp. concn on the vial). Reconstitute powder in PBS to prepare stock solutions. Store at -20oC in suitable aliquots. Stability is ~6-12 months. Do not freeze and thaw. No preservatives have been added. Sodium azide (0.1%) or merthiolate (0.05%) can be added if necessary to prevent bacterial growth.

Suggested conjugate dilutions are 1:5,000-1:30,000 ELISA. For other applications use as needed.

Cat# 20008-Fab-F, FITC-conjugate

Purified Mouse IgG Fab was coupled to FITC at F/P ratio ~4-5:1. The antibody is supplied in PBS, pH 7.4, azide in either **lyophilized** or **liquid** form (see see lot sp. concn on the vial). Reconstitute powder in PBS as desired. Store at 4oC in suitable aliquots. Stability is ~6-12 months. Do not freeze and thaw.

Suggested conjugate dilutions are 1:50-1:500 for immunofluorescence or as optimized for a given tests.

Color: Green
Absorption: 495 nm **Emission:** 528 nm

Cat# 20008-Fab-HP, HRP-conjugate

Purified Mouse IgG Fab was coupled to HRP (RZ>3.0) using periodate method. The molar enzyme to protein (E/P) ratio 4.0. The antibody is supplied in stabilizing buffer, 0.1% proclin-300 as preservative in either **lyophilized** (50 tests per vial) or **liquid** form (50 tests/100 ul or approx. 50 ug/100ul). Reconstitute powder in PBS in 0.1 ml. Store at 4oC in suitable aliquots. Stability is ~6-12 months. Do not freeze and thaw.

Suggested conjugate dilutions are 1:1,000-1:10,000 ELISA, 1:1K-1:5K for western, and 1:200-1:1000 (IHC).

All products are for In vitro research use only. PE/Cy™5 tandem conjugate. **Cy™5, Cy™5.5 and Cy™7 are trademarks of Amersham Biosciences Corp. Cy™5, Cy™5.5 and Cy™7 are for non-commercial research use only, not for therapeutic or in vivo applications

Mouse IgG Fab (IgG-3), IgA, IgG, and IgE ELISA kits

20008-Fab-con 110112A

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com