
Product Data Sheet

Normal Rabbit Serum IgG, purified (-ve control for flow cytometry)

Cat # 20009-1-100 **Size:** 100 tests (0.5 ml)
Form: Liquid Lyophilized **Storage:** Store at -20oC

Description

A healthy non-immunized population of adult Rabbit (mixed sex) was used to collect serum. IgG was prepared by ammonium sulfate, and column chromatography (Sephadex, and DEAE-ion exchange columns). The resulting preparation was judged to be >96% pure by IEP and SDS-Page.

Form and Storage

The normal Rabbit IgG is provided as liquid (100 ug/0.5 ml in PBS pH 7.4+0.05% sodium azide as preservative) or in lyophilized form (amount of protein stated on the vial). The powder product should be **stored at 4°C** for short term and -20oC for long term storage. It is stable for a minimum of 1 year. Do not store diluted solutions. Sodium azide, an inhibitor of peroxidase, must be avoided if used with antibody-HRP conjugates.

The **lyophilized products** should be dissolved in PBS, pH 7.4 to prepare desired concentration by gentle rocking or vortexing at room temperature. Centrifuge the IgG solution briefly (10,000 g, 15 min) to remove any particulate matter. It should be aliquoted and stored frozen for long term use.

Suggested Uses

This preparation of normal Rabbit IgG is suitable for flow cytometry. Recommended testing concn. is approx. 5 ul test or 1 ug/test. Optimal concentration should be tested and adjusted if necessary. This preparation may be unsuitable for some cells and may give slightly higher background. Switching to another preparation from different Rabbit or vendor may alleviate this problem.

This product is for in vitro research use only.

MSDS: A Material Safety Data Sheet is not required for this product. The product does not contain any hazardous components above 1% or any carcinogens above 0.1% as defined in 29 CFR 1910.1200, the OSHA Hazard Communication Standard. The product contains natural purified proteins in a buffer containing 0.1% or less sodium azide. MSDS can be consulted for these items. Please note that although the product is defined as not hazardous, it is still advisable to follow prudent laboratory practices when handling laboratory reagents.

This material is sold for research purposes only and is not required to appear on the TSCA inventory. It is not intended for food, diagnostic, drug, household, agricultural or cosmetic use. Its use must be supervised by a technically qualified individual experienced in handling potentially hazardous chemicals.

Related Items available from ADI

Rabbit IgG (total), IgG, IgM, IgE, IgA ELISA kits
Anti-Rabbit IgG antibody detection ELISA kits
2nd antibodies to Rabbit IgG and Isotype

20009-1-100 80522A

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com