

Product Data Sheet

**Chicken IgG (Fc) Isotype Control and Conjugates**

<input type="checkbox"/> <b>Cat#</b> 20010-3-UL	Chicken IgG (Fc), Unlabeled	<b>Size:</b> 100 ug
<input type="checkbox"/> <b>Cat#</b> 20010-3-B	Chicken IgG (Fc)-Biotin Conjugate	<b>Size:</b> 100 ug
<input type="checkbox"/> <b>Cat#</b> 20010-3-F	Chicken IgG (Fc)-FITC Conjugate	<b>Size:</b> 100 ug
<input type="checkbox"/> <b>Cat#</b> 20010-3-HP	Chicken IgG (Fc)-HRP Conjugate	<b>Size:</b> 100 ug

Chicken IgG (Fc) was purified by a combination of column chromatography (Sephadex, and DEAE-ion exchange columns) and immunoaffinity chromatography. The resulting preparation was judged to be >96% pure by immuno electrophoresis (IEP) and SDS-PAGE. It is supplied as unlabeled or coupled to biotin, HRP, FITC, etc.

Suggested conjugate dilutions are 1:20-1:500 for immunofluorescence. Or 1 test equivalent of antibody per million cells.

Absorption @495 nm Emission @528 nm

**Cat# 20010-3-Fc, unlabeled**

The antibody is supplied in PBS, pH 7.4, and 0.05% azide in either **lyophilized** or **liquid** form (1 mg/ml or see lot sp. conc on the vial). Reconstitute powder in PBS. Store at -20oC in suitable aliquots. Stability is ~6-12 months. Do not freeze and thaw.

**Cat# 20010-3-HP, HRP-conjugate**

Purified antibody was coupled to HRP (RZ>3.0) using periodate method. The molar enzyme to protein (E/P) ratio = 4.0. The antibody is supplied in stabilizing buffer, 0.1% prolcin-300 as preservative in either lyophilized (0.5 ml) or liquid form (0.5-1 mg/ml or see lot sp. conc on the vial). Reconstitute powder in PBS in 1 ml. Store at 4oC in suitable aliquots. Stability is ~6-12 months. Do not freeze and thaw.

**Cat# 20010-3-B, Biotin-conjugate**

Purified antibody was coupled to Biotin using Biotinamidocaproate N-Hydroxysuccinimide Ester (BAC) at F/P ratio ~10-20:1. The antibody is supplied in PBS, pH 7.4, 0.2% BSA and 0.05% azide in either lyophilized (0.1 mg) or liquid form (0.1 mg/0.1 ml or see lot sp. conc on the vial). Reconstitute powder in PBS in 0.1 ml to prepare 1 mg/ml solution. Store at -20oC in suitable aliquots. Stability is ~6-12 months. Do not freeze and thaw.

Suggested conjugate dilutions are 1:1,000-1:10,000 ELISA, 1:1K-1:5K for western, and 1:200-1:1000 (IHC).

Suggested conjugate dilutions are 1:5,000-1:30,000 ELISA, 1:2K-1:10K for western.

**Cat# 20010-3-F, FITC-conjugate**

Purified antibody was coupled to FITC at F/P ratio ~5. The antibody is supplied (100 tests) in PBS, pH 7.4, 1% BSA (see lot sp concn on the vial) or in powder form. Reconstitute powder in PBS. Store at -20oC in suitable aliquots. Stability is ~6-12 months. Do not freeze/thaw.

**Related Material available for ADI**

Poly-Rabbit	60419	Anti-Chicken IgG (Fc), aff pure,
Poly-Rabbit conjugate	60420	Anti-Chicken IgG (Fc)-HRP conjugate
Poly-Rabbit conjugate	60430	Anti-Chicken IgG (Fc)-FITC conjugate
Poly-Rabbit conjugate	60440	Anti-Chicken IgG (Fc)-Biotin conjugate
20365		Streptavidin-Peroxidase (HRP) conjugate
20366		Streptavidin-AP (Alk Phosphatase) conjugate
20367		Streptavidin-Fluorescein (FITC) conjugate
20368		Streptavidin-Rhodamine (TRITC) conjugate
20369		Streptavidin-Phycoerythrin (PE) conjugate

20010-3 140711A

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)