

Product Data Sheet

Camel IgG (Control, non-immune, isotype control, ELISA grade)

Cat# 20031-1 **Size:** □ 1 mg **Form:** □ Liquid □ Lyophilized **Storage:** Store at -20oC

Description

A healthy non-immunized population of Camels were used to collect serum. IgG was prepared by ammonium sulfate, and column chromatography (Sephadex, and DEAE-ion exchange columns). The resulting preparation was judged to be >80% pure by IEP.

Form and Storage

The normal Camel IgG is provided as liquid (1 mg/ml in PBS pH 7.4 + sodium azide as preservative) or in lyophilized form. The product should be **stored at 4°C** for short term and -20oC for long term storage. It is stable for a minimum of 1 year. Do not store diluted solutions.

The **lyophilized products** should be dissolved in PBS, pH 7.4 (1 mg/ml by gentle rocking or vortexing at room temperature. Centrifuge the IgG solution briefly (10,000 g, 15 min) to remove any particulate matter. It should be aliquoted and stored frozen for long term use.

Suggested Uses

This preparation of normal Camels IgG is suitable for coating the ELISA plates or as a non-immune control for ELISA, dot blot, Western or IHC. Recommended testing concn. is approx. 1-10 ug/ml. It is, however, possible that this preparation may be unsuitable for some antigens or cell/tissue extracts and may give slightly higher background. Switching to another preparation from different Camel or vendor may alleviate this problem. The purified IgG may have some residual Camel IgM or IgA.

This product is for in vitro research use only.

MSDS: A Material Safety Data Sheet is not required for this product. The product does not contain any hazardous components above 1% or any carcinogens above 0.1% as defined in 29 CFR 1910.1200, the OSHA Hazard Communication Standard. The product contains natural purified proteins in a buffer containing 0.1% or less sodium azide. MSDS can be consulted for these items. Please note that although the product is defined as not hazardous, it is still advisable to follow prudent laboratory practices when handling laboratory reagents.

This material is sold for research purposes only and is not required to appear on the TSCA inventory. It is not intended for food, diagnostic, drug, household, agricultural or cosmetic use. Its use must be supervised by a technically qualified individual experienced in handling potentially hazardous chemicals.

Related items available

Catalog#	ProdDescription
20031-1	Camel IgG (Control, non-immune, isotype control, ELISA grade)
30835-APC	Anti-Camel IgG-APC conjugate
30835-B	Anti-Camel IgG-Biotin conjugate
30835-HP	Anti-Camel IgG-HRP Conjugate
30835-RPE	Anti-Camel IgG-RPE conjugate
30835-UL	Anti-Camel IgG, aff pure, unlabeled
7820	Camel IgG ELISA Kit, 96 tests, Quantitative
RV-402110-1	Recombivirus Camel Anti-Middle East Respiratory Syndrome Coronavirus (MERS-CoV) Nucleoprotein (NP) antibody (IgG) ELISA kit, 96 tests
RV-402210-1	Recombivirus Camel Anti-Middle East Respiratory Syndrome Coronavirus (MERS-CoV) Spike protein 1 (S1) antibody (IgG) ELISA kit, 96 tests
RV-402310-1	Recombivirus Camel Anti-Middle East Respiratory Syndrome Coronavirus (MERS-CoV) Spike protein 2 (S2) antibody (IgG) ELISA kit, 96 tests
RV-402400-1	Recombivirus Camel Anti-Middle East Respiratory Syndrome Coronavirus (MERS-CoV) Spike protein 2 (S2) antibody (IgG) ELISA kit, 96 tests
RV-402410-1	Recombivirus Camel Anti-Middle East Respiratory Syndrome Coronavirus (MERS-CoV) Spike protein 2 (S2) antibody (IgG) ELISA kit, 96 tests
SP-100526-1	Ac-β- Endorphin, bovine, camel, ovine [Ac-Tyr-Gly-Gly-Phe-Met-Thr-Ser-Glu-Lys-Ser-Gln-Thr-Pro-Leu-Val-Thr-Leu-Phe-Lys-Asn-Ala-Ile-Ile-Lys-Asn-Ala-His-Lys-Lys-Gly-Gln (MW: 3480.1)]
SP-55281-1	B-Endorphin, Camel [H-Tyr-Gly-Gly-Phe-Met-Thr-Ser-Glu-Lys-Ser-Gln-Thr-Pro-Leu-Val-Thr-Leu-Phe-Lys-Asn-Ala-Ile-Ile-Lys-Asn-Ala-His-Lys-Lys-Gly-Gln-OH; MW: 3438.04]
SP-87427-1	β-Endorphin (1-27), camel, bovine, ovine [Tyr-Gly-Gly-Phe-Met-Thr-Ser-Glu-Lys-Ser-Gln-Thr-Pro-Leu-Val-Thr-Leu-Phe-Lys-Asn-Ala-Ile-Ile-Lys-Asn-Ala-His (MW: 2996.51)]
SP-89899-1	Ac-d-Endorphin (bovine, camel, mouse, ovine) [Ac-Tyr-Gly-Gly-Phe-Met-Thr-Ser-Glu-Lys-Ser-Gln-Thr-Pro-Leu-Val-Thr-Leu-Phe-Lys-Asn-Ala-Ile-Ile-Lys-Asn-Ala-His (MW: 3038.54)]
20031-1	150126A

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com