

Product Data Sheet

---

**Purified Mouse IgA control**

<b>Cat #</b> 20102-100	<b>Size:</b> 100 ug			
	<b>Form:</b> Liquid	Lyophilized		<b>Storage:</b> -20oC

---

**Description**

Mouse myeloma cells were grown in mice and ascites fluid collected. The IgA was purified by a combination salt fractionation and column chromatography (Sephadex, and DEAE-ion exchange columns). The resulting preparation was judged to be >95% pure by immuno electrophoresis (IEP) and SDS-PAGE. The purified protein is IgA Kappa.

**Form and Storage**

Mouse IgA is provided as liquid (1 mg/ml in PBS pH 7.4 + 0.05% sodium azide as preservative) or in lyophilized form. The product should be **stored at 4°C** for short term and -20oC for long term storage. It is stable for a minimum of 1 year at -20oC.

The **lyophilized products** should be dissolved in PBS, pH 7.4 to prepare desired concentration by gentle rocking or vortexing at room temperature. Centrifuge the solution briefly (10,000 g, 15 min) to remove any particulate matter. It should be aliquoted and stored frozen for long term use.

**Suggested Uses**

This preparation of mouse IgA is suitable for coating the ELISA plates or as a isotype control for ELISA, dot blot, Western or IHC or other immunoassays.

Antibody concentration must be optimized for each application under defined experimental conditions.

**MSDS:** This product is considered non-hazardous as defined by OSHA (CFR 1910.1200. Nov 25, 1983). It can be disposed of in the drain.

**Country of origin:** USA

*This product is for in vitro research use only.*

---

**Related Items**

Mouse IgA, IgG1, IgG2, IgG3, IgM, IgE Isotype controls

Mouse IgG (IgG1-3), IgM, IgE, and Isotype specific ELISA kits

ELISA kits for measuring antibodies to mouse proteins

rev 71015A

---

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)