

Product Data Sheet

**Rabbit Anti-Goat IgG Antibodies and conjugates**

**Cat# 30221-Cy2** Rabbit Anti-Goat IgG (H+L)-**Cy2** conjugate **Size:** 0.5 ml

**Cat# 30221-Cy3** Rabbit Anti-Goat IgG (H+L)- **Cy3** conjugate **Size:** 0.5 ml

Rabbits were immunized with antigen grade **Goat IgG (H+L)**. Antibodies have been isolated using ammonium sulfate, ion-exchange, and affinity chromatography to remove non-specific antibodies. Specificity of antibody has been tested using IEP, immunodiffusion, and ELISA. The product reacts with Goat IgG with no significant reactivity with other serum protein. It may also react with Goat IgA and IgM proteins. The antibody may recognize other species IgG (H+L or IgA, IgG, IgM etc). Purified anti-Goat IgG (H+L) antibody is supplied as unlabeled or conjugated to Biotin, FITC, Alk. phos., HRP, Rhodamine, phycoerythrin (PE), Texas Red, Cy3, and Cy5 to accommodate applications in various immunological techniques.

**Form and Storage**

**Cat# 30221-Cy2, Cy2-conjugate**

Purified anti-Goat IgG (H+L) antibody was coupled to Cy2 at F/P ratio ~1.7-2.0. The antibody is supplied in PBS, pH 7.4, 0.5% BSA, IgG, and 0.05% azide in either **lyophilized** (0.5 ml) or **liquid** form (0.5 ml at ~0.5 mg/1 ml). Reconstitute powder in PBS in 0.5 ml to prepare 1 mg/ml solution. Store at -20oC in suitable aliquots. Stability is ~6-12 months. Do not freeze and thaw.

**Suggested applications:** Suitable for immunomicroscopy and flow cytometry or FACS analysis as well as other antibody based fluorescent assays

**Suggested conjugate dilutions** are 1:100-1:500 for immunofluorescence. Users must optimize the dilutions for a given technique.

**Absorption Wavelength:** 492 nm

**Emission Wavelength:** 510 nm

**Cat# 30221, Cy3-conjugate**

Purified anti-Goat IgG (H+L) antibody was coupled to Cy3 at F/P ratio ~2.5-2.8. The antibody is supplied in PBS, pH 7.4, 0.5% BSA, IgG, and 0.05% azide in either **lyophilized** (0.5

ml) or **liquid** form (0.5 ml at ~0.5 mg/1 ml). Reconstitute powder in PBS in 0.5 ml to prepare 1 mg/ml solution. Store at -20oC in suitable aliquots. Stability is ~6-12 months. Do not freeze and thaw.

**Suggested applications:** Suitable for immunomicroscopy and flow cytometry or FACS analysis as well as other antibody based fluorescent assays

**Suggested conjugate dilutions** are 1:100-1:500 for immunofluorescence. Users must optimize the dilutions for a given technique.

**Absorption Wavelength:** 550 nm

**Emission Wavelength:** 570 nm

**Recommended Working Dilution for ELISA**

Working dilution for the specific application should be determined by the investigator to obtain the best conditions. Working solution should be prepared immediately before use and diluted solution should be discarded.

**Origin:** All material used in various products are derived from US sources and USDA inspected facilities.

**MSDS:** A Material Safety Data Sheet is not required for this product. The product does not contain any hazardous components above 1% or any carcinogens above 0.1% as defined in 29 CFR 1910.1200, the OSHA Hazard Communication Standard. This material is not required to appear on the TSCA inventory.

*This product is for in vitro research use only.*

**Related Material available for ADI**

**Goat IgG, IgM ELISA kits**

Anti-Goat IgG (H+L), IgG, IgG (H+L), IgG (H+L) antibodies and conjugates

30221-Cy2-Cy3

110222A

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)