

Product Data Sheet

Donkey Anti-Goat IgG Antibodies and conjugates (preadsorbed with multiple species)

Cat# 30224-200

Donkey Anti-Goat IgG-FITC conjugate (adsorbed)

Size: 200 ug/0.5 ml

Description

Donkeys were immunized with antigen grade **Goat IgG (H+L)**. Antibodies have been isolated using ammonium sulfate, ion-exchange, and affinity chromatography to remove non-specific antibodies. Specificity of antibody has been tested using IEP, immunodiffusion, and ELISA. The product reacts with Goat IgG, IgA and IgM with no significant reactivity with other Goat serum proteins. **Antibodies were further adsorbed with serum proteins from several species (bovine, Chicken, Guinea Pig, Hamster, Horse, Mouse, Rabbit, and Rat) to minimize reactivity with these species.** The antibody may recognize goat IgG (H+L or IgA, IgG, IgM etc). Purified anti-Goat IgG (H+L) antibody is supplied as unlabeled or conjugated to Biotin, FITC, Alk. phos., HRP, Rhodamine, phycoerythrin (PE), Texas Red, Cy3, Cy5, and various IRDye (700DX, 800, and 800CW) to accommodate applications in various immunological techniques.

Form and Storage

Cat# 30224-200, FITC-conjugate

Purified anti-Goat IgG (adsorbed) antibody was coupled to FITC at F/P ratio ~3:7. The antibody is supplied in PBS, pH 7.4, 0.2% BSA and 0.05% azide in either **lyophilized** (200 ug) or **liquid** form (200 ug/0.5 ml). Reconstitute powder in PBS in 0.5 ml to prepare 1 mg/ml solution. Store at -20°C in suitable aliquots. Stability is ~6-12 months. Do not freeze and thaw.

Suggested conjugate dilutions are 1:200-1:2000 for immunofluorescence.

Absorption Wavelength: 495 nm

Emission Wavelength: 528 nm

Recommended Working Dilution

Working dilution for the specific application should be determined by the investigator to obtain the best conditions. Working solution should be prepared immediately before use and diluted solution should be discarded.

Origin: All material used in various products are derived from US sources and USDA inspected facilities.

MSDS: A Material Safety Data Sheet is not required for this product. The product does not contain any hazardous components above 1% or any carcinogens above 0.1% as defined in 29 CFR 1910.1200, the OSHA Hazard Communication Standard. This material is not required to appear on the TSCA inventory.

This product is for in vitro research use only.

Related Material available for ADI

Goat IgG, IgM ELISA kits

Anti-Goat IgG (H+L), antibodies and conjugates

30224-200

90511A

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com