

Product Data Sheet

**Goat Anti-Deer (Cervid) IgG Antibodies and conjugates**

<b>Cat#</b> 30812-A	Goat Anti-Deer IgG (H+L) (Cervid family: Deer, Elk, Moose), unlabeled	<b>Size:</b> 100 ug
<b>Cat#</b> 30812-HP	Goat Anti-Deer IgG (H+L) (Cervid family: Deer, Elk, Moose)-HRP conjugate	<b>Size:</b> 100 ul
<b>Cat#</b> 30812-BT	Goat Anti-Deer IgG (H+L) (Cervid family: Deer, Elk, Moose)-Biotin conjugate	<b>Size:</b> 100 ul

Goats were immunized with antigen grade **Deer IgG (H+L)**. Antibodies have been isolated using ammonium sulfate, ion-exchange, and affinity chromatography to remove non-specific antibodies. Specificity of antibody has been tested using IEP, immunodiffusion, and ELISA. The product reacts with Deer IgG, IgA and IgM with no significant reactivity with other serum proteins. The antibody also recognize other related species IgG (Elk and moose or cervid family). Purified anti-Deer IgG (H+L) antibody is supplied as unlabeled or conjugated to Biotin, FITC, AP or HRP.

**Form and Storage**

**Cat# 30812-A, Anti-Deer IgG (Cervid), unlabeled**

The antibody is supplied in PBS, pH 7.4, and 0.05% azide in either **liquid** form or lyophilized in the same buffer. Reconstitute powder in PBS in 0.1 ml to prepare stock solution. Store at -20oC in suitable aliquots. Stability is ~6-12 months. Do not freeze and thaw.

Suggested conjugate dilutions are 1-10 for ELISA western or IHC. Use anti-Goat IgG or conjugates as secondary antibodies.

**Cat# 30812-HP, Anti-Deer IgG (Cervid) IgG-HRP-conjugate**

Purified antibody was coupled to HRP (RZ>3.0) using periodate method. The molar enzyme to protein (E/P) ratio = 4.0. The antibody is supplied in stabilizing buffer, 0.1% proclin-300 as preservative in **liquid** or **lyophilized** or form (or see lot sp. conc on the vial). Reconstitute powder in PBS to prepare stock solutions. Store at 4oC in suitable aliquots. Stability is ~6-12 months. Do not freeze and thaw.

Suggested conjugate dilutions are 1:1,000-1:10,000 ELISA, 1:1000-1:5000 for western, and 1:200-1:1000 (IHC).

**Cat# 30812-BT, Anti-Deer IgG (Cervid) IgG-Biotin-conjugate**

Purified anti-Deer (Cervid) IgG (H+L) antibody was coupled to Biotin using Biotinamidocaproate N-Hydroxysuccinimide Ester (BAC) at F/P ratio ~10-20:1. The antibody is supplied in PBS, pH 7.4, 0.2% BSA and 0.05% azide in either **lyophilized** or **liquid** form (see lot sp. conc on the vial). Reconstitute powder in PBS in 0.1 ml to prepare stock solution. Store at -20oC in suitable aliquots. Stability is ~6-12 months. Do not freeze and thaw.

Suggested conjugate dilutions are 1:5,000-1:30,000 ELISA, 1:2K-1:10K for western.

**Recommended Working Dilution for ELISA**

Working dilution for the specific application should be determined by the investigator to obtain the best conditions. Working solution should be prepared immediately before use and diluted solution should be discarded.

**Related Material available for ADI**

Catalog#	ProdDescription
30811-A	Anti-Deer IgG/Y (H+L chain) (Cervid family: Deer, Elk, Moose) unconjugated, aff pure IgG,
30811-BT	Anti-Deer IgG/Y (H+L chain) (Cervid family: Deer, Elk, Moose)-Biotin Conjugate
30811-HP	Anti-Deer IgG/Y (H+L chain) (Cervid family: Deer, Elk, Moose)-HRP Conjugate
30812-A	Anti-Deer IgG (H+L chain) (Cervid family: Deer, Elk, Moose) unconjugated, aff pure IgG,
30812-BT	Anti-Deer IgG (H+L chain) (Cervid family: Deer, Elk, Moose)-Biotin Conjugate
30812-HP	Anti-Deer IgG (H+L chain) (Cervid family: Deer, Elk, Moose)-HRP Conjugate
30813-A	Anti-Raccoon IgG (H+L chain) IgG unconjugated, aff pure
30813-HP	Anti-Raccoon IgG (H+L chain) IgG-HRP conjugate
30814-A	Anti-Alligator IgG (H+L chain) IgG unconjugated, aff pure
30814-HP	Anti-Alligator IgG (H+L chain) IgG-HRP Conjugate
30815-A	Anti-Opossum IgG (H+L chain) IgG, unconjugated, aff pure
30815-HP	Anti-Opossum IgG (H+L chain) IgG-HRP Conjugate
30816-A	Anti-Bat IgG (H+L chain) IgG, unconjugated, aff pure
30816-HP	Anti-Bat IgG (H+L chain) IgG-HRP Conjugate
30817-A	Anti-Llama IgG (H+L chain) IgG, unconjugated, aff pure
30817-HP	Anti-Llama IgG (H+L chain) IgG-HRP Conjugate
30818-A	Anti-Macaw IgG (H+L chain) IgG, unconjugated, aff pure
30818-HP	Anti-Macaw IgG (H+L chain) IgG-HRP Conjugate
30819-A	Anti-White whale IgG (H+L chain) IgG, unconjugated,
30819-HP	Anti-White whale IgG (H+L chain) IgG-HRP conjugate
30820-A	Anti-Dolphin (Bottlenose) IgG (H+L chain) IgG,
30820-HP	Anti-Dolphin (Bottlenose) IgG (H+L chain) IgG-HRP
30821-A	Anti-Dog IgG (Gray Fox/red fox/coyote) IgG, aff pure,
30821-HP	Anti-Dog IgG (Gray Fox/red fox/coyote) IgG-HRP Conjugate

**Western blot recycling kit** (Use the same blot/strip to probe with multiple antibodies) **Strip Antibodies in 5 minutes** (no heating or mercaptoethanol)

30812 121114A

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)