

Product data sheet

Anti-Mouse IgG (Fc specific)-Antibodies and conjugates

Cat# 40151-UL	Goat Anti-Mouse IgG (Fc), unlabeled	Size: 1 mg
Cat# 40151-HRP	Goat Anti-Mouse IgG (Fc)-HRP conjugate	Size: 1 ml
Cat# 40151-FITC	Goat Anti-Mouse IgG (Fc)-FITC labeled	Size: 1 ml
Cat# 40151-BT	Goat Anti-Mouse IgG (Fc)-Biotin labeled	Size: 1 ml

Goats were immunized with antigen grade mouse IgG. Antibodies have been isolated using ammonium sulfate, ion-exchange, and affinity chromatography. Antibodies were made specific for IgG-Fc by removing any antibodies crossreactive with light chains and other Ig's (IgM, IgA etc). Specificity has been tested using IEP, immunodiffusion, and ELISA. Purified antibodies react with mouse IgG with minimal reactivity with mouse IgA and IgM. The antibody may recognize other species IgG's that have common gamma chains. However, no antibody is detected to other serum proteins. Purified antibodies are supplied as unlabeled, HRP-, Biotin, FITC, and AP conjugates.

Form and Storage

Cat# 40151-UL, unlabeled

The antibody is supplied in PBS, pH 7.4, and 0.05% azide in either **lyophilized** (1.0 mg) or **liquid** form (1 mg/ml). Reconstitute powder in PBS in 1 ml to prepare 1 mg/ml solution. Store at -20oC in suitable aliquots. Stability is ~6-12 months. Do not freeze and thaw.

Cat# 40151-HRP, HRP-conjugate

Purified antibody was coupled to HRP (RZ>3.0) using periodate method. The molar enzyme to protein (E/P) ratio = 4.0. The antibody is supplied in stabilizing buffer, 0.1% prolcin-300 as preservative in either **lyophilized** (1 ml) or **liquid** form (1 ml). Reconstitute powder in PBS in 1 ml. Store at 4oC in suitable aliquots. Stability is ~6-12 months. Do not freeze and thaw.

Suggested conjugate dilutions are 1:1,000-1:10,000 ELISA, 1:1K-1:5K for western, and 1:200-1:1000 (IHC).

Cat# 40151-FITC, FITC-conjugate

Purified anti-Mouse IgG antibody was coupled to FITC at F/P ratio ~4-5:1. The antibody is supplied in PBS, pH 7.4, 0.2% BSA and 0.05% azide in either **lyophilized** or **liquid** form (or see lot sp. conc on the vial). Reconstitute powder in PBS in 1 ml to prepare 1 mg/ml solution. Store at -20oC in suitable aliquots. Stability is ~6-12 months. Do not freeze and thaw.

Suggested conjugate dilutions are 1:200-1:2000 for immunofluorescence.

Absorption Wavelength: 495 nm

Emission Wavelength: 528 nm

Cat# 40151-BTN, Biotin-conjugate

Purified antibody was coupled to Biotin using GSI Biotinylation kit # 80300 at F/P ratio ~10-20:1. The antibody is supplied in PBS, pH 7.4, 0.2% BSA and 0.05% azide in either **lyophilized** or **liquid** form (or see lot sp. conc on the vial). Reconstitute powder in PBS in 1 ml to prepare stock solution. Store at -20oC in suitable aliquots. Stability is ~6-12 months. Do not freeze and thaw.

Suggested conjugate dilutions are 1:5,000-1:30,000 ELISA.

Recommended Working Dilution for ELISA

Working dilution for the specific application should be determined by the investigator to obtain the best conditions. Working solution should be prepared immediately before use and diluted solution should be discarded.

All products are for In vitro research use only.

Related Material available for GSI

Mouse IgG (IgG1-3), IgA, IgM, and IgE ELISA kits

Anti-Mouse IgG (gamma-chain sp), IgG, IgM, IgE antibodies and conjugates

ELISA kits for the detection of mouse Antibodies

40151

130429A