

Product Specification Sheet

Tris Buffered Saline (TBS, 10X)

<input type="checkbox"/> Cat. # 80083-W	Tris Buffer pH 7.4, concentrate (TB 10X; 50 mM Tris, pH 7.4, No salts)	SIZE: 1 L
<input type="checkbox"/> Cat. # 80083	Tris Buffered Saline (TBS, pH 7.4) 10X concentrate (50 mM Tris/150mM NaCl final)	SIZE: 1 L
<input type="checkbox"/> Cat. # 80083-K	Tris Buffer-KCl (TBK, pH 7.4) 10X concentrate (50 mM Tris/150mM NaCl/5 mM KCl final)	SIZE: 1 L

#80083-W (TB, 10X)

Buffer Composition: This buffer consists of 50 mM Tris HCl, pH of 7.4 (final concn). NO SALTS.

#80083 (TBS, 10X)

Buffer Composition: This buffer consists of 50 mM Tris HCl and 150 mM Sodium Chloride, pH of 7.4 (final concn).

#80083-K (TBK, 10X)

Buffer Composition: This buffer consists of 50 mM M Tris HCl and 5mM KCl, 150 mM NaCl, pH of 7.4 (final concn).

Stabilizer or Preservatives: None. If desired, sodium azide 0.1% or other preservatives may be added.

Application(s): Suitable for ELISA or Western or molecular immunology, or other assays when applicable.

Recommended Dilution(s): This product is a 10X concentrated stock solution and should be diluted appropriately with distilled, deionized water (or equivalent) to its final working concentration. Typically, dilute 1:10 for ELISA or Western buffers (final concn 50mM Tris, 150 mM NaCl, pH 7.5).

Storage: Store container at 4oC or long term storage. Some salts may precipitate out of solution at lower temperature. Allow buffer to equilibrate to room temperature (18° to 26° C) to restore solubility of some salts. Expiration date is six (6) months from date of opening container.

Sterilization: This product is not sterile. If necessary, it can be sterile filtered through a Millipore 0.22-micron filter into clean, pre-sterilized containers.

*This product is for In vitro research use only.

Related material available from ADI

Catalog#	ProdDescription
80050	ELISA Plate Coating buffer concentrate (10X)
80060	ELISA Plate Blocking Buffer concentrate (10X) BSA-based
80062	ELISA Plate Blocking Buffer concentrate (10X) milk-based
80070	Antibody and Conjugate Diluent for ELISA concentrate (10X)
80080	Wash buffer concentrate (20X) for ELISA
80081	Phosphate Buffered Saline (PBS, pH 7.4) concentrate (20X)
80082	Tris Buffered Saline (TBS, pH 7.4) concentrate (10X)
80091	TMB substrate (1-component) for ELISA
80100	Stop solution for TMB substrate (ELISA) concentrate (10X)
80110	Single solution TMB substrate (membrane) solution for western blot

80083-Tris-Buffer

160303SV

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com